HAEMATOKRIT 210

Ⓓ	Bedienungsanleitung	5
Ⓔ	Operating Instructions	18
Ⓕ	Mode d'emploi	31
Ⓘ	Istruzioni per l'uso	44


Fig. 1

**EG-Konformitätserklärung
EC Conformity Declaration
Déclaration de conformité CE
Dichiarazione di conformità alle norme CEE**

Andreas Hettich GmbH & Co. KG • Föhrenstraße 12 • D-78532 Tuttlingen • Germany

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Contents

1	Use according to specification	20
2	Residual risks	20
3	Notes on safety	20
4	Explanations of symbols and references	21
4.1	Symbol used on the machine	21
4.2	Symbols used in the operating instructions	21
5	Delivery checklist	21
6	Space requirement	22
7	Connection to the mains	22
8	Commissioning	22
9	Opening and closing the lid	22
9.1	Opening the lid	22
9.2	Closing the lid	22
10	Control panel	23
10.1	Speed area	23
10.2	Rotation area	23
10.3	Time area	23
10.4	Key area	23
11	Adjustable parameters	24
12	Preselecting centrifuging parameters or changing them during operation	24
12.1	Speed	24
12.2	- of denser substances	24
12.3	Run time	24
12.3.1	Continuous operation	24
12.3.2	Pulsed operation	24
12.4	Brake adjustment	25
12.5	Start centrifuging run	25
12.6	End centrifuging run	25
13	Relative centrifugal force (RCF)	25
14	Haematocrit measurement	26
14.1	Haematocrit capillaries - Overview	26
14.2	Filling of the haematocrit capillaries	26
14.3	Sealing and centrifugation of the haematocrit capillaries	26
14.3.1	Standard capillary tubes, Order No.: 2074, and mylar-reinforced capillary tubes, Order No.: 1072	26
14.3.2	Self-sealing capillary tubes (Order No. 1071)	27
14.4	Loading of rotor	27
14.5	Determination of the centrifugation time	27
14.6	Reading the haematocrit value	27
15	Emergency release	27

16	Maintenance and servicing.....	28
16.1	Cleaning the centrifuge	28
16.2	Cleaning and disinfecting the haematocrit rotor	28
16.3	Procedure in case of broken haematocrit capillaries.....	28
16.4	Rotor removal / installation.....	29
17	Faults	30
18	Change mains input fuse.....	30
19	Acceptance of the centrifuges for repair	30
20	Anhang / Appendix	57
20.1	Technische Daten / Technical specification	57
20.2	Rotoren und Zubehör / Rotors and accessories.....	58

1 Use according to specification

The machine on hand is a centrifuge that is exclusively conceived for the separation of substances or respectively substance mixtures with a density of max. 1.2 kg/dm³ and which is therefore only designated for this application. A different use or application over and above this is deemed not in accordance with the specifications. The company Andreas Hettich GmbH & Co. KG undertakes no liability for damages resulting therefrom.

Belonging to the application according to specification is also the observance of all references contained in the Instruction Manual and compliance with the inspection and maintenance works.

2 Residual risks

The machine is constructed according to the state of the art and the recognized technical safety regulations. However during your application of the device hazards towards life and limb of the user or third parties or respectively impairments to the device or other material assets can arise. The machine is only to be used for the specified applications and only in an impeccable technical safety condition.

Disturbances that can interfere with the safety are to be immediately rectified.

3 Notes on safety

- **Before the initial operation of your centrifuge you should read and pay attention to the operating instructions.**
Only personnel that has read and understood the operating instructions are allowed to operate the device.
- Along with the operating instructions and the legal regulations on accident prevention, you should also follow the recognised professional regulations for working in a safe and professional manner. These operating instructions should be read in conjunction with any other instructions concerning accident prevention and environmental protection based on the national regulations of the country where the device is to be used.
- This centrifuge is a state-of-the-art piece of equipment which is extremely safe to operate.
 - However, it can lead to danger for users or others if used by untrained staff, in an inappropriate way or for a purpose other than that it was designed for.
- The centrifuge should be installed on a good, stable base.
- Before using the centrifuge absolutely check the rotor for firm placement.
- Ensure that no persons and / or dangerous materials are located within a safety zone of 300 mm around the equipment when a centrifuge is running.
- The centrifuge must not be moved or knocked during operation.
- In case of fault or emergency release, never touch the rotor before it has stopped turning.
- To avoid damage due to condensate, when changing from a cold to a warm room the centrifuge must either run hot in the cold room for 30 minutes, or warm up for 3 hours in the warm room, before connecting to the mains.
- Load centrifuge rotor evenly. All positions on rotor must be filled.
- Do not fill centrifuge containers inside the centrifuge.
- Centrifuge containers must not be filled beyond the capacity specified by the manufacturer.
- Standard centrifuge containers of glass will not stand RCF values exceeding 4000 (DIN 58970, pg. 2)
- Centrifuge containers may only be centrifuged with accessories (reducing adapters, frames, suspensions, etc.) authorised by the manufacturer (see section "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories").
- When centrifuging with maxim revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm³.
- The centrifuge may only be operated when the balance is within the bounds of acceptability.
- The centrifuge must not be operated in areas subject to danger of explosions.
- The centrifuge must not be used with:
 - inflammable or explosive materials
 - materials that react with one another producing a lot of energy.
- If users have to centrifuge hazardous materials or compounds contaminated with toxic, radioactive or pathogenic micro-organisms, they must take appropriate measures.
Without additional proceedings (like an additional bioseal between bucket and lid of bucket or angle rotor with a special bioseal between rotor and lid) a centrifuge is not a biosafety system in accordance to the regulation EN 61010-2-020. In the case of material belonging to risk group II (see the World Health Organisation's "Laboratory Biosafety Manual") they should employ a biosafety system. Under this system small drips and aerosols are prevented from escaping by a bioseal (packing ring) located between the hanger and the lid. Centrifuge containers with special screw caps, as obtainable through trade suppliers, can also be used for hazardous substances.
In the case of materials from the higher risk groups greater safety provision is required than the arrangements described above. In a biosafety system, centrifuge containers with special screw caps must be used.
- For further details of available biosafety systems see section "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories". If in doubt, you should obtain relevant information from the manufacturer.
- The centrifuge must not be operated with highly corrosive substances which could impair the mechanical integrity of rotors, hangers and accessories.

- Any rotors, hangers or accessories showing clear signs of corrosion or mechanical defects must not be used for centrifuging.
- Repairs must only be carried out by personnel authorised to do so by the manufacturer.
- Only original spare parts and original accessories licensed by the Hettich company are allowed to be utilised.
- This centrifuge is classified in Germany as a Group 3 device according to the *Medizinische Geräteverordnung MedGV* (the regulations on medical equipment).
- It conforms to safety regulations based on: IEC 61010-1, IEC 61010-2-020
- The safe operation and reliability of the centrifuge can only be guaranteed if:
 - the centrifuge is operated in accordance with the operating instructions,
 - the electrical installation on the site where the centrifuge is installed conforms to the demands of IEC stipulations,
 - prescribed tests to BGV A1, BGR 261 are carried out by an expert.
- With centrifuges for robotic use please pay attention to the notes of the key operated switch.

No claim under guarantee will be considered by the manufacturer unless the above instructions have been adhered to.

4 Explanations of symbols and references

4.1 Symbol used on the machine

Attention, general hazard area.

Before using the centrifuge implicitly read the operating instructions and pay attention to the safety relevant references!

4.2 Symbols used in the operating instructions

This symbol refers to important circumstances.

This symbol refers to safety relevant warnings and indicates possibly dangerous situations.

The non-adherence to these warnings can lead to material damage and injury to personal.

5 Delivery checklist

The following items and accessories are delivered with the centrifuge:

- | | |
|---|--------------------------------------|
| 1 | Connecting cable |
| 2 | Fuses |
| 1 | Notes on moving the equipment safely |
| 1 | Operating instructions |
| 1 | Hexagonal pin-type spanner |

The rotor(s) and associated accessories are included in the delivery in the quantity.

6 Space requirement

- According to the laboratory instrument standards EN 61010-2-20 an emergency switch to separate power supply in the event of a failure must be installed in the building electrical system.
This switch has to be placed remote from the centrifuge, preferred outside of the room in which the centrifuge is installed or near by the exit of this room.
- The necessary space requirement can be found under Dimensions in the Technical data chapter.
- The centrifuge must be set up in a suitable place, so that it is stable. During set up the required safety area of 300 mm around the centrifuge, in accordance with IEC 1010-2-2, must be observed.

Persons and hazardous materials must not be located in the safety area whilst the centrifuge is in operation.

- Substances can be discharged through the ventilation opening on the rear side of the centrifuge. The device is to be placed in such a way that the air stream is not directed towards people.
- Do not place any object in front of the ventiduct.
 - Keep a ventilation area of 300 mm around the ventiduct.

7 Connection to the mains

- Check whether the supply voltage, supply frequency and on-site mains fuse agree with the specification on the nameplate. The nameplate is located on the back of the centrifuge.
- Make sure that the mains switch is in the “0” position.
- The centrifuge must be connected to a standard mains socket using the power supply cable provided.

8 Commissioning

- Check that the centrifuge has been properly set up and that the electrical connections are correct (see Connection to the Mains and Space Requirement).
- Switch the mains switch “ON”, switch position “I”.
After a short time the control panel will switch on and the last set parameters will appear in the displays.
- When the  symbol has lit up, the lid can be opened.
- Open the lid and remove the transport protection (see Transport Protection information sheet).

9 Opening and closing the lid

9.1 Opening the lid

- When the  symbol has lit up, the lid can be unlocked and opened.
- Unlock the lid by pressing the  key and swivel the lid upwards.
- The  symbol appears.

The lid can only be opened if the centrifuge is switched on and the rotor is stationary.
If this is not possible, see chapter “Emergency release”.

9.2 Closing the lid

- Close the lid by swinging the lid down and pressing lightly. If the symbol  lights up in the rotation display the lid has been correctly locked.


Do not bang the lid shut.

10 Control panel

The control panel is divided into four areas.

- Speed
- Rotation display
- Time
- Keys

10.1 Speed area

The speed area contains the speed display which indicates the preselected speed or the ACTUAL speed during operation.

The speed can be preselected or changed during operation with the ▲ and ▼ arrow keys. At the end of operation, the display flashes in seconds until the lid is opened or a key is activated.

10.2 Rotation area

The rotation area contains the rotation display, which illuminates and rotates anti-clockwise whilst the rotor is turning.

When the ⏪ symbol lights up, the rotor is stationary and the lid can be opened.

When the lid is unlocked, the ⏪ symbol lights up.

10.3 Time area

This area contains the running time display and the ▲ and ▼ arrows, with which the switch-on time can be preselected or changed during operation.

When stationary, the preselected switch-on time is displayed, as well as the remaining run time during operation. The decimal point flashes each second during operation.

At the end of operation, the display flashes each second until the lid is opened or a key is activated.

10.4 Key area

Start key 	To start a run with the preselected parameters.
Stop key 	The STOP key allows a centrifuging run to be stopped at any time.
Pulse key 	Brief centrifuging whilst the key is pressed down. The run time is displayed in seconds in the run time display. After 60 seconds, the run time is displayed in minutes.

11 Adjustable parameters

Speed	<ul style="list-style-type: none"> • min. speed 500 rpm • max. speed 13000 rpm • in steps of 100 rpm Adjustable up to max. speed of the rotor used.
Time	<ul style="list-style-type: none"> • min. preset time 1 min • max. preset time 99 min • in steps of 1 min or
"—" 	<ul style="list-style-type: none"> • continuous operation or
Braking force	<ul style="list-style-type: none"> • full braking effect • low braking effect

12 Preselecting centrifuging parameters or changing them during operation

12.1 Speed

Each time the  and  arrow keys are activated, the speed is changed by 100 rpm. If an arrow key is kept pressed down for longer, the speed is changed at an increasing rate. If the speed is changed during operation, this becomes effective immediately, i.e. the speed is adjusted.

12.2 - of denser substances

The rotors are designed to centrifuge substances up to a maximum mean homogenous density of 1.2 kg/dm³ when rotating at the stated speed. Denser substances must be centrifuged at lower speed.

The permissible speed can be calculated using the following formula:

$$\text{Reduced speed (n}_{\text{red}}) = \sqrt{\frac{1.2}{\text{Greater density}}} \times \text{Rated speed}$$

e.g.: RPM 4000, density 1.6 kg/dm³

$$n_{\text{red}} = \sqrt{\frac{1.2}{1.6}} \times 4000 = 3464 \text{ RPM}$$

If in doubt you should obtain clarification from the manufacturer.


12.3 Run time

Each time the  and  arrow keys are activated, the run time is changed by 1 minute. If an arrow key is kept pressed down for longer, the run time is changed in minute steps, at an increasing rate. If the run time is changed during operation, this becomes effective immediately, i.e. the remaining run time is adjusted.

If an arrow key is activated during operation, the run time or speed is adjusted accordingly.

12.3.1 Continuous operation

The centrifuge can be operated in continuous operation. Continuous operation is set by pressing the  arrow key until "—" appears in the run time display. Operation is only stopped by activating the  key.



12.3.2 Pulsed operation


For brief centrifuging. The rotor turns at the preselected speed as long as the  key is pressed down.

12.4 Brake adjustment

With this centrifuge, the braking effect can be set to normal or low.

This can be set before a run, as follows:

1. Switch off the mains switch
2. Simultaneously press the pulse key  and the arrow key  in the speed area and switch on the mains switch. Release the keys again.
3. If required press the arrow button  in the speed field until the depicted display appears. In the speed field the machine version as set by the factory and in the time field the brake setting.

4. Set the required value, "1" or "0", in the time area with the  and  arrow keys.
5. Confirm the set value with the **STOP** key.
See Rotors and Accessories in the Braking Times chapter.

12.5 Start centrifuging run

After setting the centrifuging parameters, start the centrifuge by pressing the **START** key.

12.6 End centrifuging run

A centrifuging run is ended by pressing the **STOP** key or after expiry of the set run time. Once the rotor has come to a stop the lid lock is briefly triggered and the display flashes. The lid can now be opened.

13 Relative centrifugal force (RCF)

The relative centrifugal force (RCF) is given as a multiple of the acceleration of gravity (g). It is a unit-free value and serves to compare the separation and sedimentation performance.

These values are calculated using the formula below:

$$\text{RCF} = \left(\frac{\text{RPM}}{1000} \right)^2 \times r \times 1,118 \quad \Rightarrow \quad \text{RPM} = \sqrt{\frac{\text{RCF}}{r \times 1,118}} \times 1000$$

RCF = relative centrifugal force

RPM = rotational speed (revolutions per minute)

r = centrifugal radius in mm = distance from the centre of the turning axis to the bottom of the centrifuge.

For more on the centrifugal radius see the chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

The relative centrifugal force (RCF) stands in relation to the revolutions per minute and the centrifugal radius.

14 Haematocrit measurement

14.1 Haematocrit capillaries - Overview

Standard capillary tubes, Order No.: 2074

The classical heparinized haematocrit capillary tube made of glass. Tried and tested for years.

Mylar-reinforced capillary tubes, Order No.: 1072

The heparinized glass capillary is coated with plastic. Should such a capillary break the glass splinters are retained in this coating. This means that the risk of injury for the user is reduced.

Self-sealing capillary tubes, Order No. 1071

These capillary tubes, which are also heparinized, no longer need to be sealed. They are already equipped with a seal. At first, the seal material is air-permeable so that the capillary can be filled. The seal will only become impermeable when it comes into contact with blood. Therefore, there will be no leaks during centrifugation.

14.2 Filling of the haematocrit capillaries

The capillaries have different volumes according to type.

The following are the exact data:

Capillary:	Standard	Mylar-reinforced	Self-sealing
Order No.:	2074	1072	1071
Total volume without seal (µl):	49.6	74.9	74.9
Internal diameter (mm):	0.92	1.12	1.12
Overall length (mm)	75	75	75

Therefore, more blood is needed for the capillaries No. 1072 and 1071 to be filled. The dimensions, however, are still in compliance with the requirements of DIN 12846 (standard for haematocrit capillaries) and the accuracy of the haematocrit values is not affected.

For the haematocrit measurement the capillaries are always filled up to $\approx 75\%$:

Capillary:	Standard	Mylar-reinforced	Self-sealing
Order No.:	2074	1072	1071
Volume when filled up to 75%(µl)	37.2	56.2	56.2

14.3 Sealing and centrifugation of the haematocrit capillaries

14.3.1 Standard capillary tubes, Order No.: 2074, and mylar-reinforced capillary tubes, Order No.: 1072

1. When filling the capillary ensure that the end of the capillary opposite the filler opening remains dry.
2. Seal at the dry end of the capillary:
Insert the capillary tube vertically into the seal up to the point where the end of the capillary touches the bottom of the seal plate.
Slightly tilt the capillary tube to the side and remove it from the seal.
3. Place the capillary tubes with the lute side facing outwards in the haematocrit rotor and apply the rotor lid and click into place.
4. Centrifuge.
5. Determine the haematocrit value.

14.3.2 Self-sealing capillary tubes (Order No. 1071)

1. Fill the capillary tube.
2. Let the column of blood settle in the direction of the seal.

The blood must touch the seal!

3. Tap the seal end of the capillary three times on the laboratory table!
This serves to intensify the contact with the seal material and to enhance the sealing process. Only in this manner the capillaries will be sealed off!
4. Place the capillary tubes with the lute side facing outwards in the haematocrit rotor and apply the rotor lid and click into place.
5. Centrifuge.
6. Determine the haematocrit value.

14.4 Loading of rotor

- Load opposite positions equally.
 - Put lid onto rotor and click into place.
 - Place the lid and lightly press down the front edge of the lid.
 - Turn the handle on the front panel to the right.
- If the symbol  is illuminated in the rotation display, the lid has been closed properly.

14.5 Determination of the centrifugation time

The following approximation formula is used for the determination of the centrifugation time:

$$\frac{100000}{\text{RCF}} = \text{time [min]}$$

Example:

Given an RCF of 16060 $\frac{100000}{16060} \approx 7 \text{ min}$ of centrifugation time.

For the calculation of the RCF see chapter "Relative Centrifugal Force (RCF)".

14.6 Reading the haematocrit value

Perform the set-up as follows:

1. hold the haematocrit rotor.
2. move the evaluation disk (lid of the haematocrit rotor) until the 100%-line of the evaluation disk is in line with the upper margin of the liquid column.
3. Read the haematocrit value at the upper margin of the erythrocyte column.

Repeat steps 1 to 3 until all capillaries have been analysed.

15 Emergency release

The lid cannot be opened in the event of a power failure. An emergency release has to be executed by hand.

For emergency release disconnect the centrifuge from the mains.
Open the lid only during rotor standstill.

Please turn to the rear side of the cover sheet.

- Position the centrifuge on the table sufficiently far forward so that the hole (Fig. 1, A) for the emergency release on the base of the casing can be reached.
- Insert the hexagonal pin driver into this hole, push it upwards and at the same time rotate the lid-lock handle to the left.
- Open the lid.

16 Maintenance and servicing

Pull the mains plug before cleaning.

Before any other cleaning or decontamination process other than that recommended by the manufacturer is applied, the user has to check with the manufacturer that the planned process does not damage the device.

- Cleaning agents and disinfectants which lie in the pH range 5 – 8 are to be utilised. Alkaline cleaning agents with a pH value > 8 are to be avoided.
- In order to prevent appearances of corrosion through cleaning agents or disinfectants, the application guide from the manufacturer of the cleaning agent or disinfectant are absolutely to be heeded.

Haematocrit rotor, rotor lid, padded shells, and sealing rings are not autoclavable.

16.1 Cleaning the centrifuge

- Regularly clean the centrifuge housing and the centrifugal chamber; and if necessary, clean with soap or a mild cleaning agent and water. For one thing, this services purposes of hygiene, and it also prevents corrosion through adhering impurities.
- In the event of condensation water formation, dry the centrifugal chamber by wiping out with an absorbent cloth.
- If infectious materials penetrates into the centrifugal chamber this is to be disinfected immediately. For surface disinfection we recommend Bacillol Plus from the company Bode Chemicals in Hamburg. After use the disinfection medium Bacillol Plus must be thoroughly removed with water.
- Lightly grease the rubber seal of the centrifugal chamber after every cleaning.

16.2 Cleaning and disinfecting the haematocrit rotor

If despite all precautions a leakage of haematocrit capillaries during centrifugation should occur, rotor and lid must be cleaned according to the following instructions:

1. Carefully remove padded shells or sealing ring from the haematocrit rotor.
2. Soak haematocrit rotor, lid and padded shells or sealing ring in cold water, until the blood residues are completely dissolved.
3. Perform a suitable disinfection. For disinfection we recommend Helipur H plus N from the company B. Braun Melsungen. For the removal of radioactive material we recommend deco neutracon from the company Decon Laboratories Limited. The application of disinfection substances on the lid is to be avoided as the inscription can dissolve.
4. Rinse the components in cold water and dry them.
5. They can now be replaced into the centrifuge.

16.3 Procedure in case of broken haematocrit capillaries

Risk of injury!


Best work with two pairs of gloves worn on top of each other. Wear safety goggles and mask, as glass splinters could be blown off!

1. Carefully remove lid of haematocrit rotor.
2. Remove larger pieces of broken capillaries using pincers.
3. Remove rotor. Carefully and slowly remove the affected padded shells or the sealing ring with pincers and exchange them.
4. Clean and disinfect as described above.

16.4 Rotor removal / installation

Removal of haematocrit rotor:

1. Open the lid and disconnect the centrifuge from the power supply.
2. Remove the rotor lid by pressing the press lock.
3. Unscrew the fixing nut with the supplied spanner.
4. Pull the haematocrit rotor off the motor shaft.

Installation of the haematocrit rotor:

1. Insert the haematocrit rotor vertically onto the motor shaft. The motor shaft dog has to fit in the rotor slot.
2. Tighten the fixing nut with the supplied spanner.

The fixing nut must always be firmly tightened to the motor shaft.

17 Faults

If the fault cannot be rectified according to the faults table, Hettich customer services must be informed.

Please state the type of centrifuge and the factory serial number. Both values are visible on the centrifuge type plate.


MAINS RESET: - Mains switch OFF for longer than 10 secs.
- Mains switch ON.

Fault	Display	Cause of fault	Remedy
No display	---	No voltage Triggering of the excess current cut-out.	- Check distribution voltage. - Check fuses. - Mains switch ON.
System reset	- 2 -	Failure of power supply during a run Unbraked runout.	- When stationary, open lid and press (START) key.
Balance error	- 3 -	Balance error on the motor axis due to weight differences in the rotor loading.	- Open lid when rotor is stationary. - Eliminate balance error.
Kommunication	- 4 -	Fault in control unit or power unit. Unbraked runout.	- Perform mains reset when stationary.
Overload	- 5 -	Motor or motor control defective.	- Perform mains reset when stationary.
Overvoltage Undervoltage	- 6 - - 8 -	Supply voltage outside tolerance (see Technical Data) Unbraked runout.	- Perform mains reset when stationary. - Check supply voltage.
Excess temperature	- 9 -	Excess temperature switch in motor has triggered. Unbraked runout.	- When rotor is stationary, open lid using emergency unlocking (see Emergency Unlocking chapter). - Allow motor to cool down.
Controller watchdog	- C -	Fault in control unit Unbraked runout.	- Perform mains reset when stationary.
Lid error	- d -	Unbraked rundown, after standstill lid release	- Perform mains reset after standstill.
Short circuit	- E -	Short circuit in control unit / power unit.	- Perform mains reset when stationary.
Configuration	- F -	Incorrect machine version	- Inform Customer Service.

18 Change mains input fuse

Switch off the mains switch and separate the centrifuge from the mains!

The fuse holder (A) with the mains input fuses is located next to the mains switch.

- Remove the connecting cable from the machine plug socket.
- Press the snap-fit (B) against the fuse holder (A) and remove.
- Exchange defective mains input fuse.
- Reinsert the fuse holder until the snap-fit clicks shut.
- Reconnect the centrifuge to the mains supply.

19 Acceptance of the centrifuges for repair

If the centrifuge is returned to the manufacturer for repair, it must be decontaminated and cleaned to protect persons, environment and material.

We reserve the right to accept contaminated centrifuges.

Costs incurred for cleaning and disinfection are to be charged to the customer.

We ask for your understanding in this matter.

You will find further information on our homepage: <http://www.HettichLab.com>

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All rights reserved. No part of this publication may be reproduced without the written prior permission of the copyright owner.