
Labofuge 200 Biofuge 13

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PRECAUTIONS AND HAZARDS

- Please read and observe the operating instructions very carefully!
- The centrifuges must not be operated by unqualified personnel!

ATTENTION-DANGER!

- Operation of the centrifuges in explosion-hazardous locations is prohibited!
- The emergency lid release must not be operated while the motor is running!
- Never transport the centrifuges with a rotor installed on its shaft!

The centrifuges should never be operated, if

- unqualified personnel have tampered with the mechanical and electrical components
- a breakdown exists or defective accessories are inserted
- the rotor is not installed exactly
- the rotor is not loaded symmetrically
- the lid is not closed correctly on both sides

Labofuge 200

- the expiration date of the rotor has passed (see imprint on rotor upper part) or 10000 rotor cycles are reached
- the rotor is not correctly installed after being disassembled (see 4.3 Assembly and 4.4 Installation of the rotor)
- the maximum weight of a tube including the sample exceeds 24 g

Biofuge 13

- the rotor cover is not placed correctly
- the density of the sample material exceeds $1,2 \text{ g} \cdot \text{cm}^{-3}$ for top speed operation. Calculate the admissible top speed dependent upon density by the following formula:

$$N_{\text{adm}} = 13000 \text{ rpm} \sqrt{\frac{1,2}{\text{density of sample}}}$$

This is the maximum allowable setting of the speed.

ATTENTION – WARNING!

Old rotors # 3743 without a through-hole for adaption must not be used in the Biofuge 13 with the new motor shaft (long, one-slotted, cylindrical adaptor with rubber insert)!

Warning!

In case of imbalance (e.g. tube breakage or non-symmetrical loading) it is necessary to press the "stop"-key otherwise the centrifuge could be damaged.

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1. Description of unit

1.1 Definition

Both units are microprocessor controlled laboratory centrifuges with maintenance-free induction drives for many separation tasks.

They are very sturdy, conform to protective class 3 (MedGV in German) and comply with the German safety regulations VDE 0700 and UVV (VGB 7z).

Labofuge 200

This centrifuge is delivered with an integrated angle rotor.

1.2 Range of Application

Labofuge 200

The centrifuge is mostly used in medical offices and small-scale laboratories for routine tests of blood, urine etc. Also very suitable as stand-by unit in large-scale laboratories for separation tasks in research and production.

Biofuge 13

The centrifuge is designed for the use in biotechnological laboratories.

Typical applications

- Pelleting of protein precipitates
- Mini preps for Plasmid DNA isolation
- Enzymatic and other biochemical analyses
- Pelleting of antibodies
- Concentration of algae
- Separation and concentration of cells and fragments
- Fractionation of bacteria and yeast cells

1.3 Mechanical Construction

The frame of the Labofuge 200 is a solid steel construction with an integrated armour chamber and a suspension for the drive and lid mechanics. The casing is made of impact resistant plastic material and is mounted to the frame. This has the effect of reducing noise and vibration. The lid is an all-metal construction with an integrated window shield and links on both sides for the automatic lid interlock.

1.4 Electrical Construction

The microprocessor with the digital displays for speed, time and operating conditions (LED's) and the keys to adjust the parameters and to control the centri-fuge are completely installed together with the power electronic behind the front panel.

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2. Technical Data

2.1 Performance Features

Drive:	maintenance-free, frequency-controlled induction motor of long durability
Labofuge 200	
Max. speed:	5300 rpm
Max. RCF:	3030 × g (dependent on max. radius)
Max. capacity:	12 × 15 ml
Speed display:	digital, 4 digits
Speed selection:	1600 rpm to 5300 rpm in stages of 100 rpm
Biofuge 13	
Max. speed:	13000 rpm
Max. RCF:	13793 × g (dependent on max. radius)
Max. capacity:	18 × 2.0 ml
Speed display:	digital, 5 digits
Speed selection:	1600 / (2000, #3645) rpm to 13000 rpm in stages of 100 rpm
Time display:	digital, 2 digits – seconds or minutes are displayed automatically
Time selection in 2 modes:	
"hd"-mode:	continuous operation (hold)
Min-mode:	1 min to 99 min adjustable in steps of 1 min
Key for lid opening:	only activated by stationary rotor (safety circuit)
Quick-run-key:	centrifugation with max. acceleration and deceleration
Re-starting:	during braking phase
Diagnostic message:	lid position
Data storage:	speed and time adjustment
Data memory:	with NV-RAM quasi unlimited
Dimensions:	(H × W × D) 240 × 284 × 375 mm
Weight:	10 kg
Range of ambient temperature:	
– during operation	4 °C – 35 °C
– in case of storage or shipment	–10 °C – 50 °C
Max. sample-overtemperature:	15 K above ambient temperature

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3. Installation Procedures

3.1 Location

After removing the special packing, the centrifuge must be put onto a sturdy, resonance-free table.

3.2 Electrical Connection

Before connecting the unit to the mains supply, make sure that

1. the mains voltage imprinted on the identification plate is equal with the one available.
2. the mains connection is correctly fused (max. 16 A, type K used for instruments).

After pushing in the power-supply plug, the unit is ready for operation.

3.3 Lid Opening

After the centrifuge is supplied with voltage, an initial control of the electronics is performed (all display elements are illuminated). When the system check stops, the yellow lid-led goes on. The lid can be opened now. By pressing the "lid"-key the lid will be unlocked and can be opened slightly by hand.

3.4 Rotor Installation Labofuge 200

The rotor is already installed and secured for shipping. Before putting the unit into operation remove the transport packing, such as the cardboard pieces inserted between rotor and casing.

Important Notice!

The rotor must turn freely and the locking nut must be tight! During use the locking nut (central winged nut) may work slightly loose and in this case will generate abnormal noise. If this occurs the locking nut must be re-tightened. Prolonged usage of the centrifuge with a loose locking nut may also result in damage to the drive pin, the motor bearings and motor mounts. Please check the locking nut regularly to ensure that it is properly tight.

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Biofuge 13

Before installing the rotor #3757, make sure that the rotor chamber is free of contamination (dust, glass splinters etc.).

The specially formed adaptor hub is suitable for mounting the rotor only in the slot position. The rotor has to be mounted in the appropriate position of the hub (using both hands) and then has to be pressed from the top into the slot (Fig. 2).

Important!

The rotor is locked in place only when it has been firmly pushed and locked into the hub slot (Fig. 3). Then snap on the rotor cover (Fig. 4).

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Biofuge 13

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4. Rotor Exchange

No tools are required to install and remove rotors.

Labofuge 200

The rotor of the Labofuge 200 must be changed not later than 5 years (see expiration date stamped on the upper part of the rotor), or after 10000 working cycles.

4.1 Removal

Labofuge 200

When turning anti-clockwise, the central winged nut (see Fig. 1) can be removed from the motor shaft. By full load, take out 2 opposing tubes, put your index fingers into the holes and pull the rotor vertically upwards from the motor shaft.

Biofuge 13

Grasp the rotor on opposing points and detach from the cone by pulling vertically to the top.

Important!

Never pull up the rotor non-vertically otherwise the drive shaft will be damaged.

4.2 Disassembly (only Labofuge 200)

With a slight thumb pressure to the center piece (with support of the index fingers) the rotor parts can be separated from each other. Removing the upper part, the tubes slip through the borings and remain properly on the lower part.

4.3 Assembly (only Labofuge 200)

When assembling both rotor parts make sure that the pin of the upper part fits exactly into the slot positioned on the lower part. When pressing slightly, both parts snap into each other.

4.4 Installation

Labofuge 200

Before installing the rotor make sure that the drive pin on the motor shaft fits exactly into the notch provided in the rotor for this purpose. In this case the cylindrical part of the motor shaft is exactly leveled with the end of the rotor boring.

Note!

When installing, turn the rotor in appropriate position and tighten the central winged screw by hand. The pin on the upper part indicates the direction of the notch.

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5. Loading

Warning!

To avoid imbalance, jerky operation and resultant damage, make sure that the centrifuge is always loaded symmetrically.

Important! Make sure that the opposing rotor positions are loaded with equally filled and identical centrifuge tubes.

Note: Glass tubes have different wall thickness and therefore hold different volumes. The eye balancing of opposing tubes is not always sufficiently accurate.

The opposing tubes can be tared with a pharmaceutical balance.

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5.1 Usage of Commercial Sample Vessels (only for Labofuge 200)

Dimension of sample vessels	accessories		mounting instruction
	Cat.-No.	designation purpose	
17 x 100 mm 16 x 100 mm			Tube can be inserted directly into the rotor borings ATTENTION – WARNING When the rubber cushions for Monovettes are used, no more than 6 tubes can be inserted!
16 x 100 mm	# 3763	protective container prevents contamination in case of glass breakage	ATTENTION – WARNING! Do not misuse the protective containers as sample vessels! Put the tubes into the protective containers and insert them completely into the rotor borings
Vacutainer 13 x 75 mm with or without screw cap	# 3761	reducing adaptor with distance plug	to hold short 13 mm tubes 1. Disassemble the rotor (see 4.2) and remove the Monovette adaptors if present 2. Mount the reducing adaptors with inserted distance plug (with pins placed onto the conical surface of the rotor lower part) 3. Re-assemble the rotor and install (see 4.3 and 4.4) 4. Insert the short 13 mm tubes through the rotor borings into the premounted adaptors
Vacutainer 13 x 100 mm with or without screw cap	# 3761	reducing adaptor without distance plug	to hold long 13 mm tubes 1. Disassemble the rotor (see 4.2) and remove the Monovette adaptors if present 2. Push out the distance plugs using a long 13 mm tube (see Fig.5) 3. Mount the reducing adaptors by placing the distance pins on the conical surface of the rotor lower part 4. Re-assemble the rotor and install (see 4.3 and 4.4) 5. Insert the long 13 mm tubes through the rotor borings into the premounted adaptors
Sarstedt-Monovettes 10 ml with screw cap	# 3762	Monovette adaptor (rubber cushion)	for correct positioning of Monovettes ATTENTION – WARNING! Operation without Monovette adaptors provides rotor damage! 1. Disassemble the rotor (see 4.2) 2. Press the rubber cushions with the chamfer edges facing down into the bottom rim of the rotor lower part 3. Re-assemble the rotor and install (see 4.3 and 4.4) 4. Load the rotor symmetrically with maximally 6 Monovette tubes

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ATTENTION!

Do not use vessels made of Polystyrole! Use only tubes made of Borosilicate glass being resistant to fracture! The safety against fracture of conical glasses (index No. 1031) is only given generally up to 4000 min⁻¹! On using centrifugal glasses with Monovette adapters maximally 6 glasses are allowed to be arranged central symmetrically in the rotor! Otherwise the adapters must be removed to avoid rising risk of glass fracture.

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6.1 Controlling

6.1.1 Lid Opening

The lid can only be unlocked when the centrifuge is supplied with electricity and the rotor is no longer turning (**yellow** LED must be illuminated) and no error indication is displayed.

After pressing the key, lift the lid slightly by hand and open it.

6.1.2 Centrifugation Start

When start is released (green LED illuminated) then the centrifugation can be started with the "start/stop" key. After starting the lid key is blocked (**yellow** LED is switched off).

6.1.3 Centrifugation Stop

After each start, the stop is released (red LED illuminated). Pressing the "start/stop" key, the centrifugation will be terminated. After each stop, the green LED illuminates again. This means that the interrupted centrifugation can be continued again at any time.

6.1.4 Short Time Centrifugation ("quick run")

Centrifugation can also be started with the "quick run" key if the green LED is illuminated.

Quick run takes place when the key is pressed and held continuously.

ATTENTION!

As long as the key is pressed down, the rotor accelerates to its highest speed. Hereby, the pre-set speed (see 6.2.1.) will be ignored.

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6.2 Programming

6.2.1 Speed Setting

By pressing the "set" key of the speed field, the speed setting mode is activated (setting position flashes). The shown speed value can now be changed in steps of 100 rpm with the +/- counting keys (for quick input keep the key pressed). By pressing the "set" key again, the shown value will be stored as new set speed and speed setting will be terminated.

6.2.2 Time Setting

The time setting can be activated with the "set" key of the time field. The flashing position can be changed in min-steps with the +/- keys. Thus, centrifugation runs can be programmed from 1 min to 99 min. If time is not preselected (hd indication) then continuous operation is possible. The "set" key must be pressed a second time to store the new value.

6.3 Display

6.3.1 Speed

During the acceleration and centrifugation phases the value of the set speed (see 6.2.1) is displayed. When the centrifuge brakes, "br" will be indicated by the microprocessor.

6.3.2 Time

Up to 1 min, the time display takes place rhythmically in seconds. Beyond 1 min only full minutes are indicated. If a time value is preselected, the last minute will be counted down in seconds as well. Using "quick run" the time will be indicated completely in seconds.

6.3.3 Standby Mode

If the centrifuge is not used for some time, the microprocessor will switch the entire digital display off to avoid unnecessary power consumption. This will be indicated by a flashing point in the time display. By pressing one of the keys the display will be activated again, and the key function will be performed.

6.3.4 Cycle Counter Checkup (only important for Labofuge 200)

The number of cycles recorded may be checked at any time by closing the lid and keeping both "set"-keys pressed at the same time. The accumulated rotor cycles will be displayed in the speed field. **Old variant:** All decimal points of the display elements are illuminated when the last 100 cycles are reached. Upon being stepped down to zero the cycle counter will be stepped up now with each cycle and all decimal points will remain illuminated.

Important!

When 10000 cycles are terminated at latest, contact the nearest Heraeus Service Department to overhaul the Labofuge 200.

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7. Maintenance

Both units are equipped with a cycle counter which the microprocessor counts down or upwards by each starting process (see 6.3.4). Upon determination of 10000 working cycles a general check of the Labofuge 200 is required (see 7.2 Service Kit).

Biofuge 13

The carrier hub has to be handled carefully (smear slightly with silicon-grease once per week).

7.1 Cleaning

To clean the plastic surfaces a cloth moistened with non-corrosive liquid should be used.

If a glass or plastic tube breaks during centrifugation, the contents will be spilled.

In such a case, remove the glass splinters and the spilled liquid carefully. For this purpose remove and clean the rotor.

Labofuge 200

Loosen the central nut, separate both rotor bowls by pressing slightly with the thumbs to the center piece and do the cleaning (rotor is suitable to be dish washed).

Important!

Do not use alkaline cleaning reagents as they will affect the durability of the rotor! The rotor can be autoclaved maximally 20 times at 121°C. When assembling and inserting the rotor, pay attention to section 4.3 and 4.4.

Biofuge 13

For the cleaning of the aluminium rotor only neutral cleaners (pH 6.5 to 7.5) must be used.

We recommend rinsing the rotor with distilled water and drying it in a drying chamber at 50°C.

CAUTION!

Autoclaving is not allowed!

If infectious material leaks out, a disinfection of the rotor will be necessary. For this purpose use only an aldehyde-based disinfectant.

7.2 Service Kit (only for Labofuge 200)

The service kit includes all parts which have to be replaced by Heraeus service personnel or adequate technical personnel after termination of 10000 working cycles.

Spare parts (Cat. No. 75 901 001):
Rotor, motor supporting device (3 rubber feet, 6 sandpaper discs, mounting screws and washers), spin chamber sealing, motor rubber flange, variants of cycle counter components (NV-RAMs):
Labo 8 only for processor 436
1052 V01 only for processor 465
1052 V03 only for processor 474.

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8. Procedure in the Event of a Defect

Should the centrifuge operate jerkily, cause abnormal noise or vibrations, this may point to an unbalanced load of the rotor. Stop the centrifuge immediately with the "stop" key and eliminate the imbalance.

If the defect cannot be determined or if other error messages (e.g. E-03) are shown on the speed display, contact your nearest Heraeus Service Dept.

CAUTION!

Repairs must be performed only by qualified personnel!
Only original spare parts must be used!

8.1 Emergency Lid Opening

To open the electrically locked cabinet lid in case of power failure or defects, put a cylindrical pin into the borings located on both sides in the panel (see Fig.8) and unlock the lid.

ATTENTION – DANGER!

Never open the lid mechanically as long as the rotor is still running. You can check this through the round window in the middle of the lid.

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Fig.8 Emergency Lid Opening – Use Only by Stationary Rotor!

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8.2 Instrument Fuses


The laboratory centrifuges are protected against short circuits within the unit by 6.3 AT (slow blow). The fuses may only be replaced by our customer service representatives.

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Performance Features

Max. speed: 5300 rpm
 Max. RCF: 3030 × g (dependent on max. radius)
 Max. capacity: 12 × 15 ml
 Speed display: digital, 4 digits
 Speed selection: 1600 rpm to 5300 rpm in stages of
 100 rpm

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
Performance Features

Max. speed: 13000 rpm
 Max. RCF: 13793 × g (dependent on max. radius)
 Max. capacity: 18 × 2.0 ml
 Speed display: digital, 5 digits
 Speed selection: 1600 / (2000, #3645) rpm to 13000 rpm
 in stages of 100 rpm

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Accessories of Rotor #3757

#3750 Adapters 0.2ml PCR
 #3758 Adapters 0.5/0.6ml
 #3759 Adapters 0.25/0.4ml

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