**LABOFUGE 400e**  
**LABOFUGE 400**  
**LABOFUGE 400R**

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**PRECAUTIONS AND HAZARDS**

Before attempting to operate the centrifuge **read** all parts of this manual **carefully** to insure smooth operation and avoid damage to the unit or its accessories.

If a malfunction occurs, consult the Troubleshooting Guide at the end of this manual. If problems persist, call your nearest Heraeus Service Department or its accredited Agency.

- The centrifuge **must not** be operated by **unqualified personnel!**

- Avoid causing damage to the unit or its accessories through **incorrect operation!**

- **Never transport** the centrifuge with a **rotor installed on its shaft!**

- Use **only original adapters, tubes and spare parts!**

**ATTENTION - DANGER!**

For personal and environmental safety pay attention to the following features:

- Do not operate the centrifuge in **explosion-hazardous locations!** The LABOFUGES are neither explosion-proof nor inert gas shielded.

- **Never stay** in the safety zone of 30cm around the centrifuge nor **deposit dangerous goods** inside this zone during centrifugation!

- Do not spin samples which can create **flammable or explosive mixtures** when exposed to air, unless adequate safety precautions have been taken (e.g. shelter)!

- **Never spin toxic or pathological material** without or with **defective hermetic sealing**, unless adequate safety precautions have been taken (e.g. shelter)! Adequate **disinfecting procedures** must be performed, if dangerous goods have contaminated the centrifuge or its accessories.

- **Never open** the lid manually while the rotor is **spinning!**

- It is forbidden to program a rotor or bucket type, which is not **physically installed** (see 7.2.1)!

The following rules must be strictly observed:

- **Never operate** the centrifuge if it is **not correctly installed** or its **drive shaft is damaged**.

- Do not operate the centrifuge with any **parts or covers removed**.

- Never run the centrifuge when the electrical or mechanical equipment have been tampered with by **unauthorized or unskilled personnel**.

- Never attempt to run the centrifuge with an **incorrectly installed or unauthorized rotor**.

- Do not operate the centrifuge with **incorrectly installed or partial removed buckets**.

- The **maximum rotor load and speed** must be observed!

- Do not spin **corrosive samples** which may impair the material strength of the rotors and buckets **without** taking all necessary precautions (e.g. sealed tubes)

- **Never use accessories** (rotors or buckets) which show **noticeable traces of corrosion** or mechanical damage

- Do not use accessories which are **not approved** by **HERAEUS INSTRUMENTS**, except commercially available centrifuge tubes of glass or plastic.

The manufacturer is only responsible for the **security, reliability and performance** of the unit, if

- the unit is operated in accordance with the operating instructions

- installations, expansions, new adjustments, changes or repairs are performed only by **HERAEUS authorized persons**.

**NOTE:**

When performing centrifugation and other functions which may expose workers to splashed blood or body fluids, all laboratory personnel must follow universal laboratory precautions. For details contact your local laboratory safety officer

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## Setting up the Centrifuge

1. Transport the unit on both sides and set it on a flat sturdy, resonance-free surface in a well ventilated area (without exposure to heat emitting objects) and mark a safety zone of 30 cm around it (Details in Section 4)
2. Insure the mains supply available is equal (in voltage and frequency) with the required on identification plate (see Electrical Requirements, Section 4.2)
3. Plug the mains cable into the adequate wall socket and, if optional available, turn on the main power switch located under the bottom plate right side.
4. Wait for the illumination of the yellow LED next to the "lid" key. Open the lid by pressing the "lid" keypad and remove all packing materials and accessories (rotor locking nut and socket wrench for rotor fixing)
5. Remove the rotor packaging (4.3), mount it on the drive shaft, aligning the groove on the rotor's underside with the cross pin in the drive shaft. Screw locking nut in place (see Installing Rotors 5.1)
6. Rotor must move freely! When installing the swing-out rotor all trunnions must be greased and all 4 buckets of the same type and weight must be occupied (6.1)!
7. Load all adaptors (for maximum load see Section 13.) and tubes making sure the rotor is symmetrically balanced (see Section 6.2)
8. Close the lid firmly on both sides, but do not slam it. The "Open" message will disappear, if the lid is closed correctly

## Programming the First Run

9. To set the Cat. No. of the installed rotor or buckets press the "set"-key of the "time/rotor" field once and then press again and keep it pressed for the whole setting procedure. The rotor sign "ro" will appear and the flashing Cat. No. shown in the "speed/rcf" field can now be changed by using the yellow "+" or "-" keys. If the flashing Cat. No. is equal with the installed one (verify with No. imprinted on inserted rotor or bucket type), rotor setting is terminated on releasing the "set"-key (7.2.1)  
**WARNING!** Always program the Cat. No. of the really installed rotor or buckets, otherwise the maximum allowed bucket speed can be overwritten and damage to persons and property may result
10. To set speed for a timed run, press the "set" key of the "speed/rcf" field. Change the flashing digit by using the "+" or "-" keys in steps of 100 rpm and for termination press the "set" key a second time (for setting RCF instead of SPEED values see 7.2.3)
11. To set time, press the "set" keypad of "time/rotor" field. Change the flashing time value (in minutes) analogously to speed setting with "+" / "-" keys and terminate the procedure by pressing "set" key a second time (7.2.5)
12. LABOFUGE 400R only To set temperature, press the "set" keypad of the temp field and change the temperature value (in °C) analogously to the speed and time setting (7.2.6)
13. To start the centrifugation, press the dual function "start/stop" key. Slow and fast acceleration can be toggled with the "brake" key at any time (7.1.5). After termination of the programmed time, unit will shut off automatically. For premature termination of the centrifugal run press the "start/stop" key again

## Preparing the next run

14. Wait for release of "lid" key (yellow LED indicates end of cycle), open the lid and check the condition of all bags or tubes before removing them from rotor
15. Check for broken glass or leakage into rotor, buckets and/or chamber (for ceaning see 10.1, for disinfecting 10.2)

For more detailed information, operation and guidelines, consult the other marked sections of this manual.

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## TABLE OF CONTENTS

READ ME FIRST .....	2
FAST START UP AND REFERENCE GUIDE .....	3
1. DESCRIPTION .....	6
1.1 Classification .....	6
1.2 Definition .....	6
1.3 Conformity of Standards and Regulations .....	6
1.4 Warranty .....	7
1.5 Standard Accessories .....	7
2. TECHNICAL DATA .....	8
2.1 Common Features .....	8
2.2 Special Data of LABOFUGE 400e / 400 .....	8
2.3 Special Data of LABOFUGE 400R .....	9
2.4 "MEGACONTROL" - Features .....	9
3. SAFETY SYSTEMS .....	10
3.1 Rotor Chamber .....	10
3.2 Lid Lock .....	10
3.3 Thermal Motor and Electronics Protection .....	10
3.4 Imbalance Switch .....	10
3.5 Check of Acceleration Rate .....	10
3.6 Emergency Lid Release .....	11
3.7 Overtemperature Protection .....	11
4. INSTALLATION .....	12
4.1 Location .....	12
4.2 Electrical Requirements .....	12
4.3 Energizing .....	12
4.4 Opening the lid .....	12
5. ROTOR INSTALLATION AND REMOVAL .....	13
5.1 Rotor Mounting .....	13
5.2 Greasing of Rotor Trunnions .....	14
5.3 Rotor Cover .....	14
5.4 Rotor Removal .....	14
6. ROTOR LOADING .....	15
6.1 Bucket Positions .....	15
6.2 Symmetrical Charge .....	15
6.3 Tube Filling .....	16
6.4 Permissible Imbalance .....	16
6.5 Partial Load .....	17
6.6 Loading Limits .....	17
7. "MEGACONTROL" PROGRAMMING .....	18
7.1 Control Keys .....	19
7.1.1 Lid .....	19
7.1.2 Start .....	19
7.1.3 Stop .....	19
7.1.4 Quick Run .....	19
7.1.5 Smooth Acceleration and Deceleration ("brake") .....	19
7.2 Programming .....	20
7.2.1 Rotor or Bucket Selection ("rotor") .....	20
7.2.2 Speed Selection .....	20
7.2.3 RCF Selection .....	20
7.2.4 Brake Cut-out Speed ("set" + "brake") .....	21
7.2.5 Time Selection .....	21
7.2.6 Temperature Selection ("temp" only LABOFUGE 400R) .....	21

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## TABLE OF CONTENTS

7.3	Displaying Functions .....	22
7.3.1	Speed .....	
7.3.2	Relatively Centrifugal Force ("rcf") .....	
7.3.3	Time .....	
7.3.4	Rotor or Bucket Cat. No. ("rotor") .....	
7.3.5	Sample Temperature ("temp" only LABOFUGE 400R) .....	
7.3.6	Standby-Operation .....	
8.	ACCELERATION AND BRAKING PROFILES .....	23
8.1	Profile Principles .....	
8.2	Selection of Profiles .....	
8.3	Turbulence during Deceleration (Swing-out rotors) .....	
8.4	Brake Settings for Delicate Samples .....	
9.	TEMPERATURE FEATURES .....	24
9.1	LABOFUGE 400e / 400 with Direct or Indirect Air Cooling .....	
9.1.1	Preparing a Speed - Temperature Chart .....	
9.1.2	Thermo-delicate Samples .....	
9.2	Temperature Control LABOFUGE 400R .....	25
9.2.1	Pre-cooling or Pre-heating of Rotors .....	
9.2.2	Pre-cooling or Pre-heating with Samples .....	
10.	MAINTENANCE AND CARE .....	26
10.1	Cleaning .....	
10.2	Disinfection .....	
10.3	Greasing .....	27
10.4	Glass Breakage .....	
10.5	Condensed Water .....	
10.6	Servicing Schedule .....	28
11.	TROUBLE SHOOTING .....	29
11.1	Measures to Do-it-yourself .....	
11.2	Servicing .....	32
12.	IDENTIFICATION OF UNITS .....	
13.	ROTOR / BUCKET TABLE .....	

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## 1. DESCRIPTION

### 1.1 Classification

The **LABOFUGES 400e / 400 / 400R** are designed to perform a variety of separation tasks in small-scale laboratories. All units are constructed as bench-top models with an extensive range of accessories, microprocessor-controlled and equipped with a maintenance-free induction drive.

### 1.2 Definition

The **LABOFUGES 400 / 400R** have a low vibration solid polymer bottom plate (be wasted on **LABOFUGE 400e**) The housing is build of a two-pieced sheet steel and surrounds the bottom plate completely. The front panel is made of a resistant plastic material on which an easy cleanable foil is glued. Housing and front panel are fixed tightly together by screws.

The rotors of all three units are started directly with a frequency controlled and maintenance-free induction drive (brushless) with high acceleration power. This drive ensures noiseless, low vibration running at high speed with high reliability.

The drive is mounted with shock absorbers to the bottom plate and sealed with a rubber flange to the rotor chamber.

The central processor manages the controlling, the functions for speed measuring and regulation, the temperature measuring and regulation (only **LABOFUGE 400R**), safety control and error coding.

An additional microprocessor is responsible for the key and display functions.

The centrifuge is controlled by the user friendly "MEGACONTROL" system (see section 7.)

#### LABOFUGE 400e

The air flow, drawn in through the motor cover's vent holes, cools the rotor and is discharged from the centrifuge via the outlet of the lid sealing.

#### LABOFUGE 400

The cooling air flow is performed by 2 powerful fans mounted on the rear casing.

#### LABOFUGE 400R

This units is equipped with a compressor plant filled with CFC-free refrigerant for temperature control (see section 9.2)

### 1.3 Conformity of Standards and Regulations

These products are manufactured and tested in accordance with the following current standards:

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## DESCRIPTION

### Regulations

- Accident prevention regulations for electrical plants and supplies UVV VBG 4.
- Accident prevention regulations for cooling units, UVV VBG 20 (only LABOFUGE 400R).

Safety inspections according to GS or CSA, NRTL/C.

### 1.4 Warranty

The LABOFUGES carry a one-year manufacturer's warranty against defects in materials and workmanship. The CENTRIFUGES will be repaired or replaced at no charge for parts and labor for period of one year after shipping to the customer.


The warranty card is enclosed with the printed documentation. Please complete and return this card after the delivery of your CENTRIFUGE.

### 1.5 Standard Accessories

The centrifuges are supplied with a special rotor locking nut (Spare Part No. #20054407) and a 6mm socket wrench (#20360047) to fasten the rotors (see Fig. 1). 50 ml of anti-corrosive oil (Spare Part No. #70009824) to protect the metallic surface of rotors and buckets and a small vial of special grease (Spare Part No. #70006692) to coat the trunnions of swing out rotors (see 10.3) are supplied as well.

Printed documentation supplied are:

- a prospectus/catalogue with all available accessories
- the operating instructions
- a warranty registration form
- a plastic case for keeping all documents (should be fastened aside the centrifuge)

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2. TECHNICAL DATA	
<b>2.1 Common Features</b>	
Base construction:	low vibration solid polymer bottom plate (except LABOFUGE 400e)
Housing:	made of two parts of sheet steel, fixed tightly by screw and covered with annealing color
Operational panel:	synthetic resin for incorporation of key and display elements covered with an easy cleanable membrane
Operational system:	by programming of the user-friendly MEGACONTROL (features see 2.2)
Drive:	frequency controlled, maintenance-free induction motor
Rotor chamber. Dimension (Ø x H):	made of stainless steel (377 x 187)mm (3 1/2" x 7")
Lid Opening:	key button, the lid pops up and opens by the supporting of a gas lid stay
Lid Closing:	self-acting by pressing the lid down and pushing on top of lid lock
Lid Locking at run:	double safety circuit (hard- and software controlled)
Lid unlocking in case of emergency:	by using a rip cord through casing hole in front of unit
Starting / Stopping:	dual function key button indicated by LED's (start at green, stop at red)
Re-starting:	at any time during braking phase
Quick-run/ -stop:	with "quick run" press-and-hold key / quick stop acts on key release
Imbalance stop:	automatic during acceleration, if allowed imbalance is exceeded
admissible range of ambient temperature:	5 °C - 40 °C during operation, no condensation -10 °C - 50 °C storage or shipping
relative humidity:	max. 80% upto 31°C, over there linear decreasing upto 50% at 40°C

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### 2.4 "MEGACONTROL"-Features

Digital parameter displays:	<ul style="list-style-type: none"> <li>• speed or RCF value</li> <li>• time or rotor</li> <li>• temperature (only LABOFUGE 400R)</li> </ul>
Rotor selection:	integral, arranged in ascending rotor or bucket Cat. No's (see Rotor/Bucket Table in section 13)
Speed selection:	in the range from 300 rpm upto max rotor speed (see rotor/bucket table) adjustable in stages of 100 rpm
RCF selection:	in the range from min. g-value (related to 300 rpm) upto max. RCF (see rotor/bucket table) - adjustable in stages of recalculation of 100 rpm
Time selection - "hd" mode - min mode	<ul style="list-style-type: none"> <li>• continuous operation (hold)</li> <li>• 1 minute to 99 minutes, adjustable in minutes</li> </ul>
Brake cut-out point:	adjustable in stages of 100 rpm up to the preselected speed
Temperature selection (400R):	0 °C - 40 °C adjustable in stages of 1 °C
Acceleration:	2 profiles smooth start (slow running and fast running, can be switched over at any time with the "brake" key
Deceleration:	1 profile and unbraked deceleration, can be switched over at any time with the "brake" key
Time display by "quick run"	seconds mode: 1 s - 99 s
Parameter storage: (permanet data storage with NV-RAM)	<ul style="list-style-type: none"> <li>• speed or rcf value</li> <li>• rotor or bucket Cat No.</li> <li>• running time</li> <li>• brake cut-out point</li> <li>• temperature (only LABOFUGE 400R)</li> </ul>
Diagnostic indications:	<ul style="list-style-type: none"> <li>• lid open (not correctly closed)</li> <li>• imbalanced loading</li> <li>• false rotor</li> <li>• general trouble in the unit</li> <li>• data transmitting errors</li> </ul>

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#### 3.1 Rotor Chamber

##### LABOFUGE 400e

For reason of air suction cooling the fine steel chamber is not sealed to the environment.

##### ATTENTION - WARNING!

For safety use only rotors and buckets with hermetic sealed accessories!  
In case of danger of glass breakage keep sufficient distance to the unit's rear and right side!

##### LABOFUGE 400

The fine steel chamber is sealed to the unit's interior by a motor rubber cover and to the lid (in closed position) by a central profiled rubber ring (see Fig. 1).

##### LABOFUGE 400R

The sealed rotor chamber is made of stainless steel and is surrounded by the evaporator coil and protected by a CFC-free thermo-insulation.

#### 3.2 Lid Lock

The cabinet lid can only be opened when the mains power is on and the rotor speed is at zero.

The centrifuge can only be started when the lid is correctly locked and the "open" message has turned off.

#### 3.3 Thermal Motor and Electronics Protection

The power electronics is switched off in the event of motor or additionally on LABOFUGE 400R in the event of heat sink overtemperature to protect the drive and electronics (see 11.1 - error message "lid" during operation).

#### 3.4 Imbalance Switch

The centrifuges are all equipped with an adjustable imbalance switch.

If the allowable imbalance is exceeded, the message "bal" will alternate with the speed display and the drive brakes to standstill automatically (see 11.1).

The imbalance switch responds during the acceleration of the rotor if an unacceptable imbalance exists.

##### ATTENTION - WARNING!

When using the light rotors #3765 and #8171 (only LABOFUGE 400 / 400R), the imbalance function is out of order! Do not allow more than the maximum imbalance given in Rotor/Bucket Table in section 13.).

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### 3.5 Check of Acceleration Rate

Whenever the rotor #3765 and #8171 (only LABOFUGE 400 / 400R) is set, the microprocessor will check the acceleration rate to protect the drive system or other types from damage in the event of false speed or RCF programming (see 11.1 E-18).

**NOTE!** With this rotor setting slow acceleration will not be accepted (E-18).

### 3.6 Emergency Lid Release

#### ATTENTION - DANGER!

It is very dangerous to open the lid while the rotor is spinning!

In the event of a failure of the electrical supply, the automatic lid release is rendered inoperative.

To remove the samples when the rotor has come to standstill, a manually operated release is provided.

Performance as follows:

- look through the lid window to ensure rotor is not spinning
- **turn off the power, remove mains plug**
- remove plastic stopper in the center underneath the front casing using a screwdriver, spatula or knife
- pull the stopper and its attached rip cord downwards (see Fig. 2) to unlock the lid mechanically

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### 3.7 Overtemperature Protection

The LABOFUGE 400R is equipped with overtemperature protection. If the temperature in the display exceeds 50 °C, the microprocessor will indicate an overtemperature alarm (E-07, see 11.1).

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
**4.1 Transportation and Location**

Open the cardboard box carefully, remove the transport protection parts (foam). Take notice to centrifuge weight (see 2. Technical Data), seize the instrument on both sides and transport the unit to its location. Avoid impacts and do not drop the unit to preserve it from being damaged!

Location requirements:

- 1.) stable and non-vibrating floor (e.g. concrete platform)
- 2.) sturdy, even laboratory table or a suitable carriage with lockable wheels
- 3.) keep away from nearby heat emitting element or instrument
- 4.) no direct sunlight
- 5.) clearance of 10 cm to the backwall and 15 cm to both sides to ensure unrestricted air circulation
- 6.) balance the instrument with aid of a spirit level
- 7.) keep bottom range of LABOFUGE 400R proper and free (motor fan is installed in this area)

**4.2 Electrical Requirements**

Before connecting the Labofuge to the mains supply, make certain that:

- 1.) the line voltage and frequency imprinted on the identification plate is equal with the one available,
- 2.) the line voltage circuit breaker is maximum a 16Amp, typ K, slow release characteristic for commonly used instruments,
- 3.) an emergency switch is installed outside the centrifuge room to disconnect the power supply in case of a troubled run.

**4.3 Energizing, Safety Precautions**

Connect the cord plug to the appropriate wall socket and if optional available turn on the mains switch located under the bottom casing to energize the centrifuge. All displays and diodes will light up for two seconds (system check).

**Note:** When powered off, all parameters of the last run are maintained.

During centrifugation the following safety regulations have to be observed (see also 1.3):

- 1.) Mark a safety operational cycle of 30 cm around the centrifuge.
- 2.) Persons and dangerous materials (e.g. infectious liquids) must not stay in this area.
- 3.) No objects at all must not be placed onto the unit

**4.4 Opening the lid**

The lid can only be opened with the "lid" key when the yellow LED lights up. This is only possible when the rotor is not turning and no error messages are displayed.

**ATTENTION - DANGER!**

The mechanical lid unlocking device (see 3.6) must not be used to interrupt centrifugation.

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
## 5. ROTOR INSTALLATION AND REMOVAL

### CAUTION!

Before inserting the rotor, make sure that the rotor chamber is free of contaminations (dust, glass splinters etc.)!

Condensate or residual sample liquid must be removed before centrifugation!

The temperature of rotor, drive shaft and locking nut must be in a range of 10°C - 30°C during fastening to avoid self loosening caused by thermal balance!

The rotor must be removed before transporting the centrifuge!

### 5.1 Rotor Mounting

A slot matching the cross pin on the drive shaft is grooved into the underside of all rotors and must be aligned precisely when installing the rotor in the chamber.

### WARNING!

Position the rotor so that the slot grooved into it's underside seats in exact alignment with the cross pin on the drive shaft!

Failure to do so will prevent proper seating and damage to centrifuge and accessories is possible!

### Installing procedure:

- 1.) Turn the drive shaft to place the cross pin in a known position, e.g. long bar pointing to front and rear.
- 2.) Position the grooved notch on the underside of the rotor in exact alignment with the cross pin on the drive shaft (watch the marked lines on top of rotor, eventually remove rotor cover).
- 3.) Press rotor firmly and slowly to place and check the correct seat (see Fig. 3).
- 4.) When the rotor is properly seated, secure the locking nut, by turning clockwise with the aid of the 6mm socket wrench which is delivered with each centrifuge.

### Note!

If the locking nut does not easily engage the threaded portion of the drive shaft, do not try to force the rotor further onto the drive shaft. Instead, remove the locking nut and the rotor, and re-install, making sure that the rotor is properly seated. (To insure correct seating, be sure step # 2 is properly performed.)

Fig. 3

left: Correct seated!  
threaded portion extends  
completely above rotor

right: Wrong seated!  
threaded portion below  
rotor centerpiece

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**5.2 Greasing of Rotor Trunnions**

**ATTENTION!**

The trunnions of the swing-out rotor must always be provided with a thin coat of grease!

The coat of grease ensures smooth movement of the buckets and contributes to vibration-free operation of the centrifuge, particularly during acceleration and deceleration of the rotor.

Greasing should be done regularly, that is after approx. 20 runs or at least once a week.

Before grease is applied to them, the trunnions and the bucket grooves must be thoroughly cleaned with a dry cloth.

**In no case may organic solvents, alkaline or household cleaners be used which contain scouring powder!**

The best possible grease for this use is resin and acid-free vaseline (e.g. ball bearing grease). A suitable grease is available under Cat. No. 70006692.

**Molykote and graphite lubricants are also not permissible!**

**5.3 Rotor Cover**

**ATTENTION - WARNING!**

All rotors must be firmly closed for each centrifugal run with its delivered covers!

The rotors #3765 and #8171 (only LABOFUGE 400 / 400R) are designed with removable covers to keep the tubes in safe positions (#8171), to reduce air friction and noise and allows the drive to turn at maximum given speed values (refer to §13, all values are only for placed covers).

**5.4 Rotor Removal**

To remove the rotor from the chamber proceed in reversed sequence:

- 1.) Open the lid and follow all necessary laboratory safety instructions
- 2.) Remove all tubes and check for broken glass.
- 3.) Remove rotor locking nut by turning counter-clockwise (see imprinted arrow direction), using the added 6mm socket wrench.
- 4.) Size the rotor on opposing sides and pull it straight up.

**WARNING!**

**Never pull up the rotor at an angle. Doing so could bend the drive shaft and cause permanent damage to the motor.**

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## 6. ROTOR LOADING

### ATTENTION - WARNING!

The following instructions must be observed exactly to avoid damage to centrifuge and accessories!

All adaptors must only be loaded with its predestinated tubes or vessels in form, size and material!

Dangerous samples must not be centrifugated without taking all necessary precautions (e.g. shelter, sealed rotor or buckets)!

### 6.1 Bucket Positions

#### WARNING!

In swing-out rotors all positions must be occupied with identical buckets (same Cat. No.).

The buckets may have an identification letter which states the weight class (e.g. "E" or "D").

Only equal weight classes may be inserted in opposing places!

Make sure that the tubes or adapters in opposing buckets are identical (identity of tube distribution and size).

### 6.2 Symmetrical Charge

The rotor must always be loaded symmetrically. Tubes and/or racks should be evenly distributed.

Following figures show how the tubes can be correctly or incorrectly placed in rotor.

Fig. 4  
Rotor fully equipped

Fig. 5  
Rotor partly loaded,  
but symmetrically equipped  
a = best solution

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
## ROTOR LOADING


The tubes must always be distributed within the buckets in such a way that the rotor trunnions are stressed (loaded) evenly.

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**Fig. 6**  
not allowed,  
because trunnions are  
not evenly loaded

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**Fig. 7**  
not allowed,  
because trunnions are  
also not evenly loaded

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### 6.3 Tube Filling

Centrifuge tubes up to approx. 30 mm in diameter may be filled evenly by eye (up to 5 - 10 mm below the rim, depending upon size).

When a rotor is loaded with a large number of equal tubes, small differences in weight between the individual tubes will be statistically compensated.

**ATTENTION!**

Using an angle rotor and tubes **without caps**, only partial filling is possible. The real amount of filling depends on rotor's angle position (that's about 60% to 75% of the maximum capacity).

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### 6.4 Permissible Imbalance

In order to avoid rotor wobble, opposing places must be provided with identically loaded adapters, multiple carriers, tubes, bottles etc. to obtain the best possible balance.

Imbalance causes running noises and produces a negative effect on the drive system (early wear out of the motor support).

The difference in total weight (or load) between opposing buckets is rotor dependant and must not exceed the permissible values (13. Rotor Tables).

The better the taring is done, the smoother will be the centrifugal run and the longer will be the service life of the drive system. In addition you will attain a superior level of separation, because already separated components won't be whirled up again by vibration.

To determine the bucket weights, any laboratory balance is suitable, but a taring balance is best.

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
## ROTOR LOADING

### 6.5 Partial Load

With only a partial load, the tubes must be symmetrically distributed so that an even overall loading of the rotor is achieved (see also Fig. 9).

If a small number of sample tubes is to be centrifuged (less than the full capacity of the rotor), it is advisable to observe the following rules for distribution of the tubes within the adapters.

Fig. 8  
correct loading

### WARNING!

To avoid glass breakage when using tubes greater than 100 mm in length, use the configurations shown in Fig. 9. In this case, always conduct a manual swing-out test with empty tubes before committing samples.

Fig. 9 correct loading for tube length > 100 mm

### NOTE!

Best solution for long tubes upto 132 mm length is the usage of buckets # 8172 in the LABOFUGE 400/400R (vacutainer see Rotor Table in Section 2.4).

### 6.6 Loading limits

The material strength of each rotor is limited with regard to the mass and speed-dependent forces acting on its yoke and other components during centrifuge operation.

HERAEUS rotors are designed in such a way that they still possess high strength reserves at the maximum permissible load and speed. Still the maximum permissible limits must not be exceeded, otherwise the working life will be reduced extremely.

It is the operator's responsibility to avoid overloading the rotor by observing the max. weight of the buckets.

The max. permitted load is reached when all tubes are filled with a fluid of  $1.2 \text{ g}\cdot\text{cm}^{-3}$  density.

When the density is higher than 1.2 the tubes may only be partly filled. The total weight of a bucket may not exceed the value corresponding to total filling with a liquid of  $1.2 \text{ g}\cdot\text{cm}^{-3}$  density.

As an alternative in such cases, it is possible to reduce the maximum permitted speed  $n_{\text{max}}$  in accordance with the following formula:

$$n_{\text{perm}} = n_{\text{max}} \cdot \sqrt{\frac{\text{Weight of Buckets incl. Sample with Density of } 1.2 \text{ g}\cdot\text{cm}^{-3}}{\text{Weight of Buckets incl. Sample with HIGHER Density}}}$$

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## 7. "MEGACONTROL" PROGRAMMING

"MEGACONTROL" is the trade mark of the user friendly control system of HERAEUS Instruments laboratory centrifuges for easy programming, diagnostics, controlling and monitoring of the operating functions.

All key and display elements are integrated in a pad neatly divided into fields for optimum operator guidance.

Fig. 10 MEGACONTROL - key and display board

### Fields for Function Displays

- A Display for speed, RCF value, rotor type, brake cut-off point and special indications for diagnostics
- B Display for running time, "hd" (continuous operation), "ro" (rotor) or "rc" (RCF)
- C Display for temperature (LABOFUGE 400R)

### Programming Keys

- 1 "set" key for speed or RCF selection (see 7.2.2 or 7.2.3)
- 2 "brake" key a) to set brake cut-off point (see 7.2.4)  
b) to switch between fast and slow acceleration  
c) to switch brake ON and OFF (see 7.1.5)
- 3 "set" key for running time or rotor selection (see 7.2.5 or 7.2.1)
- 4 "set" key for temperature selection (400R)
- 5 "+ " key to increase the selected value
- 6 "- " key to decrease the selected value

### Control Keys

- 7 "lid" key to open the lid (7.1.1)
- 8 "quick run" press-and-hold key for short cycles braking occurs upon release (7.1.4)
- 9 "start stop" dual function key for starting and stopping (7.1.2 and 7.1.3)

### Indication Lights

- 10 red LED a) ON indicates fast acceleration  
b) ON indicates brake activated
- 11 yellow LED indicates lid may be opened
- 12 green LED indicates "quick run" may be used
- 13 green LED indicates run may be started
- 14 red LED indicates run may be terminated

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
## 7.2 Control Keys

## 7.2.1 Lid

The lid can only be unlocked with the "lid" key when the centrifuge is powered on, and the rotor is at standstill (yellow LED must be illuminated) and error messages are not indicated.

After opening the lid, "OPEN" will be indicated in the "speed/rcf" display. The starting function is blocked now.

## 7.2.2 Start

When start is released (green LED illuminated), the centrifugation can be initiated with the "start/stop" key. After centrifugation has started, the lid key is disabled (yellow LED is unlit).

## 7.2.3 Stop

The stop function becomes available (red LED is lit) immediately after each start. Pressing the "start/stop" key now will terminate the run.

After pressing stop, the green "start" LED lights again. This means that the interrupted centrifugation can be continued at any time (it's not necessary to wait for rotor standstill, if the runs was stopped inadvertently).

## 7.1.4 Quick-Run

Each centrifugation cycle can also be started with the "quick run" key when the green LED is illuminated.

The quick-run function is designed for short time runnings. Quick-run requires pressing and holding the key.

## ATTENTION - WARNING!

As long as the key is pressed, the motor accelerates at maximum rate towards the highest speed of the preselected rotor or bucket type. That means, the set speed (see 7.2.2) will be ignored when using the "quick run" function.

Centrifugation will continue until the key is released. On release the rotor will be decelerated at maximum braking. If the "quick run"-key is pressed and held for a second time, acceleration will be re-activated.

## 7.1.5 Smooth Acceleration and Deceleration

At any time during the acceleration or deceleration phases the multifunction key "brake" may be used to switch between quick and smooth centrifugation (and back again).

Fast speed changes are generally indicated by the red control LED in the speed display section. This LED goes out on switching over to smooth centrifugation (slow acceleration or unbraked deceleration).

## ATTENTION!

The light rotors #3765 and #8171 (only available for LABOFUGE 400 /400R) are not allowed to start with slow acceleration (see 11.1 E-18)!

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## 7.2.1 Rotor or Bucket Selection



**ATTENTION!**

Rotor settings are only allowed at rotor standstill!

Press "set" key in the "time/rotor" section once and then press again and keep it pressed.

The rotor input mode is indicated with "ro" and the last stored Cat. Number (#) is shown flashing in the "speed/rcf" display section.

The rotor Cat. Numbers (#) can now be displayed in ascending mode by using the "+" key and in decending mode by using the "-" key. After release of the "set" key, the new rotor Cat. No. (#) will be stored.


**ATTENTION - DANGER!**

Several bucket types are available for use with the swing-out rotor #8179 (see 13. Rotor Tables).

Therefore, the Cat.-No. of the inserted bucket must be selected to avoid accidents caused by false programming (speed selection more than maximum admissible bucket speed).

This applies especially to the Microtiter Rotor # 8177 if operated above its maximum allowed speed (see 13. Rotor Tables).


## 7.2.2 Speed Selection

On pressing the "set" key in the "speed/rcf" display section, the speed setting mode is activated (setting position flashes). The indicated speed value can now be changed in steps of 100 rpm with the "+/-" keys (for fast setting, hold pressed).

On pressing the "set" key again, the set value is stored as the new set speed and speed setting is terminated.

## 7.2.3 RCF Selection

Press the "set" key in the "speed/rcf" display section (see above) once, release, press again and keep it pressed (like rotor setting). The display switches to RCF setting mode and "rc" will be indicated in the "time/rotor" display section.


The RCF value can now be changed in steps of 100 rpm with the "+/-" keys. After release of the "set" key, the last set value will be stored.

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## 7.2.4 Brake Cut-out Speed ("set" + "brake")

Press the "set" key in the "speed/rcf" display section first and then press and hold the "brake" key (beginning is the same as rotor or rcf setting).

The setting of the brake cut-out speed occurs analogous to 7.2.2 (Speed Selection) up to the preprogrammed value of the set speed.


If there is no use of any brake cut-out speed, set this value to zero.

After release of the "brake" key, the set brake cut-out speed will be stored.

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## 7.2.5 Time Selection

The time setting can be activated with the "set" key of the "time/rotor" display section. The flashing position can be changed in minute-steps with the "+/-" keys. Thus, centrifugation runs can be programmed from 1 minute to 99 minutes.


If no time is preselected ("hd" indication) then continuous operation is possible. The "set" key must be pressed again to store the new set value.

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## 7.2.6 Temperature Selection (LABOFUGE 400R)

Press the "set" key in the "temp" display section to activate the setting mode (setting position flashes).

The setpoint value of the sample temperature is changeable in degree Celsius -steps with the "+/-" keys in the range from 0°C to 40°C.

By pressing the "set" key again, the displayed value is stored as the new set temperature and the setting procedure is terminated.

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
## 7.3 Displaying Functions

## 7.3.1 Speed Value

The speed display is given priority confront to RCF value and is updated every 0.5 second in the physical unit of  $\text{min}^{-1}$  (rpm).

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
## 7.3.2 RCF Value

The RCF value may be displayed by pressing and holding the "set" key. On key release speed values are displayed again.

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
## 7.3.3 Time Value

During "quick run" centrifugation, the elapsed time will be indicated in seconds (without flashing point). Normally the time is shown in whole minutes with a point flashing secondly.

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
## 7.3.4 Rotor or Bucket Cat.-No.

If the "set" key in the time/rotor section is kept pressed, the set rotor Cat. Number can be checked any time. New rotor settings are only possible at rest (see 7.2.1).

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
## 7.3.5 Temperature (only LABOFUGE 400R)

The display shows the sample temperature in degree Celsius (only senseful after precooling or preheating see 9.2.2) and is activated every 0.5 second.

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## 7.3.6 Standby Operation

If the centrifuge is idle for a few minutes, the microprocessor will switch off the digital display to save energy.

This state is indicated by leaving only the decimal point in the time display flashing. Pressing any key will immediately reactivate the whole display.

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## 8. ACCELERATION AND BRAKING PROFILES

To satisfy numerous applications of the LABOFUGES, 2 different reproducible acceleration and one braking profiles are available.

Additionally the unbraked deceleration can be selected with the "brake" key.

### 8.1 Profiles Principles

The profiles include a cushioned acceleration and deceleration phase to avoid re-suspension of the sample, followed by a phase with quick acceleration for shortest total run times.

The natural resonance vibrations of the rotating assembly occur within the phases of fast acceleration and deceleration to ensure that they exert a minimal effect on the drive system.

### 8.2 Selection of Profiles

The fastest acceleration and braking time can be reached with the "quick run" key. In this case each rotor is operated with the maximum possible acceleration and braking force.

The "start/stop"-key is not only provided for the fast acceleration and deceleration profile (red control LED has to be illuminated), but also for smooth (slow) centrifugation.

The slow acceleration as well as the unbraked deceleration can be processed by the "brake" key (see 7.1.5 - red LED is not illuminated in this case).

### 8.3 Turbulence during Deceleration

With swing-out buckets, the transition from the horizontal to the vertical position occurs between 600 rpm to 800 rpm.

Based on this physical effect a turbulence occurs inside the centrifugation container and may cause a re-mixing of the sample. This turbulence is suppressed more effectively when the centrifuge stops slowly and steadily which can be achieved by using e.g. a brake cut-off point (see 8.4).

### 8.4 Brake Setting for Delicate Samples

With samples producing a firm sediment during centrifugation the maximum braking force can be used.

For samples with easily disturbed sediments you can choose either to work without the brake or to select a brake cut-out point so that the rotor then coasts without braking (smooth) e.g. from below 500, 400, 300, 200 or 100 rpm to rest.

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## 9. TEMPERATURE FEATURES

### 9.1 LABOFUGE 400e/400 with Direct/Indirect Air Cooling

During centrifugation, heat is generated by air friction between the rapidly spinning rotor and the air inside the chamber. The temperature of the rotor, buckets and samples will increase.

The temperature rise depends on the rotor size, bucket type, running time, ambient temperature, location (air circulation) and especially on the rotating speed of the rotor.

**LABOFUGE 400e air suction cooling** A certain amount of the heat will be eliminated by the air flow through the rotor chamber.

**LABOFUGE 400 ventilation cooling** The thin fine steel rotor chamber is cooled from outside by the forced air flow of 2 powerful fans mounted on the unit's rear casing.

For both centrifuges thermal balance between heat generation and heat elimination will be reached after some time.

Heræus INSTRUMENTS has arranged the cooling airflow so that the warming of the samples does not exceed the standard value of  $(38 + 2)^{\circ}\text{C}$  at ambient temperature  $(23 + 2)^{\circ}\text{C}$  when the LABOFUGE 400e/400 is running at maximum rotor speed.

By carrying out appropriate tests it is possible to make a speed-temperature chart for any rotor. But take notice that the sample temperature will be affected by the ambient temperature.

#### 9.1.1 Preparing a Speed-Temperature Chart

To avoid discrepancies caused by evaporation (e.g. of water) it is best to use glycerin in the buckets to measure the temperature. The buckets should be evenly filled to the middle with 80% glycerine (taring). The required speed is selected and the instrument started. The temperature of the glycerine in each of the buckets should then be measured at regular intervals (10 to 30 minutes) until there is no change in the reading.

**ATTENTION - important note!**

It is essential to dip the thermometer into every bucket in order to avoid false readings.

Continue starting and stopping the LABOFUGE 400e/400 with the same running program until the temperature does not change.

#### 9.1.2 Thermo-Delicate Samples

Samples, which can be destroyed by exceeding a fixed temperature (e.g.  $38^{\circ}\text{C}$ ) can be centrifugated in some cases using the LABOFUGE 400e/400 only at reduced speeds.

**Note!**

A LABOFUGE 400R with temperature control is recommended for working with thermo-sensitive samples.

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## 9.2 Temperature Control LABOFUGE 400R

A sensor located in the rotor chamber measures the temperature changes very quickly. Depending on the selected rotor (bucket), speed and temperature settings, the microprocessor determines the parameters for the temperature control and display of the sample temperature.

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Fig. 11 Schematic Diagram of Temperature Control

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In the equilibrium condition (thermal balance), the sample temperature will be indicated with a tolerance of  $\pm 2$  K.

The minimum reachable sample temperature of each rotor or bucket is dependent especially on the set speed and the ambient temperature (see 13. Rotor / Bucket Table).

**NOTE!** The temperature control is only activated if the lid is closed correctly (no "OPEN" message)! This is the premise to let the centrifuge operate at standstill cooling.

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### 9.2.1 Pre-cooling or Pre-heating of Rotor and Empty Buckets

If the sample cannot tolerate temperature differences to the empty buckets, the rotor chamber (incl. rotor) must be pre-cooled or pre-heated first.

For this purpose following procedure is recommended:

- 1.) set the Cat. Number (#) of the installed rotor or bucket type
- 2.) set the required sample temperature
- 3.) set the half or the maximum rotor or bucket speed (for temperatures > 15°C: set 90% of the maximum rotor or bucket speed - defined as pre-heating)
- 4.) set running time to 15 minutes and press the start key

**NOTE!** Pre-cooling at standstill is principally possible, but it takes a cooling period for more than ten times.

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### 9.2.2 Pre-cooling or Pre-heating of Samples

Depending on the additional heat transfer from the rotor pocket (or bucket inwall) via the racks and tubes to the original sample, a much longer running period is necessary when e.g. pre-cooling with content.

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## 10. MAINTENANCE AND CARE

Heræus INSTRUMENTS GmbH recommends to let centrifuge and accessories maintain annually either by an authorized service technician or with help of the Servicing Schedule (see 10.6) by trained service specialists. For this purpose Heræus (or its Accredited Agent, as appropriate) offers an inspection and service contract. Warranty repairs are made without charge and those which are not covered by warranty and after the warranty period are charged to the customer.

### ATTENTION - WARNING!

For protection of persons, environment and property, the user is obligated to perform adequate cleaning or disinfecting procedures!

The application of specific cleaning or disinfecting procedures which are not declared in this operating instructions must be declared first with Heræus INSTRUMENTS!

Misuse can cause damage to the centrifuge or its accessories!

### 10.1 Cleaning

Care of the centrifuge is restricted especially to the cleanliness of the rotor chamber, the rotor, the buckets and their accessories.

Cleaning is not only necessary for reasons of hygiene, but is also essential to prevent corrosion caused by adhesive impurity.

Accessories of anodized aluminium like buckets, adapters etc. are particularly susceptible to corrosion. Clean them with mild detergents of pH values ranging from 6 to 8. For this purpose we recommend to use e.g. RBS Neutral Concentrat (application to producer specifications).

After cleaning, dry the aluminium accessories with a soft cloth or put them in a drying oven (rotors and buckets must never be heated > 50°C).

To increase durability and reduce their susceptibility to corrosion, it is necessary to regularly apply an anti-corrosive oil (Cat. No.70009824) to all anodized surfaces.

### 10.2 Disinfection

#### ATTENTION - DANGER OF INFECTION!

Infectious material arrives at the centrifuge in case of glass breakage or spilling! A disinfection of chamber, rotor, buckets or racks after the run is required! Wear protective gloves!

The rotor and buckets are **not autoclavable!** The Centri-Lab racks may be autoclaved at 121°C.

A contaminated rotor and rotor chamber must be treated with a neutral, universal disinfectant (e.g. with an aldehyde base, such as rpi's "Lift Away").

Best suited for this purpose are disinfectant sprays, such as e.g. Sagromed HBV, which ensure that all rotor and accessory surfaces are covered evenly.

#### ATTENTION!

The rotor #3765 (only LABOFUGE 400/ 400R) **must not be dismantled for disinfection!**

If infectious fluids should have been entered the interior rotor part, proceed as follows:

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## MAINTENANCE AND CARE

- 1.) **Observing all safety rules** (e.g. wearing of gloves), unscrew the locking nut, grip the rotor on opposing points and detach it from the drive shaft by pulling straight-up.
- 2.) Remove all rotor tubes (and any adaptors) for disinfection or recycling.
- 3.) Dip the rotor into a **disinfectant solution** (e.g. 6g formic aldehyde + 1,8g ethanedial + 1,6g glutar aldehyde per 1 liter water) and let it soak - time according to contaminant, turn the rotor upside down and let the solution run out.

### 10.3 Greasing

#### ATTENTION - GREASE THE TRUNNIONS!

The trunnions of the swing-out rotors must always be lightly greased to avoid jerky movements of the buckets.

Greasing should be done daily, but at least once per week!

Before greasing, the trunnions and the bucket grooves should be thoroughly cleaned with a dry cloth. In no case may organic solvents, alkaline or household scouring powder cleaners be used!

We recommend to use our well tested grease which is supplied with the swing-out rotor and is available under Cat. No. 70006692.

#### ATTENTION - important note!

Molykote and graphite lubricants are not permitted!

#### MOTOR BEARINGS - DO NOT GREASE!

The motor bearings are lubricated for the whole motor's life.

#### CAUTION!

Make sure that no liquids or organic solvents penetrate the annular slot around the motor shaft which might wash the grease off the bearings.

### 10.4 Glass Breakage

Centrifuge glassware shows a correspondingly higher failure rate with rising g-forces.

**ATTENTION!** Glass splinters must be removed immediately and completely from the buckets, adapters and rotor chamber. The glass splinters can damage the surface protection or could get stuck on the trunnions of the swing-out rotors where they can impair smooth movement of the swinging buckets.

When glass splinters get into the rotor chamber and are whirled around by air circulation, a very fine (black) metal dust can develop.

Metal abrasion can contaminate the rotor chamber, rotor, buckets and samples. With a strip of vaseline or trunnion grease applied vertically on the chamber wall (from top to bottom), the glass dust can be removed with a smooth centrifugation cycle.

### 10.5 Condensed Water

Air, humidity or non-hermetically sealed samples will cause condensate inside the spin chamber depending on the chamber temperature.

The condensate must be removed from the rotor chamber regularly e.g. with a dry cloth.

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### 10.6 Servicing Schedule

We recommend to perform the following maintenance routine which should be done at least once a year by an authorized service technician only:

- 1. Electrical Installation and Safety**
  - switch OFF centrifuge and disconnect unit from power
  - check voltage supply and mains fusing
  - check condition of plug and wall socket or permanent installation - (let) replace defective parts
  - check cord condition and fixing - replace or refit it
- 2. Location and Mechanical Installation**
  - check the base (ground, adjustable feet etc.) for resonance-free and stable conditions
  - check for a well ventilated place and sufficient distances to walls or adjacent equipment
  - check the levelling of the centrifuge - use any spirit level
- 3. Lid Locking Mechanism and Safety Device**
  - connect the centrifuge to power and switch ON
  - check for easy lid closing and self-acting lid opening - if in disorder, readjust lid looped hook or hinges or move lid lock assembly and smear hook slightly with grease
  - check the emergency lid unlocking device (see 3.6).
  - open the lid and turn the rotor by hand, then close the lid and try to open it using the lid key: the lid must not be opened as long as the speed values ( $> 10 \text{ rpm}$ ) are indicated - if safety circuit is out of function, call the Service of Heraeus Instruments.
- 4. Cleanliness of Spin Chamber and Motor Cover**
  - open the lid and remove the rotor (for loosening turn the locking nut with 6mm socket wrench in arrow direction see rim of rotor chamber)
  - clean the spin chamber with a dry and absorbent cloth (remove all dust and moisture)
  - check the correct condition of the motor cover and take care of the cleanliness of the annular slot around the motor shaft: Penetrating fluids can damage the upper motor bearing!
- 5. Rotor and Bucket Condition and Sealing**
  - check the condition of rotor and buckets (especially all supporting or stressed partitions as jib arms, rim of the bucket's bearing surface etc.): the rotor and/or buckets must not be used any longer, if there are visible traces of mechanical damage or rust
  - check the condition of bucket sealings (hermetic cap) and replace them in case of malfunction
- 6. Rotor Fixing and Motor Shaft**
  - check the condition of the drive (motor) shaft: the centrifuge must not be used any longer, if the drive shaft is damaged (bend, thread is worn out, horizontal grooves etc.)

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## 11. TROUBLE SHOOTING

### 11.1 Measures to do it yourself

Trouble		Cause	Remedy
Indication	Symptom		
all displays remain dark	drive shuts off suddenly, rotor coasts to rest, lid cannot be opened at rest	power failure	check mains supply, if in order call Heraeus
displays are momentarily dark	drive shuts off suddenly and coasts to rest	temporary mains failure	correct seat of mains plug or switch (option) position and start again
unusal loud running noises perhaps with imbalance indication	drive is disconnected brakes or continues running with noise	1. sudden imbalance following tube breakage and mass shifting 2. defective drive or rotor 3. disorder of speed control 4 faulty fan or compressor (Labo 400R)	if the centrifuge does not stop on its own, press "stop" key or switch OFF the mains check for obvious faults in the rotor chamber at standstill and if possible do repair, otherwise call Heraeus
"OPEN" message is indicated	centrifuge cannot be started - lid supposedly locked	1. lid not properly locked 2 defective lid switch	close lid correctly, if "OPEN" doesn't disappear, call Heraeus
"lid" message iflashes	drive is switched off, unbraked deceleration	1. lid lock was manually opened with draw-thread during operation impermissible manipulation 2 protection circuit has released	1 close lid, then switch mains OFF/ ON, press start to continue centrifugation, else press stop 2. mains OFF, check all vent holes, clean if necessary and wait 15 min., if error happens again, call Heraeus

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## TROUBLE SHOOTING

Trouble		Cause	Remedy
Indication	Symptom		
"ba1" message flashes	centrifuge stops with max braking force	1. rotor is not properly loaded	1. distribute tubes fillings evenly
		2. bucket does not swing out properly or not at all	2. clean trunnions and bucket grooves, grease with lubricant #70006692
		3. sudden tube breakage and shift of center of gravity by escaping fluid	3. remove splinters and debris, clean and reload bucket, restart with reduced speed if necessary
		4. base (floor, table) comes into vibration or centrifuge is not leveled correctly	4. enlarge or change the base and level the centrifuge (drive) or call Heraeus
		5. mechanical or other damage to rotor or bucket	5. <b>do not use this rotor and or bucket again!</b> Call Heraeus
		6. defective drive system (e.g. damaged motor shaft)	6. <b>do not start the drive again!</b> Call Heraeus
"E-03" message flashes	rotor coasts to standstill without braking	speed measuring is disturbed during run	mains OFF/ON, open lid and check rotor locking
"E-07" message flashes	rotor coasts to standstill without braking, cooling unit operates continuously	overtempera- ture alarm, indicated temperature exceeds pre- selected limit of 50°C - pay attention to ambient temperature (35°C max.) and distance to wall (10cm)	let centrifuge cool down, remove samples after standstill and cool them, pre- cool centrifuge and start new run, do not place unit in front of heating apparatus and do not expose the unit to direct sunlight

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## TROUBLE SHOOTING

Trouble		Cause	Remedy
Indication	Symptom		
"E-17" message appears	just after pressing the "lid" key at standstill	1. lid locking mechanism is clamping and doesn't open  2. faulty lid locking assy	1. turn power OFF check locking mechanism, do corrections  2. call Heraeus
"E-18" message flashes	set rotor is #3765 or #8171: drive is cut out in speed range from (0 to 1000) rpm, if the above ment. rotors are really installed, drive will brake, all other rotors will coast to standstill	1. rotor #3765 or # 8171 is preselected but another rotor is physically installed 2. speed is set to low, (1000 rpm or lower) 3. rotors #3765 or #8171 are started with slow acceleration (not allowed for this purpose)	1. to 3: wait until rotor has come to standstill; correct the parameter settings: (rotor's Cat. No., speed or normal acceleration) or install the preselected rotor open and close the lid and restart
"E-21" message flashes	drive is cut out during acceleration and coasts to standstill	insufficient acceleration power (e.g. in case of tremendous imbalance)	wait until rotor has come to standstill, open the lid and e.g. replace all buckets or racks, switch mains OFF/ON and start again
"E-xx" message flashes in speed display possible xx-values 00 - 22	during run: drive shuts off unbraked  at standstill: drive is shut off right after start	release of a thermal or electrical safety trip  troubles in the electronics, measuring or data transmission	switch OFF the centrifuge, wait 15 min switch ON and restart, if error happens again call Heraeus  switch mains OFF/ON, if error does not disappear - call Heraeus

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## 11.2 Servicing

Before calling for service, please make a note of **all** the operating parameters as follows:

- **installed rotor or bucket type (Cat. No.) and loading procedures used**
- the **setting of rotor or bucket Cat. No.**
- the setting of **speed or RCF value**
- the setting of **brake cut-off speed**
- the setting of **time**
- the setting of **temperature (LABOFUGE 400R)**

also the **control procedures** used:

- **start or stop phase**
- **quick run** pressed or released
- normal or slow **acceleration**
- **brake on or off**
- lid **closed or opened**

and describe the **fault signals** exactly as possible e.g.:

- **illuminated or flashing diagnostic LEDs**
- **E-xx - display** (make a note of any error no.)
- **unusual parameter indication** (speed, RCF, time, rotor, temperature)

and the **fault symptoms**:

- **unusual noises** during starting, acceleration, running or braking stages
- **deceleration with or without braking action etc.**

and then call **HERAEUS**.

**ATTENTION - Important note!**

Following switch-on, the centrifuge performs a check of all its display elements, displays its software identification numbers in the "speed" and "time" windows for a short time, and then completes its internal self-checks.

Make a note of these software identification number whenever you contact the **Heræus Service Department**.

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Zentrifuge	Best. Nr.	Spannung	Frequenz	Strom	Leistungsaufnahme	Geräte Sicherungen
Centrifuge	Cat. No.	Voltage	Frequency	Current	Power Input	Unit Fusing <sup>1)</sup>
LABOFUGE 400e	#75008153	230 V	50 Hz	1 A	160 W	2 x 3,15 A
LABOFUGE 400	#75008150	230 V	50 Hz	1,2 A	250 W	2 x 3,15 A
LABOFUGE 400	#75008157	120 V	60 Hz	2,5 A	250 W	2 x 6,3A
LABOFUGE 400R	#75008162	230 V	50 Hz	3,5 A	760 W	2 x 16 A
LABOFUGE 400R	#75008164	120 V	60 Hz	7 A	760 W	2 x 15 A

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## 12. IDENTIFICATION OF UNITS

<sup>1)</sup>ATTENTION! Fuses have slow blow feature and must only be replaced by an authorized service technician!

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## 13. ROTOR / BUCKET -TABLES

<sup>1)</sup> WARNUNG! Automatic imbalance detection is in this case out of function!  
Max. admissible imbalance values must never be exceeded!

<sup>2)</sup> Haematocrit Carrier

<sup>3)</sup> Round bucket with hermetically tight lid (without additional tube sealing risk groups 1 and 2, with hermetically sealed tubes e.g. Oak-Ridge tubes of Comp. Nalge risk groups 3 and 4 can be reached)

**WARNING!** Hermetical sealing is only guaranteed by trouble-free condition of the centrifuge and its accessories!

<sup>4)</sup> Vacutainer (Tubes upto 132 mm length)

<sup>5)</sup> Minimum reachable temperature of samples: 6°C (for all other rotors or buckets of LABOFUGE 400R: < 0°C at maximum rotor/bucket speed)

<sup>6)</sup> 2 Microtiter carriers with two plates in each

**DANGER!** The max. speed of 3000 rpm must never be exceeded!

---

**Heraeus**  
instruments

Heraeus Instruments GmbH  
Laboratory Division  
Postbox 1563  
D-63405 Hanau/Germany  
Telephone: (06181) 35-300  
Fax: (06181) 35 749, Sales  
(05522) 316276, Service

Heraeus Instruments Inc.  
111A - Corporate Blvd.  
South Plainfield, N.J. 07080  
Telephone: (800) 441-2554  
(908) 754-0100  
Telefax: (908) 754-9494

Technische Änderungen  
vorbehalten.  
In the interest of continuous product  
development, we reserve the right  
to make changes without express  
notice.  
Printed in Germany.  
L4e-R\_Z1  
©HS-S/TD-RM 18. Januar 1995