EBA 21

Ⓝ DE	Bedienungsanleitung.....	5
Ⓝ EN	Operating Instructions.....	23
Ⓝ FR	Mode d'emploi	39
Ⓝ IT	Istruzioni per l'uso	57


Fig. 1


Fig. 2 EBA 21

**EG-Konformitätserklärung
EC Conformity Declaration
Déclaration de conformité CE
Dichiarazione di conformità alle norme CEE**

Andreas Hettich GmbH & Co. KG • Föhrenstraße 12 • D-78532 Tuttlingen • Germany

Das bezeichnete Gerät, inklusive Zubehör entspricht den aufgeführten EG-Richtlinien und Normen.
The denoted device, including accessories corresponds to the listed EC guidelines and standards.
L'appareil désigné, y compris les accessoires, correspond aux directives CE et aux normes énumérées.
L'apparecchio designato, compresi gli accessori, è conforme alle direttive CE e alle norme citate.

Geräteart, Type of device, Type d'appareil, Tipo di apparecchio:

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Contents

1	Use according to specification	25
2	Residual risks	25
3	Technical specifications	25
4	Notes on safety	26
5	Symbol meanings	28
6	Delivery checklist	28
7	Initial operation	29
8	Opening and closing the lid.....	29
8.1	Opening the lid	29
8.2	Closing the lid	29
9	Installation and removal of the rotor.....	29
10	Loading the rotor	30
11	Aerosol tight sealing of angle rotors.....	30
12	Control and display elements	30
12.1	Control panel symbols.....	30
12.2	Control panel pushbuttons (keys)	31
12.3	Adjustment possibilities.....	31
13	Entering centrifugation parameter	32
14	Programming.....	32
14.1	Programme input/alteration.....	32
14.2	Programme recall.....	32
15	Centrifugation	32
15.1	Centrifugation with pre-set time	33
15.2	Continuous run.....	33
15.3	Short-term centrifugation.....	33
16	Emergency Stop	33
17	Acoustic Signal	33
18	Recall hours of operation.....	33
19	Relative centrifugal force (RCF)	34
20	Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm ³	34
21	Rotor Identification.....	34
22	Emergency release.....	35
23	Maintenance and servicing	35
23.1	Centrifuge (housing, lid and centrifuging chamber).....	35
23.1.1	Surface cleaning and care	35
23.1.2	Surface disinfection.....	35
23.1.3	Removal of radioactive contaminants	35
23.2	Rotors and Attachments.....	36
23.2.1	Cleaning and care.....	36
23.2.2	Disinfection	36
23.2.3	Removal of radioactive contaminants	36
23.2.4	Rotors and accessories with limited service lives	36
23.3	Autoclaving	37

23.4	Centrifuge containers.....	37
24	Faults.....	37
25	Change mains input fuse.....	38
26	Acceptance of the centrifuges for repair.....	38
27	Disposal.....	38
28	Anhang / Appendix.....	74
28.1	Rotoren und Zubehör / Rotors and accessories.....	74

1 Use according to specification

The machine presented here is a medical product (laboratory centrifuge) according to the IVD guideline 98/79/EG. The centrifuge is used to separate substances or substance mixtures with a density of max. 1.2 kg/dm³. This also includes substances and substance mixtures of human origin. The centrifuge is only intended to be used for this purpose. A different use or application over and above this is deemed not in accordance with the specifications. The company Andreas Hettich GmbH & Co. KG undertakes no liability for damages resulting therefrom.

Belonging to the application according to specification is also the observance of all references contained in the Instruction Manual and compliance with the inspection and maintenance works.

2 Residual risks

The machine is constructed according to the state of the art and the recognized technical safety regulations. Improper use and handling can result in dangers to life and limb of the user or third parties and impairments to the machine or to other material assets. The machine is only to be used for the specified applications and only in an impeccable technical safety condition.

Disturbances that can interfere with the safety are to be immediately rectified.

3 Technical specifications

Manufacturer	Andreas Hettich GmbH & Co. KG D-78532 Tuttlingen			
Model	EBA 21			
Type	1004	1004-30	1004-01	1004-31
Mains voltage ($\pm 10\%$)	220 – 240 V 1~		110 – 127 V 1 ~	
Mains frequency	50 – 60 Hz		50 – 60 Hz	
Connected load	500 VA		550 VA	
Current consumption	2.75 A		5.0 A	
Max. capacity	6 x 50 ml			
Allowed density	1.2 kg/dm ³			
Speed (RPM)	18000			
Force (RCF)	25718			
Kinetic energy	8100 Nm			
Obligatory inspection (BGR 500)	no			
Ambient conditions (EN / IEC 61010-1)	Indoors only Up to 2000 m above sea level 2°C to 40°C Maximum relative humidity 80% for temperatures up to 31°C, linearly decreasing to 50% relative humidity at 40°C.			
– Set-up site				
– Altitude				
– Ambient temperature				
– Humidity				
– Excess-voltage category (IEC 60364-4-443)	II			
– Pollution degree	2			
Device protection class	I			
	Not suitable for use in explosion-endangered areas.			
EMC	EN / IEC 61326-1, Class B		FCC Class B	
– Emitted interference, Interference immunity				
Noise level (dependent on rotor)	≤ 66 dB(A)			
Dimensions				
– Width	275 mm			
– Depth	330 mm			
– Height	247 mm	299 mm	247 mm	299 mm
Weight	approx. 11 kg		approx. 11 kg	

4 Notes on safety

No claim of warranty will be considered by the manufacturer unless ALL instructions in this manual have been followed.

- **The centrifuge should be installed on a good, stable base.**
- **Before using the centrifuge absolutely check the rotor for firm placement.**
- **When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.**
- **Rotors, suspensions and accessories that possess traces of corrosion or mechanical damage or if their term of use has expired may not be used any longer.**
- **The centrifuge may no longer be put into operation when the centrifuging chamber has safety-related damages.**
- **With swing-out rotors the trunnions must be regularly lubricated (Hettich Lubricating Grease No. 4051) in order to ensure consistent swinging out of the hangers.**

- **Before the initial operation of your centrifuge you should read and pay attention to the operating instructions. Only personnel that has read and understood the operating instructions are allowed to operate the device.**
- Along with the operating instructions and the legal regulations on accident prevention, you should also follow the recognised professional regulations for working in a safe and professional manner. These operating instructions should be read in conjunction with any other instructions concerning accident prevention and environmental protection based on the national regulations of the country where the device is to be used.
- This centrifuge is a state-of-the-art piece of equipment which is extremely safe to operate. However, it can lead to danger for users or others if used by untrained staff, in an inappropriate way or for a purpose other than that it was designed for.
- The centrifuge must not be moved or knocked during operation.
- In case of fault or emergency release, never touch the rotor before it has stopped turning.
- To avoid damage due to condensate, when changing from a cold to a warm room the centrifuge must either heat up for at least 3 hours in the warm room before being connected to the mains, or run hot for 30 minutes in the cold room.
- Only the rotors and accessories approved by the manufacturer for this device may be used (see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories").
- The centrifuge rotor may only be loaded in accordance with the chapter "Loading the rotor".
- When centrifuging with maxim revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm³.
- The centrifuge may only be operated when the balance is within the bounds of acceptability.
- The centrifuge may not be operated in explosion-endangered areas.
- The centrifuge must not be used with:
 - inflammable or explosive materials
 - materials that react with one another producing a lot of energy.
- If users have to centrifuge hazardous materials or compounds contaminated with toxic, radioactive or pathogenic micro-organisms, they must take appropriate measures.
 For hazardous substances centrifuge containers with special screw caps must strictly be used. In addition to the screw cap centrifuge containers, for materials in hazard category 3 and 4 a biosafety system must be used (see the World Health Organisation's "Laboratory Biosafety Manual").
 In a biosafety system, droplets and aerosols are prevented from escaping by a bioseal (packing ring).
 If the hanger of a biosafety system is used without the lid, the packing ring must be removed from the hanger in order to prevent the packing ring from being damaged during the centrifugation run. Damaged packing rings must not be used to seal the biosafety system.
 Without the use of a biosafety system the centrifuge is not microbiologically sealed in the sense of the EN / IEC 610101-2-020 standard.

For further details of available biosafety systems see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories". If in doubt, you should obtain relevant information from the manufacturer.

- The centrifuge must not be operated with highly corrosive substances which could impair the mechanical integrity of rotors, hangers and accessories.
- Repairs must only be carried out by personnel authorised to do so by the manufacturer.
- Only original spare parts and original accessories licensed by the Andreas Hettich GmbH & Co. KG company are allowed to be utilised.
- The following safety regulations apply:
EN / IEC 61010-1 and EN / IEC 61010-2-020 as well as their national deviations.
- The safe operation and reliability of the centrifuge can only be guaranteed if:
 - the centrifuge is operated in accordance with the operating instructions,
 - the electrical installation on the site where the centrifuge is installed conforms to the demands of EN / IEC stipulations,
 - prescribed tests to BGV A1, BGR 500 are carried out by an expert.

5 Symbol meanings

Symbol on the machine:

Attention, general hazard area.

Before using the centrifuge implicitly read the operating instructions and pay attention to the safety relevant references!

Symbol in this document:

Attention, general hazard area.

This symbol refers to safety relevant warnings and indicates possibly dangerous situations.

The non-adherence to these warnings can lead to material damage and injury to personal.

Symbol on the machine:

Equipotential: Connector (PE connector) for potential equalization (only for centrifuge with PE connector).

Symbol in this document:

This symbol refers to important circumstances.

Symbol on the machine and in this document:

Symbol for the separate collection of electric and electronic devices according to the guideline 2002/96/EG (WEEE). The device belongs to Group 8 (medical devices).

Applies in the countries of the European Union, as well as in Norway and Switzerland.

6 Delivery checklist

- 1 Connecting cable
- 1 Hex. pin driver
- 1 Operating instructions
- 1 Notes on moving the equipment safely


The rotor(s) and associated accessories are included in the delivery in the quantity ordered.

7 Initial operation

- Remove the transportation safety device from the bottom of the housing, see sheet "Transportation safety device".
- **Position the centrifuge in a stable and level manner in a suitable place. During set-up, the required safety margin of 300 mm around the centrifuge is to be kept according to EN / IEC 61010-2-020.**

When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge. Substances can be discharged through the ventilation opening on the rear side of the centrifuge. The device is to be placed in such a way that the air stream is not directed towards people.

- Do not place any object in front of the ventiduct. Keep a ventilation area of 300 mm around the ventiduct.
- Centrifuge with PE connector:
If required, connect the PE connector at the back of the device with an additional medical potential equalization system.
- Check whether the mains voltage tallies with the statement on the type plate.
- Connect the centrifuge with the connection cable to a standard mains socket. For connection ratings refer to Chapter "Technical specifications".
- Turn on the mains switch. Switch position "I".
The machine type and program version will be displayed and the LEDs light up. **OPEN** **OEFFNEN** will be displayed after 8 seconds. The LED  lights up.
- Open the lid.
The last used centrifuge data will be displayed.
- Remove the transportation safety device in the centrifuge, see sheet "Transportation safety device".

8 Opening and closing the lid

8.1 Opening the lid


The lid can only be opened when the centrifuge is switched on and the rotor is at rest. If it cannot be opened under these circumstances, see the section on "Emergency release".

- Rotate the turning handle on the front panel to the left. The  LED will now go off.
- Open the lid.


8.2 Closing the lid

Do not bang the lid shut.

- Place the lid and lightly press down the front edge of the lid.
- Turn the handle on the front panel to the right. The  LED will come on.

9 Installation and removal of the rotor


- Clean the motor shaft (B) and the collet chuck of the rotor (A) and then lightly grease the motor shaft. Particles of dirt between the motor shaft and the collet chuck prevent perfect positioning of the rotor and cause erratic running.
- Position the rotor vertically on the motor shaft and press down as far as the stop. The rotor must audibly engage in the collet chuck.
- Tighten the rotor tension nut with the supplied wrench by turning in a clockwise direction.
- Check the rotor for firm seating.
- Loosening the rotor: Loosen the tensioning nut by turning anti-clockwise. Turn the tensioning nut until the rotor is able to be lifted from the motor shaft.

10 Loading the rotor


Standard centrifuge containers of glass will not stand RCF values exceeding 4000 (DIN 58970, pg. 2).

- Check the rotor for firm seating.
With swing-out rotors all rotor positions must be lined with **identical** hangers. Certain hangers are marked with the number of the rotor position. These hangers may only be used in the respective rotor position. Hangers that are marked with a set number (e.g. S001/4) may only be used in the set.
- The rotors and hangers may only be loaded symmetrically. The centrifuge containers have to be distributed evenly on all rotor positions. For authorised combinations see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".
In the case of angle rotors all possible rotor positions must be loaded, see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

Rotor is evenly loaded


Not permitted!

Rotor is not evenly loaded

- On certain hangers, the weight of the maximum load or the weight of the maximum load and the maximum weight of the completely loaded hanger is specified. This weight may not be exceeded. In case of exception, see chapter "Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm³". The weight specified for the maximum loading includes the total weight of adapter, frame, centrifuging container and content.
- In containers with rubber inserts, the same number of rubber inserts must always be among the centrifuge containers.
- Always fill the centrifuge containers outside of the centrifuge.
- No liquid should be allowed to enter the centrifugal chamber during filling and swinging out of the hangers.
- The maximum filling quantity for the centrifuge containers specified by the manufacturer must not be exceeded.
- In order to maintain the weight differences within the centrifuge container as marginal as possible, a consistent fill level in the containers is to be heeded.

11 Aerosol tight sealing of angle rotors

To ensure aerosol sealing the lid of an aerosol sealed angular rotor must be firmly closed.

This is performed with the help of the supplied spanner that is inserted through the borehole in the rotary handle. If the rotary handle does not have a borehole then the lid must be firmly closed by turning in clockwise direction.

For available aerosol tight angle rotors, please see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

12 Control and display elements

See figure on page 2.

Fig. 2: Display and control panel

12.1 Control panel symbols

Unbalance display. The unbalance display lights up if the permissible weight difference within the rotor loading has been exceeded.


Rotation display. The rotation display lights up during the centrifugation run as long as the rotor is turning.

Lid display. The lid display lights up if the lid is closed and locked. During the centrifugation run the lid display is extinguished.

Operating errors or any occurring disturbances are symbolised in the display (see Chapter "Faults").

12.2 Control panel pushbuttons (keys)

-  Selection control key for selection of specific parameter.
The subsequent parameter is selected by every further keystroke.
-  By pushing the arrow keys the value of a previously selected parameter is increased or decreased. By keeping an arrow key depressed the value is increased or decreased with increasing speed.
-  Start the centrifugation run. The rotation display  lights up.
Acquisition of input and changes.
-  Stop the centrifugation run.
The rotor runs down with pre-selected brake step. The LED in the key lights up until the rotor stops.
Pushing the key twice triggers the EMERGENCY STOP.
-  Selecting the RCF display.
The LED in the key lights up.
-  Short-term centrifugation.
The centrifugation run is effected as long as the key is held down.
-  Storing and recalling programmes.
The subsequent programme is displayed by every further keystroke (1 - 2 - 3 - #).

12.3 Adjustment possibilities

- t/min Running time. Adjustable from 1 - 99 min. in 1-min. intervals.
- t/sec Running time. Adjustable from 1 - 59 sec. in 1-sec. intervals.

Continuous run " ∞ ". Set parameter **t/min** and **t/sec** to zero.
- RPM Revolutions per minute. A numerical value from 500 RPM up to the maximum speed of the rotor can be set. Maximum speed of the rotor, see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories". Adjustable from 500 RPM to 10000 RPM in steps of 10 and from 10000 RPM up to the maximum speed in steps of 100.
- RAD/mm Centrifugation radius. Input in mm. For centrifugation radius see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories". The input of the radius is only possible if the RCF display is selected (LED in the key  lights up).
- RCF Relative centrifugal force. A numerical value can be set, which gives a speed between 500 RPM and the maximum speed of the rotor. Adjustable up to 100 in intervals of 1, and from 100 in intervals of 10. The RCF value is automatically rounded up or rounded down with regard to the RPM interval. The input of the RCF is only possible if the RCF display is selected (LED in the key  lights up).
-  Starting steps 1 - 9. Step 9 = shortest starting time, Step 1 = longest starting time.
-  Brake steps 0 - 9. Step 9 = shortest run-down time, Step 1 = long run-down time, Step 0 = longest run-down time (brakeless run-down).
- PROG Programme Number. 3 programmes can be stored (programme positions 1 - 2 - 3). The programme position # serves as temporary storage for altered adjustments.

13 Entering centrifugation parameter

If no key is pressed for 8 seconds long after the selection or during the input of parameters, the previous values will be shown in the display. The input of parameter then has to be executed again.

- Select the RPM or RCF display with the key \odot RCF. If the RCF display is selected, the LED in the key \odot RCF lights up.
- Select the desired parameter with the key \leftarrow , and adjust with the keys \blacktriangle \blacktriangledown . In order to set continuous operation, the parameters **t/min** and **t/sec** must be set to zero with the \blacktriangledown key. Continual running is represented in the display by the following symbol, " ∞ ".
- After input of all parameters, press the key \odot START in order to store the adjustments on the programme position #. As confirmation, ***** ok ***** will be displayed for a short period.

The data on the programme position # will be overwritten with every input of parameters and pressing of the key \odot START.

14 Programming

14.1 Programme input/alteration

- Select the RPM or RCF display with the key \odot RCF. If the RCF display is selected, the LED in the key \odot RCF lights up.
- Select the desired parameter with the key \leftarrow , and adjust with the keys \blacktriangle \blacktriangledown .
- Select the parameter **PROG** with the key \leftarrow , and set the desired programme position with the keys \blacktriangle \blacktriangledown . The LED in the key \odot PROG lights up.
- Press the key \odot PROG in order to store the adjustments to the desired programme position. As confirmation, ***** ok ***** will be displayed for a short period.

14.2 Programme recall

- Select the desired programme position by pressing the key \odot PROG. The centrifugation data of the selected programme position will be displayed.
- The parameters can be checked by pressing the key \leftarrow .

15 Centrifugation

When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.

If the permissible weight difference within the rotor loading has been exceeded, the drive shuts down during the start-up, the unbalance display lights up, and **IMBALANCE** is displayed.

If the speed in the selected program is higher than the maximum speed of the rotor, it will not be possible to start a centrifugation run. **N > ROTOR MAX** will be displayed (see Chapter "Faults").

A centrifugation run can be stopped at any time by pushing the key \odot STOP.

All parameters can be selected and altered during the centrifugation run (see Chapter "Entering centrifugation parameter").

You can switch-over at any time between the RPM and RCF display with the key \odot RCF. The input of the centrifugation radius is necessary if you are working with the RCF display.

If \blacktriangle **OPEN** \blacktriangle **OEFFNEN** is displayed, a further operation of the centrifuge is only possible after opening the lid once.

If **R xx n-max xxxxx** is displayed, then no centrifugation run has taken place as the rotor was changed beforehand, refer to Chapter "Rotor Identification".

- Turn on the mains switch. Switch position I.
- Load the rotor and close the centrifuge lid.

15.1 Centrifugation with pre-set time

- Adjusting time or recall a programme with pre-set time (see Chapter "Entering centrifugation parameter" or "Programme recall").
- Press the key **START**. The rotation display  lights up as long as the rotor is turning.
- After expiration of the time or with truncation of the centrifugation run by pushing the key **STOP**, the run-down is effected with the selected brake step. The brake step is displayed.


During the centrifugation run the rotational speed of the rotor or the subsequently resulting RCF value, the sample temperature (only in centrifuges with cooling) and the remaining time will be displayed.

15.2 Continuous run

- Adjusting the symbol ∞ or recall a continuous run programme (see Chapter "Entering centrifugation parameter" or "Programme recall").
- Press the key **START**. The rotation display  lights up as long as the rotor is turning. The time metering begins at 00:00.
- Press the key **STOP** in order to stop the centrifugation run. The run-down is effected with the selected brake step. The brake step is displayed.

During the centrifugation run the rotational speed of the rotor or the subsequently resulting RCF value, the sample temperature (only in centrifuges with cooling) and the expired time will be displayed.

15.3 Short-term centrifugation

- Hold down the key **IMPULS**. The rotation display  lights up as long as the rotor is turning.
- Let go of the key **IMPULS** again in order to stop the centrifugation run. The run-down is effected with the selected brake step. The brake step is displayed.

During the centrifugation run the rotational speed of the rotor or the subsequently resulting RCF value, the sample temperature (only in centrifuges with cooling) and the expired time will be displayed.

16 Emergency Stop

- Press the key **STOP** twice.

With Emergency Stop the run-down is effected with brake step 9 (shortest run-down time). Brake step 9 is displayed. If brake step 0 was pre-selected, the run-down time is technically longer than with brake step 9.

17 Acoustic Signal

The acoustic signal sounds:

- Upon the appearance of a disturbance in 2 second intervals.
- After completion of a centrifugation run and rotor standstill in 30 second intervals.




The acoustic signal is stopped by opening the lid or pressing any key.

The signal after completion of the centrifugation run can be activated or deactivated in the following manner, if the rotor is at standstill:

- Hold down the key  for 8 seconds.
After 8 seconds, **SOUND / BELL** appears in the display.
- Set **OFF** or **ON** with the key  or .
- Press the key **START** in order to store the setting.
As confirmation, ***** ok ***** will be displayed for a short period.

18 Recall hours of operation

Recall hours of operation is only possible during rotor standstill.

- Hold down the key  for 8 seconds.
After 8 seconds, **SOUND / BELL** appears in the display.
- Press the key  once again.
The centrifuge's hours of operation (**CONTROL:**) are displayed.
- Press the key  or  to exit the hours of operation recall.

19 Relative centrifugal force (RCF)

The relative centrifugal force (RCF) is given as a multiple of the acceleration of gravity (g). It is a unit-free value and serves to compare the separation and sedimentation performance.

These values are calculated using the formula below:

$$\text{RCF} = \left(\frac{\text{RPM}}{1000} \right)^2 \times r \times 1,118 \quad \Rightarrow \quad \text{RPM} = \sqrt{\frac{\text{RCF}}{r \times 1,118}} \times 1000$$

RCF = relative centrifugal force

RPM = rotational speed (revolutions per minute)

r = centrifugal radius in mm = distance from the centre of the turning axis to the bottom of the centrifuge.

For more on the centrifugal radius see the chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

The relative centrifugal force (RCF) stands in relation to the revolutions per minute and the centrifugal radius.

20 Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm³

When centrifuging with maximum revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm³.

The speed must be reduced for materials or mixtures of materials with a higher density.

The permissible speed can be calculated using the following formula:

$$\text{Reduced speed (n}_{\text{red}}) = \sqrt{\frac{1.2}{\text{Greater density [kg/dm}^3]}} \times \text{maximum speed [RPM]}$$

e.g.: maximum speed RPM 4000, density 1.6 kg/dm³

$$n_{\text{red}} = \sqrt{\frac{1.2 \text{ kg/dm}^3}{1.6 \text{ kg/dm}^3}} \times 4000 \text{ RPM} = 3464 \text{ RPM}$$

In the exceptional case that the maximum loading indicated on the hanger is exceeded, the speed must also be reduced.

The permissible speed can be calculated using the following formula:

$$\text{Reduced speed (n}_{\text{red}}) = \sqrt{\frac{\text{maximum load [g]}}{\text{actual load [g]}}} \times \text{maximum speed [RPM]}$$

e.g.: maximum speed RPM 4000, maximum load 300 g, actual load 350 g

$$n_{\text{red}} = \sqrt{\frac{300 \text{ g}}{350 \text{ g}}} \times 4000 \text{ RPM} = 3703 \text{ RPM}$$

If in doubt you should obtain clarification from the manufacturer.

21 Rotor Identification

After every start of a centrifugation run the rotor utilised is identified.

After a change of rotor the drive switches off and the rotor code (R xx) as well as the maximum rotational speed (n-max=xxxxx) of the rotor are displayed.

A further operation of the centrifuge is only possible after a single opening of the lid.

If, following a rotor change, the maximum speed of the rotor is less than the set speed, the speed is limited to the maximum speed of the rotor.

22 Emergency release

The lid cannot be opened in the event of a power failure. An emergency release has to be executed by hand.

For emergency release disconnect the centrifuge from the mains.
Open the lid only during rotor standstill.

See figure on page 2.

- Switch off the mains switch (switch position "0").
- Look through the window in the lid to be sure that the rotor has come to a standstill.
- Position the centrifuge on the table sufficiently far forward so that the hole (Fig. 1, A) for the emergency release on the base of the casing can be reached.
- Insert the hexagonal pin driver into this hole, push it upwards and at the same time rotate the lid-lock handle to the left.
- Open the lid.

23 Maintenance and servicing

Pull the mains plug before cleaning.

Before any other cleaning or decontamination process other than that recommended by the manufacturer is applied, the user has to check with the manufacturer that the planned process does not damage the device.

- Centrifuges, rotors and accessories must not be cleaned in rinsing machines.
- They may only be cleaned by hand and disinfected with liquids.
- The water temperature must be between 20 – 25°C.
- Only detergents/disinfectants may be used which:
 - have a pH between 5 - 8
 - do not contain caustic alkalis, peroxides, chlorine compounds, acids and alkaline solutions
- In order to prevent appearances of corrosion through cleaning agents or disinfectants, the application guide from the manufacturer of the cleaning agent or disinfectant are absolutely to be heeded.

23.1 Centrifuge (housing, lid and centrifuging chamber)

23.1.1 Surface cleaning and care

- Clean the centrifuge housing and the centrifuging chamber regularly, using soap or a mild detergent and a damp cloth if required. For one thing, this services purposes of hygiene, and it also prevents corrosion through adhering impurities.
- Ingredients of suitable detergents: soap, anionic tensides, non-ionic tensides.
- After using detergents, remove the detergent residue by wiping with a damp cloth.
- The surfaces must be dried immediately after cleaning.
- In the event of condensation water formation, dry the centrifugal chamber by wiping out with an absorbent cloth.
- Lightly rub the rubber seal of the centrifuge chamber with talcum powder or a rubber care product after each cleaning.
- The centrifuging chamber is to be checked for damage once a year.

If damage is found which is relevant to safety, the centrifuge may no longer be put into operation. In this case, notify Customer Service.

23.1.2 Surface disinfection

- If infectious materials penetrates into the centrifugal chamber this is to be disinfected immediately.
- Ingredients of suitable disinfectants: ethanol, n-propanol, isopropyl alcohol, glutardialdehyde, quaternary ammonium compounds.
- After using disinfectants, remove the disinfectant residue by wiping with a damp cloth.
- The surfaces must be dried immediately after disinfecting.

23.1.3 Removal of radioactive contaminants

- The agent must be specifically labelled as being an agent for removing radioactive contaminants.
- Ingredients of suitable agents for removing radioactive contaminants: anionic tensides, non-ionic tensides, polyhydrated ethanol.
- After removing the radioactive contaminants, remove the agent residue by wiping with a damp cloth.
- The surfaces must be dried directly after removing the radioactive contaminants.

23.2 Rotors and Attachments

23.2.1 Cleaning and care

- In order to prevent corrosion and material changes, rotors and accessories must be cleaned regularly with soap or a mild detergent and a damp cloth. Cleaning is recommended at least once a week. Contaminants must be removed immediately.
- Ingredients of suitable detergents:
soap, anionic tensides, non-ionic tensides.
- After using detergents, remove detergent residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotors and accessories must be dried directly after cleaning.
- Angle rotors, container and hanger made of aluminium are to be lightly greased after drying using acid-free grease, e.g. vaseline.
- In the case of biosafety systems (for further details of available biosafety systems see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories"), the packing rings must be checked and cleaned regularly (weekly). The sealing ring is to be replaced immediately upon indication of crack formation, embrittlement or abrasive wear. In order to prevent the packing ring from twisting when opening and closing the cover, the packing ring must be lightly rubbed with talcum powder or a rubber care product.
- In order to prevent corrosion as a result of moisture between the rotor and the motor shaft, the rotor should be disassembled and cleaned at least once a month, and the motor shaft should be lightly greased.
- The rotors and the accessory parts have to be checked on a monthly basis for wear and corrosion.

Rotors and attachments may no longer be utilised upon indication of wear and tear or corrosion.

- Check the firm seating of the rotor on a weekly basis.

23.2.2 Disinfection

- If infectious material should get on the rotors or accessories, they must be appropriately disinfected.
- Ingredients of suitable disinfectants:
glutaraldehyde, propanol, ethyl hexanol, anionic tensides, corrosion inhibitors.
- After using disinfectants, remove disinfectant residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotors and accessories must be dried directly after disinfection.

23.2.3 Removal of radioactive contaminants

- The agent must be specifically labelled as being an agent for the removal of radioactive contaminants.
- Ingredients of suitable agents for removing radioactive contaminants:
anionic tensides, non-ionic tensides, polyhydrated ethanol.
- After removing the radioactive contaminants, remove agent residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotors and accessories must be dried directly after removing the radioactive contaminants.

23.2.4 Rotors and accessories with limited service lives

The use of certain rotors, hangers and accessory parts is limited by time.

These are marked with the maximum permitted number of operating cycles or with an expiration date and the maximum permitted number of operating cycles or just with the expiration date; e.g.:

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For safety reasons, rotors, hangers and accessory parts may no longer be used if either the indicated maximum number of operating cycles or the indicated expiration date has been reached.

23.3 Autoclaving

Swing-out rotors, angle rotors made of aluminium, suspension made of metal, lids with bio-containment as well as stands and reductions can be autoclaved at 121° C / 250°F (20 mins).

Otherwise you must ask the manufacturer.

No statement can be made about the degree of sterility.

The lids of the rotors and containers must be removed prior to autoclaving.

Autoclaving accelerates the ageing process of plastics. In addition, autoclaving may discolour plastics.

We recommend that the packing rings of the bio-safety system be replaced after autoclaving.

23.4 Centrifuge containers

- With leakiness or after the breakage of centrifuging containers broken container parts and leaked centrifugation material are to be completely removed.
- The rubber inserts as well as the plastic sleeves of the rotors are to be replaced after a glass breakage.

Remaining glass splitters cause further glass breakage!

- If this concerns infectious material, a disinfection process is to be executed immediately.

24 Faults

If the fault cannot be eliminated with the help of the fault table, please inform Customer Service.

Please specify the type of centrifuge and the serial number. Both numbers can be found on the name plate of the centrifuge.

Perform a MAINS RESET:


- Switch off the mains switch (switch position "0").
- Wait at least 10 seconds and then switch on the mains switch again (switch position "1").

Message / fault		Cause	Remedy
No display	---	No voltage. Mains input fuses defective.	<ul style="list-style-type: none"> – Check supply voltage. – Check mains power input fuse, refer to Chapter "Change mains input fuse". – Mains switch ON.
TACHO - ERROR	01	Faulty speedometer.	<ul style="list-style-type: none"> – Open the cover. – Switch off the mains switch (switch position "0"). – Wait at least 10 seconds. – Turn the rotor vigorously by hand. – Switch on the mains switch again (switch position "1"). The rotor must turn during switch-on.
	02	No rotor installed. Defective motor, frequency converter or drive.	
IMBALANCE	---	The rotor is unevenly loaded.	<ul style="list-style-type: none"> – Open lid. – Check the loading of the rotor, see chapter "Loading the rotor". – Repeat the centrifugation run.
CONTROL - ERROR	04, 06 - 09	Error in lid locking or lid closure.	<ul style="list-style-type: none"> – Perform a MAINS RESET.
N > MAX	05	Rotation too fast	
N < MIN	13	Rotation too slow	
ROTORCODE	10	Incorrect rotor coding	
MAINS INTERRUPT	---	Power failure during the centrifugation run. (The centrifugation run was not finished.)	<ul style="list-style-type: none"> – Open lid. – Push START button. – Repeat the centrifugation run if necessary.
VERSIONS-ERROR	12	Mismatch between electronic components	<ul style="list-style-type: none"> – Perform a MAINS RESET.
CONTROL-ERROR	21 - 27	Error / defective control unit	
SER I/O - ERROR	30 - 38	Error / defective interface	
FU / CCI - ERROR	60 - 83	Error / defective motor control	
N > ROTOR MAX	---	Speed in the selected program greater than the maximum speed of the rotor.	<ul style="list-style-type: none"> – Check the set speed. – Reduce the set speed

25 Change mains input fuse

Switch off the mains switch and separate the centrifuge from the mains!

The fuse holder (A) with the mains input fuses is located next to the mains switch.

- Remove the connecting cable from the machine plug socket.
- Press the snap-fit (B) against the fuse holder (A) and remove.
- Exchange defective mains input fuses.

Only use fuses with the rating defined for the type. See the following table.

- Reinsert the fuse holder until the snap-fit clicks shut.
- Reconnect the centrifuge to the mains supply.

Model	Type	Fuse	Order no.
EBA 21	1004, 1004-30	T 3,15 AH/250V	E997
EBA 21	1004-01, 1004-31	T 6,3 AH/250V	2266

26 Acceptance of the centrifuges for repair

If the centrifuge is returned to the manufacturer for repair, it must be decontaminated and cleaned to protect persons, environment and material.

We reserve the right to accept contaminated centrifuges.

Costs incurred for cleaning and disinfection are to be charged to the customer.

We ask for your understanding in this matter.

27 Disposal

Before disposal, the device must be decontaminated and cleaned to protect people, the environment and property.

When you are disposing of the device, the respective statutory rules must be observed.

Pursuant to guideline 2002/96/EC (WEEE), all devices supplied after August 13, 2005 may not be disposed as part of domestic waste. The device belongs to group 8 (medical devices) and is categorized in the business-to-business field.

The icon of the crossed-out trash can shows that the device may not be disposed as part of domestic waste.

The waste disposal guidelines of the individual EC countries might vary. If necessary, contact your supplier.

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




1) Sample temp. rise during maximum speed and 1 hour running time


1416		1054-A								
Winkelrotor 12-fach / Angle rotor 12-times  35°										
		0507	0553	---	---	---	0501			
										
Kapazität / capacity	ml	15	5	1,6 – 5,0	4	6	1,1 – 1,4	2,7 - 3	2,6 – 2,9	
Maße / dimensions Ø x L	mm	17 x 100	12 x 75 13 x 75	13 x 75	12 x 60	12 x 82	8 x 66	11 x 66	13 x 65	
Anzahl p. Rotor / number p. rotor		12	12	12	12	12	12	12	12	
Drehzahl / speed	RPM	6000	6000	6000	6000	6000	6000	6000	6000	
RZB / RCF	⁴⁾	4146	3300	3300	3260	3300	3300	3300	3300	
Radius / radius	mm	103	82	82	81	82	82	82	82	
 9 (97%)	sec	18								
 9	sec	23								
Probenerwärmung/Sample temp. rise	K ¹⁾	12								

1416											
Winkelrotor 12-fach / Angle rotor 12-times  35°											
		0509							0518		
											
Kapazität / capacity	ml	15	4,9	4,5 - 5	7,5 – 8,2	9 - 10	10	4 - 7	8	8,5 - 10	15
Maße / dimensions Ø x L	mm	17 x 120	13 x 90	11 x 92	15 x 92	16 x 92	15 x 102	13 x 100	16 x 125	16 x 100	17 x 100
Anzahl p. Rotor / number p. rotor		6	12	12	12	12	12	12	12	12	
Drehzahl / speed	RPM	6000	6000	6000	6000	6000	6000	6000	6000	6000	
RZB / RCF	²⁾	4146	4146	4146	4146	4146	4146	4146	4146	4146	
Radius / radius	mm	103	103	103	103	103	103	103	103	103	
 9 (97%)	sec	18									
 9	sec	23									
Probenerwärmung/Sample temp. rise	K ¹⁾	12									

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









- 1) Sample temp. rise during maximum speed and 1 hour running time
- 2) Centrifuge containers of glass will not stand RCF values exceeding 4000

1022		3) 2031		2023		2024				
Winkelrotor 12-fach / Angle rotor 12-times   45°										
		0536	2078							
										
Kapazität / capacity	ml	2,0	1,5	0,5	0,8	0,4	0,2			
Maße / dimensions	∅ x L mm	11 x 38	8 x 30	8 x 45	6 x 45	6 x 18				
Anzahl p. Rotor / number p. rotor		12	12	12	12	12				
Drehzahl / speed	RPM	18000	18000	18000	18000	18000				
RZB / RCF	²⁾	25718	25718	25718	25718	25718				
Radius / radius	mm	71	71	71	71	71				
 9 (97%)	sec						16			
 9	sec						19			
Probenerwärmung/Sample temp. rise	K ^{1) 4)}						18			

1024		3) 2031		2023		2024				
Winkelrotor 24-fach / Angle rotor 24-times  +  6) 40										
		0536	2078							
										
Kapazität / capacity	ml	2,0	1,5	0,8	0,5	0,4	0,2			
Maße / dimensions	∅ x L mm	11 x 38	11 x 38	8 x 45	8 x 30	6 x 45	6 x 18			
Anzahl p. Rotor / number p. rotor		24	24	24	24	24	24			
Drehzahl / speed	RPM	15000	15000	15000	15000	15000	15000			
RZB / RCF	²⁾	21382	21382	21382	21382	21382	21382			
Radius / radius	mm	85	85	85	85	85	85			
 9 (97%)	sec						30			
 9	sec						31			
Probenerwärmung/Sample temp. rise	K ^{1) 4)}						17			

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




- 1) Sample temp. rise during maximum speed and 1 hour running time
- 2) If the closed rubber seal (order no. E612) is installed, the sample temperature rise will increase.
- 3) If the rubber seal with ventiduct (order no. E1877) is installed, only centrifuge containers made of plastic may be used.
- 4) recommended for high-speed centrifugation
- 5) Centrifuge containers of glass will not stand RCF values exceeding 4000
- 6) in conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".




1115		1123									
Ausschwingrotor 6-fach / Swing out rotor 6-times  $\angle 40^\circ$											
		0716								0716	
		2 x 								1 x 	
		0507	2079							0518	
											
Kapazität / capacity	ml	15	10	4,9	4 – 5,5	7,5 – 8,2	9 - 10	8,5 - 10	4 - 7	15	
Maße / dimensions	Ø x L mm	17 x 100	17 x 70	13 x 90	15 x 75	15 x 92	16 x 92	16 x 100	16 x 75	17 x 100	
Anzahl p. Rotor / number p. rotor		6	6	6	6	6	6	6	6	6	
Drehzahl / speed	RPM	5000	5000	5000	5000	5000	5000	5000	5000	5000	
RZB / RCF		2711	2376	2711	2711	2711	2711	2711	2543	2711	
Radius / radius	mm	97	85	97	97	97	97	97	91	97	
 9 (97%)	sec	8									
 9	sec	12									
Probenerwärmung/sample temp. rise	K ¹⁾	9									

1115		1123							
Ausschwingrotor 6-fach / Swing out rotor 6-times  $\angle 40^\circ$									
		1054-A							
									
		0553	0501						0701
									
Kapazität / capacity	ml	5	6	1,1 – 1,4	2,7 - 3	2,6 – 2,9	1,6 - 5	5	4
Maße / dimensions	Ø x L mm	12 x 75	12 x 82	8 x 66	11 x 66	13 x 65	13 x 75	13 x 75	12 x 60
Anzahl p. Rotor / number p. rotor		6	6	6	6	6	6	6	6
Drehzahl / speed	RPM	5000	5000	5000	5000	5000	5000	5000	5000
RZB / RCF		2068	2068	2068	2068	2068	2068	2068	1984
Radius / radius	mm	74	74	74	74	74	74	74	71
 9 (97%)	sec	8							
 9	sec	12							
Probenerwärmung/sample temp. rise	K ¹⁾	9							

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All rights reserved. No part of this publication may be reproduced without the written prior permission of the copyright owner.