ROTOLAVIT

Ⓝ	Bedienungsanleitung	6
Ⓝ	Operating Instructions	37
Ⓝ	Mode d'emploi	68
Ⓝ	Istruzioni per l'uso	99


Fig. 1


Fig. 2 ROTOLAVIT


Fig. 3


Fig. 4

Hettich AG, CH-8806 Baech, Switzerland
Seestrasse 204a
Tel. +41 44 786 80 20, Fax +41 44 786 80 21
mail@hettich-ag.ch
www.hettich-ag.ch

KONFORMITÄTSERKLÄRUNG / **D**ECLARATION OF CONFORMITY
DECLARATION DE CONFORMITE / **D**ICHIARAZIONE DI CONFORMITA

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We declare under our sole responsibility
that the in vitro diagnostic medical
device

Cell washing centrifuge Hettich Rotolavit from serial number 1000100 to ...
Cell washing centrifuge Hettich Rotolavit from serial number 2000100 to ...
Cell washing centrifuge Medion DACIII from serial number 4000100 to ...

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classified as follows according to the directive on in vitro diagnostic medical devices 98/79/EC

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meets all the provisions of the directive on in vitro diagnostic medical devices 98/79/EC which apply to it.

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Applied common technical specifications, harmonised standards, national standards or other normative documents

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A handwritten signature in blue ink, appearing to read "P. Lutz", is positioned above the printed name and title.

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www.hettich-ag.ch

© 2003 by Hettich AG

All rights reserved. No part of this publication may be reproduced without the written prior permission of the copyright owner.

Änderungen vorbehalten! , Modifications reserved! , Sous réserve de modifications ! , Con riserva di modifiche!

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Contents

1	Use according to specification	39
2	Residual risks	39
3	Technical specifications	39
4	Notes on safety	40
5	Symbol meanings	41
6	Delivery checklist	41
7	Unpacking the centrifuge	41
8	Initial operation	42
9	3-way valve (only for centrifuges with 3-way valves)	43
10	Opening and closing the lid	43
10.1	Opening the lid	43
10.2	Closing the lid	43
11	Inserting and removing the rotor	43
12	Loading the rotor	44
13	Control and display elements	44
13.1	Control panel symbols	44
13.2	Control panel buttons	44
13.3	Possible settings	45
13.3.1	Global parameters	45
13.3.2	Program parameters (parameters for the programs 1 to 5)	45
13.3.3	Program parameters for the centrifugation program Spin	46
14	Programming	47
14.1	Global parameters adjustment	47
14.2	Entering Spin centrifugation program	47
14.3	Entering program parameters	47
14.4	Programme recall	48
15	Functions	48
15.1	Description	48
15.1.1	Washing process	49
15.1.2	Suspension	50
15.1.3	Centrifugation run for acceleration of the agglutination reaction	50
15.1.4	Centrifugation program Spin	51
15.1.5	Cleaning program (Clean Program)	51
15.2	Combining functions	52
15.2.1	Performing washing and suspension procedure	52
15.2.2	Performing centrifugation run for acceleration of the agglutination reaction	54
15.2.3	Production of a suspension without previous washing procedure	54
15.2.4	Checking the samples during the washing procedure or production of the suspension	55
15.2.5	Performing centrifugation run with Spin centrifugation program	55
15.2.6	System with rinsing cleaning program (Clean Program)	55
15.3	Examples for the function selection	57
15.3.1	Washing with subsequent centrifugation with shaking	57
15.3.2	Washing with subsequent suspension without shaking the suspension	57

15.3.3	Washing with subsequent suspension with shaking the suspension	58
15.3.4	Only suspension without shaking.....	58
15.3.5	Only suspension with shaking.....	58
15.3.6	Centrifugation run with decanting	59
15.3.7	Centrifugation run with SPIN (S) program.....	59
16	Emergency stop	59
17	Adjusting the filling volume.....	59
18	Entering the rotor type.....	60
19	Audible signal.....	60
20	Querying operating hours.....	60
21	Relative centrifugal force (RCF)	61
22	Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm ³	61
23	Emergency release	61
24	Maintenance and servicing	62
24.1	Centrifuge	62
24.2	Rotor	63
24.3	Autoclaving	63
24.4	Glass breakage.....	63
24.5	Clean system with cleaning solution	64
24.5.1	Device without 3-way valve.....	64
24.5.2	Device with 3-way valve:.....	64
24.6	Removal of the drain gutter.....	64
25	Faults	65
25.1	Operating faults.....	65
25.2	Error messages.....	66
26	Changing the fuse	67
27	Acceptance of the centrifuges for repair.....	67
28	Disposal	67
29	Anhang / Appendix.....	130
29.1	Rotoren und Zubehör / Rotors and accessories.....	130

1 Use according to specification

The machine presented here is a medical product (laboratory centrifuge) according to the IVD guideline 98/79/EG. This appliance is for cleaning erythrocytes and for carrying out quick anti-human globulin tests (direct and indirect Coombs' tests) when cross-matching, looking for and differentiating anti-bodies as well as when determining features of rare blood groups.

During a washing process, which consists of several wash cycles, the physiological saline solution is added successively, the erythrocytes are sedimented by the centrifugation process, and finally the physiological saline solution is decanted again.

After manually adding the anti-human globulin serum, a further centrifugation speeds up the agglutination test.

The device is only intended for the designated purpose.

A different use or application over and above this is deemed not in accordance with the specifications. The company Hettich AG undertakes no liability for damages resulting therefrom.

Belonging to the application according to specification is also the observance of all references contained in the Instruction Manual and compliance with the inspection and maintenance works.

2 Residual risks

The machine is constructed according to the state of the art and the recognized technical safety regulations. Improper use and handling can result in dangers to life and limb of the user or third parties and impairments to the machine or to other material assets. The machine is only to be used for the specified applications and only in an impeccable technical safety condition.

Disturbances that can interfere with the safety are to be immediately rectified.

3 Technical specifications

Manufacturer	Hettich AG CH-8806 Baech	
Model	ROTOLAVIT	
Type	1006	1006-01
Mains voltage (± 10%)	200 – 240 V 1~	100 – 120 V 1 ~
Mains frequency	50 – 60 Hz	50 – 60 Hz
Connected load	180 VA	180 VA
Current consumption	0.75 A	1.5 A
Max. capacity	24 x 5 ml	
Allowed density	1.2 kg/dm ³	
Speed (RPM)	3500	
Force (RCF)	1438	
Kinetic energy	250 Nm	
Obligatory inspection (BGR 500)	no	
Ambient conditions (EN / IEC 61010-1)	Indoors only Up to 2000 m above sea level 2°C to 40°C Maximum relative humidity 80% for temperatures up to 31°C, linearly decreasing to 50% relative humidity at 40°C.	
– Set-up site		
– Altitude		
– Ambient temperature		
– Humidity		
– Excess-voltage category (IEC 60364-4-443)	II	
– Pollution degree	2	
Device protection class	I	
Not suitable for use in explosion-endangered areas.		
EMC		
– Emitted interference, Interference immunity	EN / IEC 61326-1, Class B	FCC Class B
Noise level (dependent on rotor)	≤ 44 dB(A)	
Dimensions		
– Width	333 mm	
– Depth	420 mm	
– Height	278 mm	
Weight	23.5 kg	

4 Notes on safety

No claim of warranty will be considered by the manufacturer unless ALL instructions in this manual have been followed.

- **The centrifuge should be installed on a good, stable base.**
- **When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.**
- **Rotors, suspensions and accessories that possess traces of corrosion or mechanical damage or if their term of use has expired may not be used any longer.**
- **The centrifuge may no longer be put into operation when the centrifuging chamber has safety-related damages.**

- **Before the initial operation of your centrifuge you should read and pay attention to the operating instructions. Only personnel that has read and understood the operating instructions are allowed to operate the device.**
- Along with the operating instructions and the legal regulations on accident prevention, you should also follow the recognised professional regulations for working in a safe and professional manner. These operating instructions should be read in conjunction with any other instructions concerning accident prevention and environmental protection based on the national regulations of the country where the device is to be used.
- This centrifuge is a state-of-the-art piece of equipment which is extremely safe to operate. However, it can lead to danger for users or others if used by untrained staff, in an inappropriate way or for a purpose other than that it was designed for.
- The centrifuge must not be moved or knocked during operation.
- In case of fault or emergency release, never touch the rotor before it has stopped turning.
- To avoid damage due to condensate, when changing from a cold to a warm room the centrifuge must either heat up for at least 3 hours in the warm room before being connected to the mains, or run hot for 30 minutes in the cold room.
- Only the rotors and accessories approved by the manufacturer for this device may be used (see section "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories").
- The centrifuge rotor may only be loaded in accordance with the chapter "Loading the rotor".
- When centrifuging with maxim revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm³.
- The centrifuge may only be operated when the balance is within the bounds of acceptability.
- The centrifuge may not be operated in explosion-endangered areas.
- The centrifuge must not be used with:
 - inflammable or explosive materials
 - materials that react with one another producing a lot of energy.
- If users have to centrifuge hazardous materials or compounds contaminated with toxic, radioactive or pathogenic micro-organisms, they must take appropriate measures.
- The centrifuge must not be operated with highly corrosive substances which could impair the mechanical integrity of rotors, hangers and accessories.
- Repairs must only be carried out by personnel authorised to do so by the manufacturer.
- Only original spare parts and original accessories licensed by the Hettich AG company are allowed to be utilised.
- Components contaminated with blood (e.g. rotor, centrifugation area) must after replacement be disposed of in the special waste for materials contaminated with blood.
- The following safety regulations apply:
EN / IEC 61010-1 and EN / IEC 61010-2-020 as well as their national deviations.

- The safe operation and reliability of the centrifuge can only be guaranteed if:
 - the centrifuge is operated in accordance with the operating instructions,
 - the electrical installation on the site where the centrifuge is installed conforms to the demands of EN / IEC stipulations,
 - prescribed tests to BGV A1, BGR 500 are carried out by an expert.

5 Symbol meanings

Symbol on the machine:

Attention, general hazard area.

Before using the centrifuge implicitly read the operating instructions and pay attention to the safety relevant references!

Symbol in this document:

Attention, general hazard area.

This symbol refers to safety relevant warnings and indicates possibly dangerous situations.

The non-adherence to these warnings can lead to material damage and injury to personal.

Symbol in this document:

This symbol refers to important circumstances.

Symbol on the machine and in this document:

Symbol for the separate collection of electric and electronic devices according to the guideline 2002/96/EG (WEEE). The device belongs to Group 8 (medical devices).

Applies in the countries of the European Union, as well as in Norway and Switzerland.

6 Delivery checklist

- 1 Drainage pipe (Ø 16 mm)
- 1 Supply hose (Ø 7 mm) with intake pipe; for the physiological saline solution
- 1 Connecting cable
- 1 Hex. pin driver
- 1 Release pin
- 1 Operating instructions
- 1 Notes on moving the equipment safely

The rotor(s) and associated accessories are included in the delivery in the quantity ordered.

7 Unpacking the centrifuge

Do not lift by the front panel.

Observe the weight of the centrifuge, refer to chapter "Technical specifications".

- Lift up the centrifuge on both sides and take it out of the box with an adequate amount of people to help you.

8 Initial operation

- According to the laboratory instrument standards EN / IEC 61010-2-020 an emergency switch to separate power supply in the event of a failure must be installed in the building electrical system.
This switch has to be placed remote from the centrifuge, preferred outside of the room in which the centrifuge is installed or near by the exit of this room.
- Remove the transportation safety device from the bottom of the housing, see sheet "Transportation safety device"
- **Position the centrifuge in a stable and level manner in a suitable place. During set-up, the required safety margin of 300 mm around the centrifuge is to be kept according to EN / IEC 61010-2-020.**

When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.

- Do not place any object in front of the ventiduct.
Keep a ventilation area of 300 mm around the ventiduct.
- Insert the supplied drainage hose (∅ 16 mm) firmly to the correct coupling connection on the back of the centrifuge.
Connect the free end of the drainage hose to a suitable collection container or drain.


The drainage hose must not be kinked and there must be a continual downward slope to the collecting container. The fault-free functioning of the drain must be checked regularly. If the drainage hose is kinked or blocked, the liquid accumulates in the centrifugal area, pours into the motor area and causes damage to the centrifuge!

There are physiological saline solutions, which are preserved with sodium azide. This can result in explosive azide salts in the drainage pipes. Ascertain with the manufacturer the type of saline solution before you start to put the physiological saline solution needed in to the drain.

- Device without 3-way valve:
Insert the free end of the supplied drainage hose (∅ 7 mm) firmly to the correct coupling connection on the back of the centrifuge. Place the intake pipe of the supply hose into the physiological saline solution container. If need be, pull out the intake pipe (metal part) from the supply hose and connect the supply hose with the help of adapter (plastic part) to the container of the physiological saline solution.
- Device with 3-way valve:
Insert the free ends of the two supplied supply hoses (∅ 7 mm) firmly to the coupling connections "NaCl" and "H₂O" of the 3-way valve. See the chapter on the "3-way valve" about connections of the 3-way valves.
- Place the intake pipe of the supply hose, which is connected to the "NaCl" coupling connection, into the physiological saline solution container. Place the intake pipe of the supply hose, which is connected to the "H₂O" coupling connection, into the distilled water container. If need be, pull out the intake pipe (metal part) from the supply hose and connect the supply hose with the help of adapter (plastic part) to the container. Turn the 3-way valve to the "NaCl" position so that the system is connected to the physiological saline solution; see the chapter on the "3-way valve".
- Check whether the mains voltage tallies with the statement on the type plate.
- Connect the centrifuge with the connection cable to a standard mains socket. For connection ratings refer to Chapter "Technical specification".
- Turn on the mains switch. Switch position "I".
The machine type and program version will be displayed and the LEDs light up. **⚡ OPEN LID ⚡** will be displayed after 8 seconds. The LED **⚡** lights up.
- Open the lid.
The last used centrifuge data will be displayed.
- Remove the transportation safety device in the centrifuge, see sheet "Transportation safety device".

9 3-way valve (only for centrifuges with 3-way valves)

To thoroughly rinse out the system, the 3-way valve can be switched from physiological saline solution (NaCl) to distilled water (H₂O).


Connections:

NaCl: Supply hose (∅ 7 mm) of the physiological saline solution.

H₂O: Supply hose (∅ 7 mm) of the distilled water.

ROTOLAVIT: Connecting hose to the device.

Positions of the 3-way valve:

NaCl: The system is connected to the physiological saline solution.

H₂O: The system is connected to the distilled water.

10 Opening and closing the lid

10.1 Opening the lid

The lid can only be opened when the centrifuge is switched on and the rotor is at rest. If it cannot be opened under these circumstances, see the section on "Emergency release".

- Swing handle rail on the lid upwards. The  LED will now go off.
- Open the lid.


10.2 Closing the lid

Do not bang the lid shut.

- Place the lid and swing handle rail on the lid downward. The  LED will come on.

11 Inserting and removing the rotor

- Hold the rotor down into the recessed grip (A) and place it on the hub (C). Both the driving plates (B) in the rotor must be located in the two grooves (D) of the hub (C). The alignment (E) of the driving plates (B) is indicated on the rotor.
- Removing the rotor: Hold the rotor by the recessed grip (A) and lift it off upwards.

12 Loading the rotor

- The rotors and hangers may only be loaded symmetrically. For authorised combinations see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".
In order to avoid that the physiological salt solution of the non-occupied rotor positions is sprayed into the centrifugation area during centrifugation, all rotor positions should always be loaded with centrifugation containers.
- The maximum filling quantity for the centrifuge containers specified by the manufacturer must not be exceeded.
- The maximum loading weight is stated on every rotor. This weight may not be exceeded. The weight specified for the maximum loading includes the total weight of centrifuging container and content.
- Due to cleaning reasons centrifugation containers made of plastic may only be used for a single application. With multiple use of these centrifugation containers the user is responsible for the suitable cleaning and sterilization. If centrifugation containers have to be used several times then centrifugation containers made of glass should be used.

13 Control and display elements

See figure on page 2.

Fig. 2: Display and control panel

13.1 Control panel symbols

Unbalance display. The unbalance display lights up if the permissible weight difference within the rotor loading has been exceeded.

Rotation display. The rotation display lights up during the centrifugation run as long as the rotor is turning.

Lid display. The lid display lights up if the lid is closed and locked. During the centrifugation run the lid display is extinguished.

Operating errors or any occurring disturbances are symbolised in the display (see Chapter "Faults").

13.2 Control panel buttons


- Button for selecting the individual parameter.
Every time the button is pressed, the subsequent parameter is selected.


- The value of a previously selected parameter is raised or lowered by pressing the arrow buttons. The value can be raised or lowered more quickly if the arrow button is held down.
- Calling up program.
Every time the button is pressed, the subsequent program is displayed (1 - 2 - 3 - 4 - 5 - Spin (S) - Clean Program)

- Starts the washing procedure, the production of the suspension and the cleaning program (Clean Program). The rotation indicator  lights up.
- Adopts entries and changes.

- Starts the centrifugation run for the acceleration of the agglutination reaction and the centrifugation run with the Spin centrifugation program. The rotation indicator  lights up.

- Terminates the washing procedure, the production of the suspension, the centrifugation run for the acceleration of the agglutination reaction, the centrifugation run with the Spin centrifugation program and the cleaning program (Clean Program).
The rotor finishes with a fixed braking level. The LED in the button is lit until the rotor comes to a standstill.
- Pressing the button twice activates the EMERGENCY STOP.
- Exits the parameter input and saves the settings of the global parameters.
- Exits the parameter input without saving the settings of the parameters for the programs 1 to 5.
- During the filling volume calibration terminate the injection of the saline solution.

- Activates the adjustment of the filling volume of saline solution.
The button is only activated when the lid is opened. The LED in the button lights up as the saline solution is injected.
- Refills the system with saline solution during the sequence of the cleaning program (Clean Program).

- Interrupts the washing procedure or the production of the suspension.
- Starts the injection of saline solution whilst adjusting the filling volume.

13.3 Possible settings

There are two adjustment levels:

The global parameters and the program parameters (parameters for the programs 1 to 5 and Spin).

13.3.1 Global parameters

 The adjustment of the global parameters affects the program parameters for the programs 1 to 5.

For the adjustment of the parameters refer to chapter "Global parameters adjustment".

VOLUME ADJUST XX	Correction value for the filling volume of the saline solution. Adjustable from -20 to +20, in single digit steps. 0 = no correction of the filling volume. Preset to 0.
ROTOR: XX PLACE	Rotor type used (12 x or 24 x). Can be adjusted to 12 or 24. Preset to 12.
Shakings XX	Shakings XX:
ShakeT (min) XX	Number of shaking movements during a washing cycle or before the centrifugation run for the agglutination reaction. Adjustable from 0 to 100, in single digit steps. 0 = no shaking is performed. Preset to 15 shaking movements.
	ShakeT (min) XX:
	Duration of the shaking during a washing cycle or before the centrifugation run for the agglutination reaction. Adjustable from 15 minutes to 1 minute, in 1 minute steps.
Susp.Agit	Number of shaking movements after the production of a suspension. Adjustable from 0 to 100 in single digit steps. 0 = no shaking is performed. Preset to 15 shaking movements.
D.SpinDown (s)	Activating or deactivating the adjustment facilities of the DSpinDown program parameters. Through the input of the running time for the centrifugation run to centrifuge the drops down the parameter DSpinDown can be selected. Adjustable from 0 to 20 seconds, in 1 second steps. Preset to 0. 0 = program parameter DSpinDown can not be selected.
SuspensionM	Activating or deactivating the adjustment facilities of the SalSusp/ml program parameters. Adjustable YES or NO. Preset to NO. YES = program parameter SalSusp/ml can be selected. NO = program parameter SalSusp/ml can not be selected.
Agit.SpinM	Activating or deactivating the adjustment facilities of the Agit.Spin program parameters. Adjustable YES or NO. Preset to NO. YES = program parameter Agit.Spin can be selected. NO = program parameter Agit.Spin can not be selected, that means that before the centrifugation run for the agglutination reaction shaking is always performed.
SpinDecantM	Activating or deactivating the adjustment facilities of the SpinDecant program parameters. Adjustable YES or NO. Preset to NO. YES = program parameter SpinDecant can be selected. NO = program parameter SpinDecant can not be selected.
SOUND / BELL OFF	Acoustic signal (refer to chapter "Audible signal"). Adjustable ON1 or OFF. Preset to ON1.
CONTROL: XX	Operating hours of the centrifuge (refer to chapter "Querying operating hours").
VERS 12 °C/* 00, FU/CCI – 1001, FU/CCI –S 01.00	Internal system data. No adjustments possible.

13.3.2 Program parameters (parameters for the programs 1 to 5)

For the adjustment of the parameters refer to chapter "Entering program parameters".

#Cyc/wash	Number of wash cycles. Adjustable from 0 – 9. 0 = no washing cycle is performed.
SalWash/ml	Filling volume of saline solution per fuse (in the washing cycle). Adjustable from 0 - 5 ml, in 0.1ml increments. If the decimal place is greater than zero, it is indicated by a + sign, e.g. 1.7 = 1+. 0.0 = no saline solution is filled.
Wash(rpm)	Rotating speed for the centrifugation run of a washing cycle. A numerical value of 500 RPM to 3500 RPM can be set in increments of 10.
Twash(min), Twash(sec)	Time for the centrifugation run of a washing cycle. Adjustable from 0 - 9 min 59 s, either in 1 second or 1 minute increments. 0:00 = no centrifugation run is performed, and also no decanting.

Decant(rpm)	Rotation speed for the decanting of the excess during a washing cycle. Adjustable from 100 RPM to 1500 RPM, in single digit steps. Preset to 425 RPM.
DSpinDown	Centrifugation run, to centrifuge the droplets down on the tube wall, activate or deactivate. Adjustable ON or OFF. Preset to ON. ON = a centrifugation run is performed. OFF = no centrifugation run is performed. The program parameter DSpinDown can only be selected if the global parameter D.SpinDown (s) does not have the value 0.
SalSusp/ml	Filling volume of saline solution per fuse (for the production of the suspension). Adjustable from 0 - 5 ml, in 0.1ml increments. If the decimal place is greater than zero, it is indicated by a + sign, e.g. 1.7 = 1+. 0.0 = no saline solution is filled. The program parameter SalSusp/ml can only be selected if the global parameter SuspensionM has the value "YES".
Agit.Spin	Switches the shaking that is performed before the centrifugation run for the agglutination reaction ON or OFF. Adjustable ON or OFF. Preset to ON. ON = shaking is performed. OFF = shaking is not performed. The program parameter Agit.Spin can only be selected if the global parameter Agit.SpinM has the value "YES".
Spin(rpm)	Revolutions per minute for the centrifugation cycle for the agglutination reaction. A numerical value of 500 RPM to 3500 RPM can be set in increments of 10.
Tspin(min), Tspin(sec)	Time for the centrifugation cycle for the agglutination reaction. Adjustable from 0 - 9 min 59 s, either in 1 second or 1 minute increments. 0:00 = a centrifugation run is not performed.
SpinDecant	Switches the decanting after the centrifugation run for the agglutination reaction ON or OFF. Adjustable ON or OFF. Preset to OFF. ON = decanting is performed. OFF = decanting is not performed. The program parameter SpinDecant can only be selected if the global parameter SpinDecantM has the value "YES".
PROG = X	Program position on which the program is stored. 5 programs can be stored (program positions 1 - 2 - 3 - 4 - 5).

13.3.3 Program parameters for the centrifugation program Spin

For the adjustment of the parameters refer to chapter "Entering Spin centrifugation program".

Spin(rpm)	Speed for the centrifugation run. A numerical value of 500 RPM to 3500 RPM can be set in increments of 10.
Tspin(min), Tspin(sec)	Time for the centrifugation run. Time for the centrifugation cycle for the agglutination reaction. Adjustable from 0 - 9 min 59 s, either in 1 second or 1 minute increments.
PROG = S	Program position S (Spin) to where the Spin program is saved.

14 Programming

14.1 Global parameters adjustment

The adjustment of the global parameters has an effect on the program parameters for programs 1 to 5.

If after the selection or during input of global parameters a key is not pressed for 16 seconds then the previous values will be displayed on the display again. The adjustments performed up to this moment in time are **automatically saved**. If the adjustments are not to be saved the device must be switched OFF during the 16 seconds.

- Hold the button down for 8 seconds. After 8 seconds **VOLUME ADJUST XX** appears in the display.
- Select the desired parameter with the key , and adjust with the keys .
- To select the parameter **ShakeT (min)**, set the parameter **Shakings** with the key to smaller than "0".
- Press the key or to save the adjustment. After pressing the key ***** ok ***** is briefly displayed as confirmation, after pressing the key no confirmation is given.

14.2 Entering Spin centrifugation program

On the Spin (S) program location a centrifugation program can be entered. In this program only the speed and the running time can be adjusted.

If after the selection or during input of program parameters a key is not pressed for 16 seconds then the previous values will be displayed on the display again. The adjustments performed up to this moment are **not automatically saved**. The input of parameter then has to be executed again.

The input of the program parameters can be terminated at any time by pressing the key. In this case the adjustments are not saved.

- With the or key select the **Spin** program position. The centrifugation data of the selected program position is displayed.
- Select the desired parameter with the key , and adjust with the keys .
- Press the key to save the adjustments to the **S** program position. As confirmation, ***** ok ***** will be displayed for a short period.

The previous data for the program position is overwritten when it is saved.

14.3 Entering program parameters

If after the selection or during input of program parameters a key is not pressed for 16 seconds then the previous values will be displayed on the display again. The adjustments performed up to this moment are **not automatically saved**. The input of parameter then has to be executed again.

The input of the program parameters can be terminated at any time by pressing the key. In this case the adjustments are not saved.

- For adjustment of the global parameters refer to chapter "Globale Parameter einstellen".
- Select the desired parameter with the key , and adjust with the keys .
- Select the parameter **PROG** with the key , and set the desired programme position with the keys . The LED in the key lights up.
- Press the key in order to store the adjustments to the desired programme position. As confirmation, ***** ok ***** will be displayed for a short period.


The previous data for the program position is overwritten when it is saved.

14.4 Programme recall

- Select the desired program position by pressing the ▲ or ▼ button.
The centrifugation data of the selected program position is displayed.
With selection of the cleaning program, **CLEAN PROGRAM** is displayed.
- The parameters can be checked by pressing the ◀ button.
With selection of the cleaning program (CLEAN PROGRAM) no parameters can be verified.
To quit the parameter display, press the ⓐSTOP button or don't press any button for 16 seconds.

15 Functions

15.1 Description

* the function can be selected or deselected

Display	Action
FILL	Fill the tubes with saline solution.
SPIN	Centrifugation run
DCNT	Decanting the excess.
DOWN	Centrifugation run to centrifuge the remaining droplets on the side of the tube to the bottom of the tube.
AGIT	Shaking
SUSP	Production of a suspension


During a washing process, which consists of several wash cycles, the physiological saline solution is added successively, the erythrocytes are sedimented by the centrifugation process, and finally the physiological saline solution is decanted again.

After manually adding the anti-human globulin serum, a further centrifugation speeds up the agglutination test.

The following functions can be selected or deselected:

- Centrifugation run to centrifuge the remaining droplets on the side of the tube to the bottom of the tube (DOWN).
- The shaking of the tube at the end of a washing cycle (AGIT)
- The production of a suspension (SUSP) and/or the shaking of the tube after the production of a suspension (AGIT).
- The shaking of the tubes before a centrifugation run that serves to accelerate the agglutination reaction (AGIT).
- The decanting of the excess after the centrifugation run that serves to accelerate the agglutination reaction (DCNT).

15.1.1 Washing process

* the function can be selected or deselected

A washing process consists of several successive automatic current wash cycles. The number of wash cycles is variable.

A washing cycle consists of 3 to 5 steps:

1. FILL

Filling in the physiological saline solution.

The tubes are filled with physiological saline solution at a rotation speed of 1100 RPM. The number of revolutions per minute is not adjustable. The physiological saline solution is injected directly into the tubes, in order to receive a good resuspension of the cells.

2. SPIN

Sedimentation.

The erythrocytes are sedimented by a selectable rotation speed. The time only begins to count after the set rotation speed is reached. After the time has elapsed, a quick braking follows to prevent a resuspension of the pellet.

3. DCNT

Decanting.

The excess is decanted at a selectable rotation speed. During decanting, the rotor turns in the opposite direction of rotation.

4. DOWN

Centrifugation run to centrifuge the remaining droplets on the side of the tube to the bottom of the tube. The speed is 2000 RPM and can not be adjusted.

This centrifugation run can be selected and deselected.

Selecting centrifugation run:

In the global parameters **D.SpinDown (s)** adjust the running time and set the program parameter **DSpinDown** to "ON".

Deselect centrifugation run:

Either set the global parameter **D.SpinDown (s)** to "0" or in the global parameter **D.SpinDown (s)** adjust the running time and set the program parameter **DSpinDown** to "OFF".

5. AGIT

Shaking.

The pellet is loosened again for the subsequent wash cycle by a fast, brief movement of the rotor.

In order not to loosen the pellet, no shaking occurs during a washing process with only one wash cycle or during the last wash cycle of a washing process.


Select shaking:

Either in the global parameters **Shakings** select the number of shaking movements, or in the global parameters **ShakeT (min)** adjust the duration of the shaking.

Deselect shaking:

Set the global parameter **Shakings** to "0".

15.1.2 Suspension

* the function can be selected or deselected

A suspension can be produced with or without a previous washing process. At a speed of 1100 RPM the tubes are filled with a selectable amount of saline solution. The production of a suspension can be selected or deselected.

Select suspension:

Set the global parameter **SuspensionM.** to "YES", and in the program parameter **SalSusp/ml** adjust the filling volume.

Deselect suspension:

Either set the global parameter **SuspensionM.** to "NO" or set the program parameter **SalSusp/ml** to "0".

The shaking after the production of the suspension can be selected or deselected.


Select shaking:

In the global parameter **Susp.Agit.** adjust the number of shaking movements.

Deselect shaking:

Set the global parameter **Susp.Agit.** to "0".

15.1.3 Centrifugation run for acceleration of the agglutination reaction

* the function can be selected or deselected

After manually adding the anti-human globulin serum, a further centrifugation speeds up the agglutination test.

The shaking before the centrifugation run can be selected or deselected.

Select shaking:

Either set the global parameter **Agit.SpinM.** to "NO" or set the global parameter **Agit.SpinM.** to "YES" and the program parameter **Agit.Spin** to "ON".

Deselect shaking:

Set the global parameter **Agit.SpinM.** to "YES" and the program parameter **Agit.Spin** to "OFF".

The decanting of the excess after the centrifugation run can be selected or deselected.


Select decanting:

Set the global parameter **SpinDecantM.** to "YES" and the program parameter **SpinDecant** to "ON".

Deselect decanting:

Either set the global parameter **SpinDecantM.** to "NO" or set the global parameter **SpinDecantM.** to "YES" and the program parameter **SpinDecant** to "OFF".

15.1.4 Centrifugation program Spin


On the Spin (S) program position a centrifugation program can be saved. In this program only the speed and the running time can be adjusted.

15.1.5 Cleaning program (Clean Program)

To avoid the formation of salt crystals the system must be rinsed through daily after use with distilled water. Before the next use the system must be refilled again with saline solution.

This is performed with the help of the cleaning program (Clean Program).

In this program no adjustments can be made.

The program consists of the following steps:

1. FILL
At a speed of 1100 RPM 10 ml of distilled water is filled into each tube (120 ml with the 12 x rotor and 240 ml with the 24 x rotor).
2. SPIN
Centrifugation run with a running time of 10 seconds and a speed of 1500 RPM.
3. DCNT
Decanting.
The distilled water is decanted at a speed of 600 RPM. For decanting the rotor turns in the reverse direction.
4. AGIT
Shaking is performed 5 times.
5. DCNT
Decanting.
The distilled water is decanted at a speed of 600 RPM. For decanting the rotor turns in the reverse direction.
6. Refill process
The tubes are refilled with saline solution again. 60 ml of saline solution is pumped through the tubes.

A container must be held in front of the injection pipe to catch the distilled water.

15.2 Combining functions

Through the careful operation of the centrifuge the user determines the precision of the test results.

 If the permissible weight difference within the rotor loading has been exceeded, the drive shuts down during the start-up, the unbalance display lights up, and **IMBALANCE** is displayed.
A function can be terminated at any time by pressing the **STOP** key.
If **OPEN LID** or **WASH COMPLETE** or **SUSP. COMPLETE** or **man. SUSP. DONE** are displayed, then the further operation of the centrifuge is only possible after a single opening of the lid.

 Be sure that the system is filled with physiological saline solution before working with test material.
In order not to damage the rotor a washing procedure may only be started with the **START WASH** key or a centrifugation run may only be started with the **SPIN** key only after the rotor has come to a standstill.

- Turn on the mains switch. Switch position I.

15.2.1 Performing washing and suspension procedure

- Load the rotor and close the centrifuge lid.
- Set the parameters or call up a programme (see the "Programming" chapter).
- Press the **START WASH** button. The rotation indicator  lights up as long as the rotor is turning.

The washing procedure and the production of a suspension are executed:

1. Filling in the physiological saline solution.

The LED in the **SALINE** key illuminates whilst the saline solution is being pumped in and **FILL** is displayed.

PROG	CYC	ml	RPM	t/ min:s
1	3	3	1100	FILL

If the program parameter **SalWash/ml** has the value "0.0", then saline solution is not filled.

2. Sedimentation.

The time only begins to count after the set rotation speed is reached.

Briefly **SPIN** is displayed.

PROG	CYC	ml	RPM	t/ min:s
1	3	3	3500	SPIN

If the program parameter **Twash(min)** and **Twash(sec)** have the value "0:00", then a centrifugation run is not performed.

3. Decanting.

DCNT is displayed.

PROG	CYC	ml	RPM	t/ min:s
1	3	3	425	DCNT

If the program parameters **Twash(min)** and **Twash(sec)** have the value "0:00", then decanting is not performed.

4. Centrifugation run to centrifuge the remaining droplets on the side of the tube to the bottom of the tube.

DOWN is displayed.

PROG	CYC	ml	RPM	t/ min:s
1	3	3	2000	DOWN

The centrifugation run is only performed if the function is activated (global parameter **D.SpinDown (s)** > "0" and program parameter **D.SpinDown** = "ON").

5. Shaking.
AGIT is displayed.


The steps 1 to 5 are performed until all washing cycles have been performed.

Shaking is only performed if the function is activated. (Global parameter **Shakings** > "0" or global parameter **ShakeT (min)** > "0").

In order not to loosen the pellet, no shaking occurs during a washing process with only one wash cycle or during the last wash cycle of a washing process.

6. Suspension.
SUSP is displayed.

The suspension is only produced if the function is activated (global parameter **SuspensionM.** = "YES" and program parameter **SalSusp/ml** > "0").

7. Shaking.
AGIT is displayed.

Shaking is only performed if the function is activated. (Global parameter **Susp.Agit.** > "0").

- After completion of the washing procedure or the suspension and with termination of the washing procedure or the suspension, by pressing the **STOP** key the run-out is performed at a preset brake level. The brake level is displayed.

Display after production of a suspension when rotor is stationary:

Display after a washing cycle when rotor is stationary:

After a washing procedure a centrifugation run for acceleration of the agglutination reaction must be performed. For this the **SPIN** key must be pressed.

The **START WASH** key is locked. By pressing the **START WASH** key the following display appears:

If a centrifugation run is not desired the **STOP** key must be pressed twice.

After the first pressing of the **STOP** key the following display appears:

After the second pressing of the **STOP** key the locking of the **START WASH** key is released and the centrifugation data of the selected program position is displayed.

During the run the remaining number of washing cycles, filling volume of the saline solution, the speed of the rotor and the remaining time of the centrifugation run are displayed.

15.2.2 Performing centrifugation run for acceleration of the agglutination reaction

- Put the anti-human globulin serum in the tubes and close the centrifuge lid.
- Press the **SPIN** button. The rotation indicator  lights up as long as the rotor is turning.

1. Shaking.

Spin and **AGIT** are displayed.


Shaking is only performed if the function is activated (global parameter **Agit.SpinM.** = "NO", or global parameter **Agit.SpinM.** = "YES" and program parameter **Agit.Spin** = "ON").

2. Sedimentation.


The time only begins to count after the set rotation speed is reached.

Spin and **SPIN** are briefly displayed.


3. Decanting.

Spin and **DCNT** are displayed.

Decanting is only performed if the function is activated. (Global parameter **SpinDecantM.** = "YES" and program parameter **SpinDecant** = "ON").

- After the time has elapsed or after the centrifugation cycle has been interrupted by pressing the **STOP** button, it ends with a fixed braking level. The braking level is displayed.
Display when rotor is stationary:

The rotation speed of the rotor and the time remaining is displayed during the centrifugation cycle.

15.2.3 Production of a suspension without previous washing procedure

It is possible to produce a suspension without previously performing a washing process.

- Load the rotor and close the centrifuge lid.
- Set the parameters or call up a programme (see the "Programming" chapter).


The adjustment of the following values is absolutely necessary.

Program parameter **#Cyc/wash** = "0", **Twash(min)**, **Twash(sec)** = "0:00", **SalWash/ml** = "0.0";
Global parameter **D.SpinDown (s)** = "0" or program parameter **DSpinDown** = "OFF";
Global parameter **SuspensionM.** = "YES" and program parameter **SalSusp/ml** > "0".

- Press the **START WASH** button. The rotation indicator  lights up as long as the rotor is turning.

1. Suspension.


SUSP is displayed.

The suspension is only produced if the function is activated (global parameter **SuspensionM.** = "YES" and program parameter **SalSusp/ml** > "0").

2. Shaking.

AGIT is displayed.

Shaking is only performed if the function is activated. (Global parameter **Susp.Agit.** > "0").

- After production of the suspension termination or through pressing the **STOP** key the run-out is performed with a fixed braking level. The braking level is displayed.
Display when rotor is stationary:

15.2.4 Checking the samples during the washing procedure or production of the suspension

The washing procedure and the production of the suspension can be interrupted for checking the samples.

- Press the **CHECK** button. The LED in the **CHECK** key illuminates.
Display when rotor is stationary:

- Open the lid and check the tubules.
- Close the lid again and press the **START WASH** key to continue with the washing procedure or the production of the suspension.

15.2.5 Performing centrifugation run with Spin centrifugation program

- Load the rotor and close the centrifuge lid.
- With the **▲** or **▼** key select the **Spin** program position.
The centrifugation data of the **Spin** program position is displayed.
- Press the **SPIN** button. The rotation indicator **⦿** lights up as long as the rotor is turning.
The time only begins to count after the set rotation speed is reached.
Spin and briefly **SPIN** are displayed.

- After the time has elapsed or after the centrifugation cycle has been interrupted by pressing the **STOP** button, it ends with a fixed braking level. The braking level is displayed.
Display when rotor is stationary:


The rotation speed of the rotor and the time remaining is displayed during the centrifugation cycle.

15.2.6 System with rinsing cleaning program (Clean Program)


Rinsing out the system:

- Device without 3-way valve:
Take off the supply hose from the physiological saline solution container and connect it to a container which is filled with distilled water.
- Device with 3-way valve:
Turn the 3-way valve to the "H₂O" position so that the system is connected to the distilled water; see the chapter on the "3-way valve".
- Load every second rotor position with an empty tube and close the centrifuge lid.


- With the **▲** or **▼** key select the **Clean Program** program position.
- Press the **START WASH** key. The rotation display **⦿** illuminates as long as the rotor is turning.

1. Filling the distilled water.
The LED in the **SALINE** key illuminates whilst the distilled water is injected in.
Clean and **FILL** is displayed.

2. Centrifugation run.
Clean and **SPIN** is displayed.

- 3. Decanting.
Clean and **DCNT** is displayed.

- 4. Shaking.
Clean and **AGIT** is displayed.

- 5. Decanting.
Clean and **DCNT** is displayed.

- After completion or after termination of the cleaning program by pressing the **STOP** key the run-out is performed with the fixed braking level. The braking level is displayed.
Display with stationary rotor:

- Open the lid.
Display after completion of the cleaning program:

If the cleaning program is terminated with the **STOP** key before the distilled water was decanted the following display appears:

 The cleaning program must be started again.

- Switch the unit OFF.

Filling the system:

- Device without 3-way valve:
Take off the supply hose from the distilled water container and connect it again to the physiological saline solution container.
- Device with 3-way valve:
Turn the 3-way valve to the "NaCl" position so that the system is connected to the physiological saline solution; see the chapter on the "3-way valve".
- Switch the unit ON

Display:

 If **OPEN LID** is displayed open the lid.

- Hold a container with a capacity of at least 80 ml in front of the injection pipe in order to catch the distilled water. Then press and hold the **SALINE** key for 3 seconds.
If a different key than the **SALINE** key is pressed the following display appears:

The system is filled with saline solution.

Display:

After the filling the following display appears briefly:

Then the centrifugation data of program position 1 is displayed.

15.3 Examples for the function selection

 The adjustment of the bold marked values is absolutely required for the sequence of the different functions. The adjustment of the global parameters affects the program parameters for the programs 1 to 5. In the following there are respectively 2 possibilities displayed for adjustment of the different functions.

15.3.1 Washing with subsequent centrifugation with shaking

Basic program for cross samples

Adjustments			
Global Parameters		Program Parameters	
Shakings	15	#Cyc/wash	3
Susp.Agit.	15	SalWash/ml	3.0
		Wash(rpm)	3500
		Twash(min), Twash(sec)	0:35
		Decant(rpm)	425
D.SpinDown (s)	3	DSpinDown	OFF
SuspensionM.	YES	SalSusp/ml	0.0
Agit.SpinM.	YES	Agit.Spin	ON
		Spin(rpm)	1100
		Tspin(min), Tspin(sec)	0:20
SpinDecantM.	YES	SpinDecant	OFF
		PROG =	1

Adjustments			
Global Parameters		Program Parameters	
Shakings	15	#Cyc/wash	3
Susp.Agit.	15	SalWash/ml	3.0
		Wash(rpm)	3500
		Twash(min), Twash(sec)	0:35
		Decant(rpm)	425
D.SpinDown (s)	0		
SuspensionM.	NO		
Agit.SpinM.	NO		
		Spin(rpm)	1100
		Tspin(min), Tspin(sec)	0:20
SpinDecantM.	NO		
		PROG =	1

Press the **(START WASH)** button. Press the **(SPIN)** for centrifugation run.

15.3.2 Washing with subsequent suspension without shaking the suspension

Adjustments			
Global Parameters		Program Parameters	
Shakings	15	#Cyc/wash	3
Susp.Agit.	0	SalWash/ml	3.0
		Wash(rpm)	3500
		Twash(min), Twash(sec)	0:35
		Decant(rpm)	425
D.SpinDown (s)	3	DSpinDown	OFF
SuspensionM.	YES	SalSusp/ml	4.0
Agit.SpinM.	YES	Agit.Spin	ON
		Spin(rpm)	1100
		Tspin(min), Tspin(sec)	0:20
SpinDecantM.	YES	SpinDecant	OFF
		PROG =	2

Adjustments			
Global Parameters		Program Parameters	
Shakings	15	#Cyc/wash	3
Susp.Agit.	0	SalWash/ml	3.0
		Wash(rpm)	3500
		Twash(min), Twash(sec)	0:35
		Decant(rpm)	425
D.SpinDown (s)	0		
SuspensionM.	YES	SalSusp/ml	4.0
Agit.SpinM.	NO		
		Spin(rpm)	1100
		Tspin(min), Tspin(sec)	0:20
SpinDecantM.	NO		
		PROG =	2

Press the **(START WASH)** button.

15.3.3 Washing with subsequent suspension with shaking the suspension

Adjustments			
Global Parameters		Program Parameters	
Shakings	15	#Cyc/wash	3
Susp.Agit.	15	SalWash/ml	3.0
		Wash(rpm)	3500
		Twash(min), Twash(sec)	0:35
		Decant(rpm)	425
D.SpinDown (s)	3	DSpinDown	OFF
SuspensionM.	YES	SalSusp/ml	4.0
Agit.SpinM.	YES	Agit.Spin	ON
		Spin(rpm)	1100
		Tspin(min), Tspin(sec)	0:20
SpinDecantM.	YES	SpinDecant	OFF
		PROG =	2

Adjustments			
Global Parameters		Program Parameters	
Shakings	15	#Cyc/wash	3
Susp.Agit.	15	SalWash/ml	3.0
		Wash(rpm)	3500
		Twash(min), Twash(sec)	0:35
		Decant(rpm)	425
D.SpinDown (s)	0	SalSusp/ml	4.0
SuspensionM.	YES	Spin(rpm)	1100
Agit.SpinM.	NO	Tspin(min), Tspin(sec)	0:20
		SpinDecantM.	NO
		PROG =	2

Press the **(START WASH)** button.

15.3.4 Only suspension without shaking

Adjustments			
Global Parameters		Program Parameters	
Shakings	15	#Cyc/wash	0
Susp.Agit.	0	SalWash/ml	0.0
		Wash(rpm)	3500
		Twash(min), Twash(sec)	0:00
		Decant(rpm)	425
D.SpinDown (s)	3	DSpinDown	OFF
SuspensionM.	YES	SalSusp/ml	4.0
Agit.SpinM.	YES	Agit.Spin	ON
		Spin(rpm)	1100
		Tspin(min), Tspin(sec)	0:20
SpinDecantM.	YES	SpinDecant	OFF
		PROG =	3

Adjustments			
Global Parameters		Program Parameters	
Shakings	15	#Cyc/wash	0
Susp.Agit.	0	SalWash/ml	0.0
		Wash(rpm)	3500
		Twash(min), Twash(sec)	0:00
		Decant(rpm)	425
D.SpinDown (s)	0	SalSusp/ml	4.0
SuspensionM.	YES	Spin(rpm)	1100
Agit.SpinM.	NO	Tspin(min), Tspin(sec)	0:20
		SpinDecantM.	NO
		PROG =	3

Press the **(START WASH)** button.

15.3.5 Only suspension with shaking

Adjustments			
Global Parameters		Program Parameters	
Shakings	15	#Cyc/wash	0
Susp.Agit.	15	SalWash/ml	0.0
		Wash(rpm)	3500
		Twash(min), Twash(sec)	0:00
		Decant(rpm)	425
D.SpinDown (s)	3	DSpinDown	OFF
SuspensionM.	YES	SalSusp/ml	4.0
Agit.SpinM.	YES	Agit.Spin	ON
		Spin(rpm)	1100
		Tspin(min), Tspin(sec)	0:20
SpinDecantM.	YES	SpinDecant	OFF
		PROG =	3

Adjustments			
Global Parameters		Program Parameters	
Shakings	15	#Cyc/wash	0
Susp.Agit.	15	SalWash/ml	0.0
		Wash(rpm)	3500
		Twash(min), Twash(sec)	0:00
		Decant(rpm)	425
D.SpinDown (s)	0	SalSusp/ml	4.0
SuspensionM.	YES	Spin(rpm)	1100
Agit.SpinM.	NO	Tspin(min), Tspin(sec)	0:20
		SpinDecantM.	NO
		PROG =	3

Press the **(START WASH)** button.

15.3.6 Centrifugation run with decanting

Adjustments			
Global Parameters		Program Parameters	
Shakings	15	#Cyc/wash	3
Susp.Agit.	15	SalWash/ml	3.0
		Wash(rpm)	3500
		Twash(min), Twash(sec)	0:35
		Decant(rpm)	425
D.SpinDown (s)	3	DSpinDown	OFF
SuspensionM.	YES	SalSusp/ml	0.0
Agit.SpinM.	YES	Agit.Spin	OFF
		Spin(rpm)	1100
		Tspin(min), Tspin(sec)	0:20
SpinDecantM.	YES	SpinDecant	ON
		PROG =	4

Adjustments			
Global Parameters		Program Parameters	
Shakings	15	#Cyc/wash	3
Susp.Agit.	15	SalWash/ml	3.0
		Wash(rpm)	3500
		Twash(min), Twash(sec)	0:35
		Decant(rpm)	425
D.SpinDown (s)	0		
SuspensionM.	NO		
Agit.SpinM.	NO		
		Spin(rpm)	1100
		Tspin(min), Tspin(sec)	0:20
SpinDecantM.	YES	SpinDecant	ON
		PROG =	4

Press the **[SPIN]** button.

15.3.7 Centrifugation run with SPIN (S) program

Adjustments			
Global Parameters		Program Parameters	
No adjustments required		Spin(rpm)	3500
		Tspin(min), Tspin(sec)	0:20
		PROG =	S

Press the **[SPIN]** button.

16 Emergency stop

- Press the **[STOP]** button twice.

In the event of an emergency stop, it finishes with a braking level of 9 (shortest run-down time). Braking level 9 is displayed.

17 Adjusting the filling volume

It is recommended to check the filling volume of the physiological saline solution every day and, if required, adjust it.

 In order to prevent the physiological saline solution from spraying or spilling, it is imperative to place a vessel in front of the injection pipe before adjusting the filling volume.

 If no key is pressed for 16 seconds or the **[STOP]** key is pressed then the previous values are displayed again on the display. The adjustments performed up to this moment in time are **automatically saved**. If the adjustments are not to be saved the device must be switched OFF during the 16 seconds.

It is only possible to adjust the filling volume when the rotor has stopped and the lid is opened.

- Hold the **[SALINE]** button down for 4 seconds.
After 4 seconds **CALIBRATE xx.xml** appears in the display.

 The filling volume displayed corresponds to the filling volume of a fully loaded rotor (12-fold or 24-fold). The calculation is performed with the parameter set in **Sal/Wash/ml**.

- Hold a vessel in front of the injection pipe and press the **[CHECK]** button. The vessel is filled with the displayed filling volume of physiological saline solution.
- If the filling volume does not correspond with the display, hold the **[◀]** button down for 8 seconds.
After 8 seconds **VOLUME ADJUST XX** appears in the display.

- Set the correction factor using the ▲▼ buttons. A value of 20 to –20 can be set. An adjustment of 1 corresponds to an adjustment in filling volume of 0,5% of the measured filling volume.

$$\text{Correction} = \text{deviation (ml)} \times \frac{200}{\text{measured filling volume (ml)}}$$

Deviation (ml) = desired filling volume (ml) – measured filling volume (ml)

- Press the **START WASH** button to store the setting.
*** ok *** is briefly displayed as a confirmation.

18 Entering the rotor type

It is necessary to enter the rotor type used (12-fold or 24-fold) in order to calculate the filling volume. It is only possible to enter the rotor type when the rotor has stopped.

 If no key is pressed for 16 seconds or the **STOP** key is pressed then the previous values are displayed again on the display. The adjustments performed up to this moment in time are **automatically saved**. If the adjustments are not to be saved the device must be switched OFF during the 16 seconds.

- Hold the **ENTER** button down for 8 seconds.
After 8 seconds **VOLUME ADJUST XX** appears in the display.
- Press the **ENTER** button once more.
ROTOR: XX PLACE is displayed.
- Set the rotor type (12-fold or 24-fold) with the ▲ or ▼ button
- Press the **START WASH** button to store the setting.
*** ok *** is briefly displayed as a confirmation.

19 Audible signal

The audible signal sounds:

- in 2 second intervals if a fault occurs.
- in 30 second intervals after the centrifugation cycle has ended and the rotor has stopped.

The audible signal can be stopped by opening the lid or by pressing any button.

The signal after completion of the centrifugation run can be activated or deactivated in the following manner, if the rotor is at standstill:

 If no key is pressed for 16 seconds or the **STOP** key is pressed then the previous values are displayed again on the display. The adjustments performed up to this moment in time are **automatically saved**. If the adjustments are not to be saved the device must be switched OFF during the 16 seconds.

- Hold the **ENTER** button down for 8 seconds.
After 8 seconds **VOLUME ADJUST XX** appears in the display.
- Repeat pressing the **ENTER** button until **SOUND / BELL** is displayed.
- Set **OFF** (off) or **ON1** (on) using the ▲ or ▼ button.
- Press the **START WASH** button to store the setting.
*** ok *** is briefly displayed as a confirmation.

20 Querying operating hours

It is only possible to query the operating hours when the rotor has stopped.

 If no key is pressed for 16 seconds or the **STOP** key is pressed then the previous values are displayed again on the display. The adjustments performed up to this moment in time are **automatically saved**. If the adjustments are not to be saved the device must be switched OFF during the 16 seconds.

- Hold the **ENTER** button down for 8 seconds.
After 8 seconds **VOLUME ADJUST XX** appears in the display.
- Repeat pressing the **ENTER** button until the operating hours (**CONTROL:**) of the centrifuge are displayed.
- To quit the operating hours display, press the **STOP** button or don't press any button for 16 seconds.

21 Relative centrifugal force (RCF)

The relative centrifugal force (RCF) is given as a multiple of the acceleration of gravity (g). It is a unit-free value and serves to compare the separation and sedimentation performance.

These values are calculated using the formula below:

$$\text{RCF} = \left(\frac{\text{RPM}}{1000} \right)^2 \times r \times 1,118 \quad \Rightarrow \quad \text{RPM} = \sqrt{\frac{\text{RCF}}{r \times 1,118}} \times 1000$$

RCF = relative centrifugal force

RPM = rotational speed (revolutions per minute)

r = centrifugal radius in mm = distance from the centre of the turning axis to the bottom of the centrifuge.
For more on the centrifugal radius see the chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

The relative centrifugal force (RCF) stands in relation to the revolutions per minute and the centrifugal radius.

22 Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm³

When centrifuging with maximum revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm³.

The speed must be reduced for materials or mixtures of materials with a higher density.

The permissible speed can be calculated using the following formula:

$$\text{Reduced speed } (n_{\text{red}}) = \sqrt{\frac{1.2}{\text{Greater density [kg/dm}^3]}} \times \text{maximum speed [RPM]}$$

e.g.: maximum speed RPM 4000, density 1.6 kg/dm³

$$n_{\text{red}} = \sqrt{\frac{1.2 \text{ kg/dm}^3}{1.6 \text{ kg/dm}^3}} \times 4000 \text{ RPM} = 3464 \text{ RPM}$$

In the exceptional case that the maximum loading indicated on the hanger is exceeded, the speed must also be reduced.

The permissible speed can be calculated using the following formula:

$$\text{Reduced speed } (n_{\text{red}}) = \sqrt{\frac{\text{maximum load [g]}}{\text{actual load [g]}}} \times \text{maximum speed [RPM]}$$

e.g.: maximum speed RPM 4000, maximum load 300 g, actual load 350 g

$$n_{\text{red}} = \sqrt{\frac{300 \text{ g}}{350 \text{ g}}} \times 4000 \text{ RPM} = 3703 \text{ RPM}$$

If in doubt you should obtain clarification from the manufacturer.

23 Emergency release

The lid cannot be opened during power failure. An emergency release has to be executed by hand.

For emergency release disconnect the centrifuge from the mains.

Open the lid only during rotor standstill.

Only the plastic release pin provided may be used for emergency release.

See figure on page 2.

- Switch off the mains switch (switch position "0").
- Look through the window in the lid to be sure that the rotor has come to a standstill.
- Insert the release pin (see scope of supply) horizontally into the hole (Fig. 1, A). Push the unlocking pin in until the handle can be lifted when the pin is pressed down.
- Open the lid.

24 Maintenance and servicing

Pull the mains plug before cleaning.

For safety reasons you must wear gloves and a respiratory mask when cleaning the equipment used for blood processing.

Before any other cleaning or decontamination process other than that recommended by the manufacturer is applied, the user has to check with the manufacturer that the planned process does not damage the device.

- Centrifuges, rotors and accessories must not be cleaned in rinsing machines.
- They may only be cleaned by hand and disinfected with liquids.
- The water temperature must be between 20 – 25°C.
- Only detergents/disinfectants may be used which:
 - have a pH between 5 - 8
 - do not contain caustic alkalis, peroxides, chlorine compounds, acids and alkaline solutions
- In order to prevent appearances of corrosion through cleaning agents or disinfectants, the application guide from the manufacturer of the cleaning agent or disinfectant are absolutely to be heeded.
- Certain preservation agents, in azide-free saline solutions can cause long term damage effects to the plastic parts in the device. Regular cleaning prevents salt deposits and lengthens the life span of these parts.

24.1 Centrifuge

- The following must be carried out daily:
 - Check the hoses and their attachments. The hoses must not be cracked or blocked and must be securely attached. The used saline solution must be allowed to flow off freely through the drainage pipe.
 - The centrifugal area must be clean and free from dried up salt crystals and other deposits. Clean the centrifugation area, the drain gutter and the covering ring with a moist cloth or sponge. The drain gutter and the covering ring can be removed from the centrifugation area for cleaning, refer to Chapter "Removal of the drain gutter".
 - Check the filling volume of the saline solution (see the "Adjusting the filling volume" chapter).
 - The system must be rinsed through with distilled water to avoid formation of salt crystals, refer to Chapter "System with rinsing cleaning program (Clean Program)".
- It is important that the hoses are kept clean and free from dried up salt crystals and other deposits.
- The system must be regularly cleaned, refer to Chapter "Clean system with cleaning solution". It is recommended to clean the appliance at least once a week.
- Clean the centrifuge housing and the centrifuging chamber regularly, using soap or a mild detergent and a damp cloth if required. This serves as hygiene protection and prevents corrosion caused by impurities.
- Ingredients of suitable detergents:
 - soap, anionic tensides, non-ionic tensides.
- After using detergents, remove the detergent residue by wiping with a damp cloth.
- The surfaces must be dried immediately after cleaning.
- Lightly rub the rubber seal of the centrifuge chamber with talcum powder or a rubber care product after each cleaning.
- Surface disinfection:
 - If infectious materials penetrates into the centrifugal chamber this is to be disinfected immediately.
 - Ingredients of suitable disinfectants:
 - ethanol, n-propanol, isopropyl alcohol, glutardialdehyde, quaternary ammonium compounds.
 - After using disinfectants, remove the disinfectant residue by wiping with a damp cloth.
 - The surfaces must be dried immediately after disinfecting.
- Removal of radioactive contaminants:
 - The agent must be specifically labelled as being an agent for removing radioactive contaminants.
 - Ingredients of suitable agents for removing radioactive contaminants:
 - anionic tensides, non-ionic tensides, polyhydrated ethanol.
 - After removing the radioactive contaminants, remove the agent residue by wiping with a damp cloth.
 - The surfaces must be dried directly after removing the radioactive contaminants.
- The centrifuging chamber is to be checked for damage once a year.

If damage is found which is relevant to safety, the centrifuge may no longer be put into operation. In this case, notify Customer Service.

24.2 Rotor

- It is important that the rotor is kept clean and free from dried up salt crystals and other deposits.
- Either soak the rotor in warm, distilled water or pour the water directly down into the rotor for a few minutes. The water must flow out of all the injection nozzles.
- If the injection nozzles are blocked, insert the included plastic pin into the injection nozzles and carefully slide it in and out until the nozzles become clear again.
- In order to prevent corrosion and material changes, rotors and accessories must be cleaned regularly with soap or a mild detergent and a damp cloth. Cleaning is recommended at least once a week. Contaminants must be removed immediately.
- Ingredients of suitable detergents: soap, anionic tensides, non-ionic tensides.
- After using detergents, remove detergent residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotors and accessories must be dried directly after cleaning.
- Disinfection:
 - If infectious material should get on the rotors or accessories, they must be appropriately disinfected.
 - Ingredients of suitable disinfectants: glutaraldehyde, propanol, ethyl hexanol, anionic tensides, corrosion inhibitors.
 - After using disinfectants, remove disinfectant residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
 - The rotors and accessories must be dried directly after disinfection.
- Removal of radioactive contaminants:
 - The agent must be specifically labelled as being an agent for the removal of radioactive contaminants.
 - Ingredients of suitable agents for removing radioactive contaminants: anionic tensides, non-ionic tensides, polyhydrated ethanol.
 - After removing the radioactive contaminants, remove agent residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
 - The rotors and accessories must be dried directly after removing the radioactive contaminants.
- The rotor is to be checked for corrosion damage every month.

If there are signs of wearing or corrosion, the rotors and accessories must no longer be used.

24.3 Autoclaving

The system must be regularly cleaned, refer to Chapter "Cleaning and rinsing system".

Centrifugation containers made of glass can be autoclaved at 121°C (20 min).

24.4 Glass breakage

- In the case of glass breakage, the pieces of glass and leaked centrifugal liquid must be carefully removed from the centrifugal area and tube holders.
- If the leakage involves infectious materials, the area must be disinfected immediately.

24.5 Clean system with cleaning solution

24.5.1 Device without 3-way valve

- Make about 400 ml of 0.5% sodium hypochlorite cleaning solution.
- Take off the supply hose from the physiological saline solution container and connect it to the cleaning solution container.
- Rinse and fill the system with cleaning solution using the clean program; see the chapter on "Rinsing the system with the clean program".

Also carry out the filling of the system with cleaning solution.

- Take off the supply hose from the cleaning solution container and connect it to a container which is filled with distilled water.
- Rinse the system with distilled water using the clean program; see the paragraph on "Rinsing the system" in the chapter on "Rinsing the system with the clean program".
- Take off the supply hose from the distilled water container and connect it again to the physiological saline solution container.
- Fill the system with physiological saline solution using the clean program; see the paragraph on "Filling the system" in the chapter on "Rinsing the system with the clean program".
- Dry the centrifugal area and rotor.
- Adjust the filling volume (see the "Adjusting the filling volume" chapter).

24.5.2 Device with 3-way valve:

- Make about 400 ml of 0.5% sodium hypochlorite cleaning solution.
- Take off the supply hose from the distilled water container and connect it to the cleaning solution container.
- Turn the 3-way valve to the "H₂O" position so that the system is connected to the cleaning solution; see the chapter on the "3-way valve".
- Rinse and fill the system with cleaning solution using the clean program; see the chapter on "Rinsing the system with the clean program".

Also carry out the filling of the system with cleaning solution.

- Take off the supply hose from the cleaning solution container and connect it again to the distilled water container.
- Rinse the system with distilled water using the clean program; see the paragraph on "Rinsing the system" in the chapter on "Rinsing the system with the clean program".
- Turn the 3-way valve to the "NaCl" position so that the system is connected to the physiological saline solution; see the chapter on the "3-way valve".
- Fill the system with physiological saline solution using the clean program; see the paragraph on "Filling the system" in the chapter on "Rinsing the system with the clean program".
- Dry the centrifugal area and rotor.
- Adjust the filling volume (see the "Adjusting the filling volume" chapter).

24.6 Removal of the drain gutter

The drain gutter and the covering ring can be removed from the centrifugation area for cleaning.

See illustrations on page 2.

Removal of the grain gutter and the covering ring:

- Remove the covering ring (Fig. 3, c) from the centrifugation area.
- Carefully flap the inner sealing ring (Fig. 3, a) upwards and remove the drain gutter (Fig. 3, b) out of the centrifugation area.

Insertion of the drain gutter and covering ring:

- At the back of the centrifugation area carefully flap the sealing ring (Fig. 4, a) upwards and push the drain gutter (Fig. 4, b) under the sealing ring (Fig. 4, a) refer to Fig. 4.
The drain hole of the drain gutter (Fig. 4, b) must be located above the drain hole (Fig. 4, c) in the centrifugation area.
- Carefully flap the sealing ring (Fig. 4, a) upwards around the drain gutter on the inside and press the drain gutter (Fig. 4, b) downwards. The drain gutter (Fig. 3, b) must be located below the sealing ring (Fig. 3, a), refer to Fig. 3.
- Place the covering ring (Fig. 3, c) on the drain gutter in such a way (Fig. 3, b) that the inscription "**this side up**" can be read, refer to Fig. 3.

25 Faults

If the fault cannot be eliminated with the help of the fault table, please inform Customer Service.

Please specify the type of centrifuge and the serial number. Both numbers can be found on the name plate of the centrifuge.

25.1 Operating faults

Fault	Cause	Remedy
Incomplete wash cycle	<ul style="list-style-type: none"> – A 24-fold rotor is being used, even though a 12-fold rotor has been set as the rotor type. – The filling volume is set too low. – The injection nozzles are blocked. 	<ul style="list-style-type: none"> – Check the ROTOR: parameter in the programme. – Check the SALINE (ml) parameter in the programme. – Clean the injection nozzles.
Pellet does not form on the base of the centrifuge container	<ul style="list-style-type: none"> – The revolutions per minute during the centrifugation run for the agglutination test are too low. – The tube holders are left hanging in the decantation position. 	<ul style="list-style-type: none"> – Check the rotor function.
No pellet/pellet too small	<ul style="list-style-type: none"> – A 12-fold rotor is being used, even though a 24-fold rotor has been set as the rotor type. – The filling volume is set too high. 	<ul style="list-style-type: none"> – Check the ROTOR: parameter in the program. – Check the SALINE (ml) parameter in the program.
The liquid is not decanted.	<ul style="list-style-type: none"> – The rotor mechanism is defective. 	<ul style="list-style-type: none"> – Check the rotor function.

25.2 Error messages

 Perform a MAINS RESET:



- Switch off the mains switch (switch position "0").
- Wait at least 10 seconds and then switch on the mains switch again (switch position "1").

Message / fault		Cause	Remedy
No display	---	<ul style="list-style-type: none"> - No voltage. - Overvoltage protection tripped out. - Fuse blown. 	<ul style="list-style-type: none"> - Check supply voltage. - Check the fuse, see chapter "Changing the fuse". - Mains switch ON.
TACHO - ERROR	01	- Faulty speedometer.	<ul style="list-style-type: none"> - Open the lid. - Switch off the mains switch (switch position "0"). - Turn the rotor vigorously by hand. - Switch on the mains switch again (switch position "1"). The rotor must turn during switch-on.
	02	<ul style="list-style-type: none"> - No rotor installed. - Defective motor, frequency converter or drive. 	
IMBALANCE	---	Imbalance about motor axis through weight differential in rotor assembly.	<ul style="list-style-type: none"> - Open lid. - Correct imbalance.
CONTROL - ERROR	04, 06 - 09	Error in lid locking or lid closure.	<ul style="list-style-type: none"> - Open lid. - Perform a MAINS RESET.
N > MAX	05	Rotation too fast	
N < MIN	13	Rotation too slow	
ROTORCODE	10	Incorrect rotor coding	
POWER INTERRUPT	---	Power failure, centrifuging not properly completed	<ul style="list-style-type: none"> - Open lid. - Push (START WASH) button.
VERSIONS-ERROR	12	Mismatch between electronic components	<ul style="list-style-type: none"> - Open lid. - Perform a MAINS RESET.
CONTROL-ERROR	21 - 27	Error / defective control unit	
SER I/O - ERROR	30 - 38	Error / defective interface	
FU / CCI - ERROR	60 - 83	Error / defective motor control	
LOW SALINE	---	<ul style="list-style-type: none"> - The physiological saline solution container is empty. - The input hose for the physiological saline solution is kinked, bent or broken. - Hose pump, flow rate measurer defective. 	<ul style="list-style-type: none"> - Open the lid or MAINS RESET - Replace the physiological saline solution container. - Clean the input pipe for the saline solution and straighten any kinks.
N > ROTOR-MAX	---	Speed in the selected program greater than the maximum speed of the rotor.	<ul style="list-style-type: none"> - Check the set speed. - Reduce the set speed

26 Changing the fuse


Switch off the mains switch and separate the centrifuge from the mains!

The fuse holder (A) with the fuse is next to the mains switch.

- Remove the connecting cable from the machine plug socket.
- Screw the screw cap (B) of the fuse holder (A) off by turning it counter clockwise and pull it out with the fuse.
- Exchange the defective fuse.

Only T2.5 AH/250V, 6.3 x 32 mm fuses with - and -approval (order no. E2268) may be used.

- Tightly screw the screw cap (B) of the fuse holder (A) back on by turning in the clockwise direction.
- Reconnect the centrifuge to the mains supply.

27 Acceptance of the centrifuges for repair

If the centrifuge is returned to the manufacturer for repair, it must be decontaminated and cleaned to protect persons, environment and material.

We reserve the right to accept contaminated centrifuges.

Costs incurred for cleaning and disinfection are to be charged to the customer.

We ask for your understanding in this matter.

28 Disposal

Before disposal, the device must be decontaminated and cleaned to protect people, the environment and property.

When you are disposing of the device, the respective statutory rules must be observed.

Pursuant to guideline 2002/96/EC (WEEE), all devices supplied after August 13, 2005 may not be disposed as part of domestic waste. The device belongs to group 8 (medical devices) and is categorized in the business-to-business field.

The icon of the crossed-out trash can shows that the device may not be disposed as part of domestic waste.

The waste disposal guidelines of the individual EC countries might vary. If necessary, contact your supplier.

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