

**Heraeus®  
Megafuge® 11  
Centrifuge Series**



**Outstanding Capacity, Convenience  
and Value for Routine Applications**

Introducing the Heraeus Megafuge 11 Centrifuge Series – a new Thermo Scientific centrifuge offering unsurpassed convenience and exceptional versatility for routine applications. With its quiet performance and compact footprint, the Heraeus Megafuge 11 Centrifuge Series is the ideal choice for high-capacity, high-throughput processing in the clinical lab.

*Featuring Fast Rotor  
Exchanges with AUTO-LOCK®*

**DJB Labcare - The UK's Leading Centrifuge Specialist**  
info@djblabcare.co.uk www.djblabcare.co.uk  
01908 612598

# Heraeus Megafuge 11 Centrifuge Series



## Capacity, Convenience and Value for Routine Applications

The Thermo Scientific Heraeus Megafuge 11 Centrifuge Series is an excellent choice for high-capacity, high-throughput processing in the lab. Offering outstanding reliability, the Heraeus Megafuge 11 Centrifuge Series delivers consistent separations quickly, with maximum convenience and comfort.

### Exceptional Capacity

The Heraeus Megafuge 11 Series provides high capacity in a compact footprint, so you can process large numbers of samples without giving up precious bench space.

- Process 60 x 5/7 ml vacutainer or 16 x 50 ml conical tubes with the versatile T41 swing-out rotor.
- Process 28 x 5/7 ml vacutainer or 4 x 50 ml conical tubes in sealed conditions with the S41 rotor.
- The AC15.4 fixed-angle rotor accommodates 30 x 15 ml tubes.
- The T20 microplate rotor spins 2 x 3 microplates in sealed buckets.

- The FIBERLite® 48 x 2 ml microtube rotor enables high-speed microvolume centrifugation at > 23,000 x g.

### Convenient Operation

- The patented AUTO-LOCK® rotor exchange system on the Heraeus Megafuge 11 Series lets you swap rotors in seconds. You can change protocols quickly and perform daily cleaning tasks with ease.
- The user-friendly control panel gives you one-touch access to frequently-used programs.

### Proven Reliability

- The Heraeus Megafuge 11 Series' brushless induction motor provides years of service with minimal maintenance.
- Whisper-quiet operation eliminates distractions.

### Proven Safety

- The Heraeus Megafuge 11 Series supports a protected laboratory environment with a full set of sealed accessories, tested by CAMR\* in Porton Down, UK, the gold standard in biological safety testing.

\*CAMR is now incorporated within the UK Health Protection Agency (HPA)



### Quick Recall Keys

provide single keystroke access to frequently used programs.

### G-Force Control

allows protocol transfer from other centrifuges without cross-reference to RPM.

### Bright LED Displays

for easy, at-a-glance monitoring of set and run conditions across the lab.

### Pre-Cooling Function

on the refrigerated Megafuge 11R for rotor and chamber preparation.

DJB Labcare - The UK's Leading Centrifuge Specialist  
info@djblabcare.co.uk www.djblabcare.co.uk 0800 612598

## Heraeus Megafuge 11 Centrifuge Series Rotors and Accessories

### T41 High-Capacity Swing-Out Rotor



**T41 High-Capacity Swing-Out Rotor (11210435)**  
280 ml standard bucket (11210436) with sealing lid (11175712) and variety of inserts to meet your specific application requirements

Description	Rotor Capacity	Max Tube ø (mm)	Radius (mm)	Max Speed (rpm)	Max RCF (xg)	Cat. No.
<b>T41 Swing-Out Rotor, 4 X 280 ml</b>						<b>11210435</b>
Standard buckets, 280 ml (set of 4)		77	164	4100	3082	11210436
Buckets, 4 X 50 ml conical (set of 4)	<b>16</b>	33	169	4100	3176	11175752
Buckets, 3 X 50 ml conical skirted (set of 4)	<b>12</b>	33	169	4100	3176	11175759
Buckets, 15 X 5/7 ml vacutainer (set of 4)	<b>60</b>	13.5	165	4100	3101	11175763
Sealing lids for 11210436 (set of 4)						11175712
<b>Adapters for T41 Swing-Out Rotor (set of 4)</b>						
1 X 280 ml round	4	64	161	4100	3026	11175713
1 X 250 ml flat	4	61	161	4100	3026	11175714
1 X 200 ml flat	4	55	156	4100	2932	11175715
1 X 175 ml conical	4	60	166	4100	3120	11175716
1 X 100 ml	4	38	161	4100	3026	11175719
2 X 65 ml round	8	34	156	4100	2932	11175720
4 X 50 ml round	16	28.5	132	4100	2481	11175721
1 X 50 ml conical skirted	4	29	163	4100	3063	85540175
<b>2 X 50 ml conical</b>	<b>8</b>	<b>29</b>	<b>134</b>	<b>4100</b>	<b>2518</b>	<b>11175722</b>
<b>3 X 50 ml conical*</b>	<b>12</b>	<b>29</b>	<b>134</b>	<b>4100</b>	<b>2518</b>	<b>11210348</b>
1 X Centriprep	4	28	160	4100	3007	11175742
5 X 30 ml round	20	25.5	157	4100	2951	11175723
3 x 25 ml Universal	12	26	153	4100	2875	11204125
5 X 25 ml	20	24	155	4100	2913	11175724
4 X 25 ml Corex	16	24	155	4100	2913	11175725
9 X 16 ml round	36	18	160	4100	3007	11175726
<b>12 X 10 ml vacutainer</b>	<b>48</b>	<b>17</b>	<b>160</b>	<b>4100</b>	<b>3007</b>	<b>11175727</b>
<b>4 X 15 ml vacutainer *</b>	<b>16</b>	<b>16</b>	<b>163</b>	<b>4100</b>	<b>3063</b>	<b>11175738</b>
7 X 15 ml urine	28	17	161	4100	3026	11175729
<b>4 X 15 ml conical</b>	<b>16</b>	<b>17</b>	<b>163</b>	<b>4100</b>	<b>3063</b>	<b>11175730</b>
<b>5 X 15 ml conical *</b>	<b>20</b>	<b>17</b>	<b>163</b>	<b>4100</b>	<b>3063</b>	<b>11175731</b>
4 X 14 ml Corning round	16	15.5	163	4100	3063	11175747
12 X 13 ml	48	16	160	4100	3007	11175732
9 X 10 ml	36	16	160	4100	3007	11175728
19 X 8 ml	76	12	161	4100	3026	11175733
<b>12 X 5/7 ml vacutainer</b>	<b>48</b>	<b>12.5</b>	<b>160</b>	<b>4100</b>	<b>3007</b>	<b>11175740</b>
19 X 5 ml RIA	76	13	161	4100	3026	11175734
19 X 3 ml RIA	76	11	161	4100	3026	11175735
13 X 2 ml	52	10	159	4100	2988	11175736

\* In open buckets only.

### S41 Swing-Out Rotor



**S41 Swing-Out Rotor (11210246)**  
200 ml standard bucket (11175611) with sealing lid (11175612) and adapters for variety of applications

Description	Rotor Capacity	Max Tube ø (mm)	Radius (mm)	Max Speed (rpm)	Max RCF (xg)	Cat. No.
<b>S41 Swing-Out Rotor, 4 X 200 ml</b>						<b>11210246</b>
Standard buckets, 200 ml (set of 4)		57	161	4100	3026	11175611
Sealing lids for 11175611 (set of 4)						11175612
<b>Adapters for S41 Swing-Out Rotor (set of 4)</b>						
200 ml PP bottle p/n 41193271	4	56.5	161	4100	3026	direct
1 X 100 ml	4	45	159	4100	2990	11174181
1 X 50 ml round	4	35	158	4100	2970	11174591
<b>1 X 50 ml conical</b>	<b>4</b>	<b>29</b>	<b>159</b>	<b>4100</b>	<b>2990</b>	<b>11174528</b>
1 X 25 ml universal flat/skirt	4	25	158	4100	2970	11174527
3 X 20 ml	12	22	159	4100	2990	11174502
<b>3 X 15 ml conical</b>	<b>12</b>	<b>17.5</b>	<b>156</b>	<b>4100</b>	<b>2930</b>	<b>11174537</b>
7 X 15 ml Nalgene	28	17	156	4100	2930	11174183
3 x 15 ml urine	12	17	159	4100	2990	11210419
<b>3 x 15 ml vacutainer</b>	<b>12</b>	<b>17.5</b>	<b>159</b>	<b>4100</b>	<b>2990</b>	<b>11210523</b>
3 X 12 ml Greiner flat	12	17.5	159	4100	2990	11174509
<b>5 X 10 ml vacutainer</b>	<b>20</b>	<b>16</b>	<b>159</b>	<b>4100</b>	<b>2990</b>	<b>11174503</b>
5 X 5/10 ml flat	20	17	159	4100	2990	11174534
<b>7 X 5/7 ml vacutainer</b>	<b>28</b>	<b>13</b>	<b>159</b>	<b>4100</b>	<b>2990</b>	<b>11210604</b>
<b>9 X 5 ml vacutainer</b>	<b>36</b>	<b>13.5</b>	<b>159</b>	<b>4100</b>	<b>2990</b>	<b>11174505</b>
5 X 5 ml flat	20	17	159	4100	2990	11174533
12 X 3-5 ml	48	12	158	4100	2970	11174506
6 X 2 ml	24	11	158	4100	2970	11174508



Indicates that the centrifuge accessory has been successfully tested for biological containment by C.A.M.R. (HPA) Porton Down, UK.

DJJ Labcare - The UK's Leading Centrifuge Specialist  
 info@djlabcare.co.uk www.djlabcare.co.uk 01908 625598



Sealed plate carrier  
(11175631)



### T20 Microplate Rotor

Description	Rotor Capacity	Radius (mm)	Max Speed (rpm)	Max RCF (xg)	Cat. No.
<b>T20 Microplate Rotor</b>	<b>2x3 Microplates / 2x1 DeepWell</b>	<b>115</b>	<b>4100</b>	<b>2160</b>	<b>11175750</b>
Sealed plate carriers (set of 2)	2x3 Microplates / 2x1 DeepWell				11175631
Rubber cushion for microplates					11174207
Flask adapters (set of 2)	2x1 T-25 or T-75 Falcon® type flasks				11210762

### AC15.4 Fixed Angle Rotor

Description	Rotor Capacity	Max Tube ø (mm)	Radius (mm)	Max Speed (rpm)	Max RCF (xg)	Cat. No.
<b>AC15.4 Fixed Angle Rotor 37°</b>	<b>30 x 15 ml</b>	<b>17.5</b>	<b>135</b>	<b>4100</b>	<b>2537</b>	<b>11175755</b>
<b>Adapters (each)</b>						
1 x 10/15 ml	30		135	4100	2537	direct
1 X 13.5 ml	30		130	4100	2440	11203666
1 X 8/7 ml vacutainer	30		130	4100	2440	11172596
1 X 5 ml Hemogard	30		116	4100	2180	11172595
1 X 5 ml	30		116	4100	2180	11172402
1 X 3 ml	30		112	4100	2100	11172287
1 X 2/1.5 ml	30		105	4100	1970	11172288

### AC100.10A Fixed Angle Rotor

Description	Rotor Capacity	Max Tube ø (mm)	Radius (mm)	Max Speed (rpm)	Max RCF (xg)	Cat. No.
<b>AC100.10A Fixed Angle Rotor 25°</b>	<b>6 x 100 ml</b>	<b>38</b>	<b>99</b>	<b>10000</b>	<b>11066</b>	<b>11210437</b>
<b>Adapters (set of 6)</b>						
1 X 50 ml round	6	29	93	10000	10397	11174713
1 X 50 ml conical	6	30	93	10000	10397	11174714
1 X 38 ml round	6	25.5	88	10000	9838	11174715
1 X 15 ml conical	6	17	93	10000	10397	11174716
2 X 10 ml round	12	16	83	10000	9279	11174717

### AC50.10 Fixed Angle Rotor

Description	Rotor Capacity	Max Tube ø (mm)	Radius (mm)	Max Speed (rpm)	Max RCF (xg)	Cat. No.
<b>AC50.10 Fixed Angle Rotor 40°</b>	<b>6 x 50 ml con</b>	<b>29.5</b>	<b>112</b>	<b>10000</b>	<b>12520</b>	<b>11175754</b>
<b>Adapters</b>						
1 X 50 ml conical (each)	6	30	112	10000	12520	direct
1 X 50 ml round (set of 6)	6	28.5	112	10000	12520	11210577
1 X 30/38 ml (set of 8)	6	26	91	10000	10174	11174599
1 X 15 ml conical (set of 4)	6	17.5	105	10000	11738	11174606
1 X 10 ml (set of 4)	6	17	112	10000	12520	11177378

### AC2.14 Sealed Microtube Rotor

Description	Rotor Capacity	Max Tube ø (mm)	Radius (mm)	Max Speed (rpm)	Max RCF (xg)	Cat. No.
<b>AC2.14 Sealed Rotor 45°</b>	<b>24 x 2 ml</b>	<b>11</b>	<b>84</b>	<b>14600</b>	<b>20016</b>	<b>11175741</b>
Spare lid for AC2.14 Rotor						11175642
<b>Adapters: see adapters for FIBERLite High-Throughput Microtube Rotor</b>						

### FIBERLite® High-Throughput Microtube Rotor

Description	Rotor Capacity	Max Tube ø (mm)	Radius (mm)	Max Speed (rpm)	Max RCF (xg)	Cat. No.
<b>FIBERLite Microtube Rotor 45°</b>	<b>48 x 2 ml</b>	<b>11</b>	<b>97</b>	<b>14600</b>	<b>23113</b>	<b>11210351</b>
<b>Adapters (set of 24)</b>						
1 x 2/1.5 ml	24 / 48	10				direct
1 x 0.5 / 0.8 ml	24 / 48	8				41174928
1 x 0.5 ml	24 / 48	6				41174938
1 x 0.2 ml PCR	24 / 48	6				11174631

### DC6.11 Drum Rotor

Description	Rotor Capacity	Max Tube ø (mm)	Radius (mm)	Max Speed (rpm)	Max RCF (xg)	Cat. No.
<b>DC6.11 Drum Rotor 90°</b>		<b>N/A</b>	<b>72</b>	<b>14600</b>	<b>17158</b>	<b>11175743</b>
<b>Adapters (one rack)</b>						
10 x 1.5-2 ml	60	11				11174561
20 x 0.5-0.8 ml Microtainer	120	8				11174573
20 x 0.7 ml Rhesus	120	6				11174574
21 x 0.6 ml Viatron	126	6				11174563
32 x 0.25-0.4 ml Beckman	192	6				11174562

DJB Labcare - The UK's Leading Centrifuge Specialist  
 info@djblabcare.co.uk www.djblabcare.co.uk  
 01908 612598



1. Snap in your swing-out rotor and lock it in place.



2. Instantly unfasten and remove it.



3. Snap in your microplate rotor, and you're ready to spin again – in less than 5 seconds!

## Maximum Safety and Greater Convenience with AUTO-LOCK

Our patented AUTO-LOCK rotor locking system makes routine processing safer, faster and easier:

- **Worry-free spins** – AUTO-LOCK ensures that your rotor is properly locked in place with minimum effort

- **Fast protocol changes** – You can swap rotors in seconds, without having to remove or tighten lug nuts

- **Easy cleaning** – AUTO-LOCK makes routine cleaning straightforward



Adapters and accessories are available for all commonly used sample tubes and vessels. Bucket sealing lids provide certified protection against aerosols, bacteria and viruses.



### High Capacity, Advanced Technology

The Heraeus Megafuge 11 Centrifuge Series offers a high-throughput, 48 x 2 ml carbon fiber rotor. This advanced material is easier to manage and offers a safer alternative to conventional metallic rotors, which can crack with age. Capable of supporting higher capacities and higher G-forces, the 48-place carbon fiber rotor also provides a high level of chemical resistance against corrosive agents, making it a truly cost-effective investment for your lab.



With its generous sample capacity and wide selection of rotors, the Heraeus Megafuge 11 Centrifuge Series is perfect for routine processing applications in the clinical lab.

DJB Labcare - The UK's Leading Centrifuge Specialist  
 info@djblabcare.co.uk www.djblabcare.co.uk 020 2008 612598

## Heraeus Megafuge 11 Centrifuge Series Specifications



Specifications	Megafuge 11	Megafuge 11R
<b>Max Capacity</b>		
Swing-out	4 x 280 ml	4 x 280 ml
Fixed Angle	6 x 100 ml	6 x 100 ml
<b>Max Density</b>	1.2	1.2
<b>Max Load</b>	1.34 kg	1.34 kg
<b>Max Speed</b>		
Swing-out	4100 rpm	4100 rpm
Fixed Angle	14600 rpm	14600 rpm
<b>Max RCF</b>		
Swing-out	3176 x g	3176 x g
Fixed Angle	23113 x g	23113 x g
<b>Control System</b>	Microprocessor	Microprocessor
<b>Programs</b>	4 + pre-cooling with direct recall	4 + pre-cooling with direct recall
<b>Speed Set/Display</b>		
Range	500 - 14600 rpm	500 - 14600 rpm
Step	10 rpm at less than 1000 rpm 100 rpm at above 1000 rpm	10 rpm at less than 1000 rpm 100 rpm at above 1000 rpm
Accuracy	±20 rpm (display)	±20 rpm (display)
<b>RCF Set/Display</b>		
Range	10 - 23100	10 - 23100
Step	10 x g at less than 1000 x g 100 x g at above 1000 x g	10 x g at less than 1000 x g 100 x g at above 1000 x g
Accuracy	±20 x g (display)	±20 x g (display)
<b>Radius Control Accuracy</b>	1 mm	1 mm
<b>Timer Set/Display</b>	30 sec - 99 min + continuous	30 sec - 99 min + continuous
<b>Temperature Set/Display</b>	N/A	-9°C to 40°C
<b>Acceleration Rates</b>	5	5
<b>Braking Rates</b>	5 + inertial	5 + inertial
<b>Drive System</b>	Direct	Direct
<b>Motor Type</b>	Brushless induction	Brushless induction
<b>Power (average / max)</b>	350 / 500W	550 / 800 W
<b>Sound level (6x50 ml at 10000 rpm)</b>	61 dBA	54 dBA
<b>Sound level (4x280 ml at 4100 rpm)</b>	62 dBA	59 dBA
<b>Standards</b>	IEC 1010-1, IEC 1010-2-020, CE marked	IEC 1010-1, IEC 1010-2-020, CE marked

Ordering Information	Megafuge 11	Megafuge 11R
<b>Catalog Numbers</b>		
230V - 50/60 Hz	75004410	75004412
120 V - 60 Hz	75004411	75004413
<b>Product Dimensions</b>		
Height / lid closed [cm / in]	37 / 14.6	37 / 14.6
Width [cm / in]	40 / 15.7	57 / 22.4
Depth [cm / in]	50 / 19.7	60 / 23.6
Weight [kg / lb]	40 / 88.1	72 / 158.7
<b>Packing Dimensions</b>		
Height [cm / in]	60 / 23.6	61 / 24.0
Width [cm / in]	61 / 24.0	76 / 29.9
Depth [cm / in]	54 / 21.2	74 / 29.1
Weight [kg / lb]	52 / 114.6	85 / 187.3

### Unrivalled Quality and Reliability

Thermo Scientific Heraeus centrifuges are renowned for their superior quality and outstanding reliability. From compact tabletop instruments to high-capacity floor models, Heraeus centrifuges provide years of safe, user-friendly operation, and deliver consistent, high-throughput separations.

### Worldwide Service and Support

We are committed to keeping your lab equipment working at peak performance levels. Our goal is to help you lower ownership costs, manage labs more effectively, and increase productivity. Contact your local account representative to learn more about our service offerings,

including service agreements, preventive maintenance, on-site field repair, depot repair, compliance services and educational services.

© 2007 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.