

# Cyto Accessories

## Cells – fundamental to life

Cytology makes a valuable contribution to medical diagnostics. High cell yield and good cell presentation are essential for the correct evaluation of cytological preparations.

The Hettich cyto system provides cytologists with a wide selection of accessories for preparing slides for the cytological examination of various body fluids.

## The assembly of a complete cyto insert is really easy



suspension 1660 for cyto chambers

suspension 1680 for angle chambers




two-tier adapter **5280** illustrated with slide carriers with fastening ring **1662** and cyto chambers **1663**



filter cards



one-funnel chambers



angle chambers



multi-funnel chambers



slide carrier **1662**  
slide carrier for 2 chambers **1670**  
slide carrier with fastening slider **1470**

**DJB Labcare - The UK's Leading Centrifuge Specialist**  
[info@djblabcare.co.uk](mailto:info@djblabcare.co.uk) [www.djblabcare.co.uk](http://www.djblabcare.co.uk) 01908 612598

## Cyto Chambers

### cyto chambers for slide carrier 1470

fitting into suspension 1660 and two-tier adapter 5280

cyto chambers		
capacity in ml	8	2 x 8
area in mm <sup>2</sup>	240	240
diameter in mm	17.5	17.5
<b>Cat. No.</b>	<b>1471</b>	<b>1475</b>

### cyto chambers for slide carrier 1662 and slide carrier for 2 chambers 1670

fitting into suspension 1660 and two-tier adapter 5280

cyto chambers						
capacity in ml	1	2	4	8	3 x 2	4 x 1
area in mm <sup>2</sup>	30	60	120	240	3 x 60	4 x 30
diameter in mm	6.2	8.7	12.4	17.5	3 x 8.7	4 x 6.2
<b>Cat. No.</b>	<b>1663</b>	<b>1664</b>	<b>1665</b>	<b>1666</b>	<b>1667</b>	<b>1668</b>
filter cards						
<b>Cat. No.</b>		<b>1675</b>		<b>1676</b>	<b>1677</b>	<b>1678</b>

## Angle Chambers

### angle chambers for slide carrier 1662

fitting into suspension 1680

angle chambers			
area in mm <sup>2</sup>	30	60	120
diameter in mm	6.2	8.7	12.4
<b>Cat. No.</b>	<b>1671</b>	<b>1672</b>	<b>1673</b>
filter cards			
<b>Cat. No.</b>	<b>1696</b>	<b>1697</b>	<b>1698</b>

You can find out how to sediment **liquor cerebrospinalis, urine, bronchial fluids, pleural fluid** and **ascites** on slides, and which accessories to use for obtaining best preparations by consulting Hettich's cyto application notes. These notes can be ordered free of charge or downloaded from our homepage [www.hettichlab.com](http://www.hettichlab.com).

### slide carriers



slide carrier with fastening slider  
**Cat. No. 1470**



slide carrier with fastening ring  
**Cat. No. 1662**



slide carrier for 2 chambers  
**Cat. No. 1670**

### suspensions



lid  
**Cat. No. 1661**



suspension  
**Cat. No. 1660**




suspension  
**Cat. No. 1680**



suspension  
**Cat. No. 1452**

### two-tier adapter



**Cat. No. 5280**

DJB Labcare - The UK's Leading Centrifuge Specialist  
[info@djblabcare.co.uk](mailto:info@djblabcare.co.uk) [www.djblabcare.co.uk](http://www.djblabcare.co.uk) 01908 612598

## Use of cyto chambers

slide carriers		suspension/adapter			rotor		centrifuge
Cat. No.	qty.	Cat. No.	lid	qty.	Cat. No.	places	
1470 1662 1670	4	1660	1661	4	1624	4	ROTOFIX 32 A UNIVERSAL 320/320 R
	6	1660	1661	6	1626	6	
	4	1452	-	4	1494	4	UNIVERSAL 320/320 R
	4	1452	-	4	1724	4	ROTINA 38/38 R
	6	1660	1661	6	1726	6	
	8	1660	1661	8	1748	8	
	8	5091 / 5280	5095	4	1798	4	
	8	5092 / 5280 / 2 x 5220	5095	4	5624	4	ROTANTA 460/460 R ROTANTA 460 S/460 RS
	8	5091 / 5280	5095	4	5694	4	
	12	5091 / 5280	5095	6	4446	6	
	12	4396/5280	-	6	4256	6	
12	5091 / 5280	5095	6	4296	6	ROTIXA 50 S/50 RS	

## Use of angle chambers

slide carriers		suspension/adapter			rotor		centrifuge
Cat. No.	qty.	Cat. No.	lid	qty.	Cat. No.	places	
1662	4	1680	incl.	4	1048	4	MIKRO 22 R
	4			1624	4	ROTOFIX 32 A UNIVERSAL 320/320 R	
	6			1626	6		
	8			1648	8		
	6			1726	6	ROTINA 38/38 R	
	8			1748	8		

Disposable Pasteur pipettes for removing the supernatant are available under Cat. No. **0552**.

Please note that RCF values in excess of 1,100 may lead to glass breakage.