

# Product Specifications

# ACCULAB

sartorius group



=====  
Balance meets requirements of Council  
Directive 89/336/EEC (CE marking; EMC)  
=====

Recommended accessories:

- Calibration weight with DKD\*  
certificate YCW6548-00
- Lock-down capability (for anti-theft  
locking device) LC-1
- Data printer, adjustable for use in  
legal metrology, with GLP/GMP  
function YDP03-0CE

\*DKD = German Calibration Service,  
recognized throughout Europe  
=====

## ATL-6201-I

### ACCULAB ATILON

Electronic precision balance with  
internal calibration at the press of a  
key and backlit amber display segments  
and symbols; choice of 4 filter levels  
for optimal weighing results even under  
adverse ambient conditions

=====  
Weighing capacity 6200 g  
Readability 0.1 g  
Tare range (subtractive) 6200 g  
Response time, typical 1.7 s  
Repeatability <=± 0.1 g  
Linearity <=± 0.1 g  
=====

Weighing pan 180x180 mm  
Balance housing (WxDxH) 215x305x97 mm  
Net weight approx. 3.7 kg  
=====

Standard features/equipment supplied:

- Built-in motorized calibration weight
- 25-contact bidirectional RS-232  
data interface
- ISO/GLP-compliant data records
- Below-balance weighing port
- Overload protection
- Anti-theft locking device
- Leveling feet; level indicator near  
display
- AC adapter
- Built-in application programs  
(mass unit conversion, counting  
with automatic reference sample  
updating, weighing in percent, net-  
totalizing, animal weighing,  
calculation by a factor  
(division/multiplication)
- Built-in clock