

Article Specification



SARTORIUS M-Power AZ6101

Entry level to professional weighing for applications in education. High quality weighing systems ensure stable and reliable measurement values. Precision balance with 7-segment display, stability indicator and 4 selectable filter levels, plus optimized display update rate

```
=====
Weighing capacity          6100 g
Readability                0.1 g
Tare range (subtractive)  6100 g
Repeatability              <=± 0.1 g
Linearity                  <=± 0.3 g
Calibration weight        5000 g
=====
```

```
=====
Pan size                   174x143 mm
Bal. housing (WxDxH)       200x275x75 mm
Net weight                 approx. 1.4 kg
=====
```

Standard features/equipment supplied:

- Bidirectional RS-232 data interface port
- AC adapter
- Antitheft locking device
- Overload protection
- Built-in application programs:
Net-total formulation, weighing in percent, counting, mass unit conversion by toggling, dynamic weighing (averaging weights based on a selected number of individual weight measurements: 5/10/20/50/100)

Optional accessories:

```
=====
Dust cover                 6960TE02
Calibration weight        YCW653-AC-00
Printer                    YDP20-0CE
External rechargeable  
battery set               YRB08Z
=====
```