

Article Specification

USB data interface

YADAP-USB



SARTORIUS M-Prove

AY5101

Smart balance in stackable design for daily use in schools. Meets the demands of ambitious school experiments, supported by high resolution, large weighing ranges and digital filtering.

```
=====
Weighing capacity          5100 g
Readability                 0.1 g
Tare range (subtractive)   5100 g
Repeatability              <=± 0.1 g
Linearity                   <=± 0.2 g
=====
```

```
=====
Weighing pan                142x130 mm
Balance housing (WxDxH)     178x248x80 mm
Net weight                   1,1 kg
=====
```

Standard equipment/features:

- Lock-down capability (anti-theft locking device)
- Overload protection
- Polished stainless steel weighing pan
- Removable Flip-up cover
- Built-in application programs:
Counting with variable reference sample, weighing in percent, display hold, totaling, specific gravity, under-pan weighing, mass unit conversion
- Prepared for data interface (YADAP-RS or YADAP-USB)
- AC adapter

Optional accessories:

```
=====
RS-232 data interface      YADAP-RS
```