

Certificate of Calibration

Customer: DJB Labcare Limited

Certificate Number: 0914001

Instrument Type: SSA-20 CO2 Analyser

Serial Number: 141002

Calibration Date: 16/09/2014

Calibration Next Due: 15/09/2015

Calibration Environment: 20°C +/- 3°C

Calibration & Testing Gases	Analyser Reading CO ₂
Atmospheric CO ₂	0.04%
Calibration Certified 5.00% CO ₂	4.99%
Calibration Certified 10.07% CO ₂	10.05%
Nitrogen 99.999%	0.00%

Comments:

The instrument was calibrated and tested using certified gases which are traceable to UK National Standards.

Calibration Technician: Shaun Dawson

Calibration Technician Signature: *Shaun Dawson*

