

Certificate number: 12082020

12 Howard Way Cromwell Business Centre Newport Pagnell Bucks MK16 9QS

> Tel: 01908 612598 Fax: 01908 217974

E-mail: service@djblabcare.co.uk Website: www.djblabcare.co.uk

Certificate of Calibration

Equipment: Fastime 3 Stopwatch Serial number: 018

Date: 12/08/2020 Calibration frequency: Every 12 months

Calibrated by: Jude Huckle

Results:

Set Value	Measured Value
05:00.00	05:00.03
10:00.00	10:00.00
60:00.00	60:00.04

Method of test:

All calibrations are carried out using an in house calibration work instruction WI40

Expanded uncertainty is 0.48s

Approved by: C Bradnam

