

# CERTIFICATE OF CALIBRATION

Certificate number: IN03049\_3/8863  
Date of calibration: 04 April 2012  
Product: G100  
Serial number: IN03049  
Calibration checked at: 22.1°C to 31.2°C

## Primary Gas Channels

Carbon Dioxide (CO2)	
Certified Gas (%)	Reading (%)
10.15	10.1
2.59	2.5
5.00	4.9

Barometer	
Certified (mb)	Reading (mb)
1006	1011

Approved by: Dawn Hemings

(Name)



(Signature)

All gases are traceable to certified National Standards

This unit must be serviced at regular 12 monthly intervals by a Thermo Fisher approved service facility.