

Spectrafuge™ Mini Centrifuges and Slide Spinner

Compact in design and personal sized, the *Spectrafuge Mini* is ideal for microfiltration and quick spin downs from the walls and caps of microcentrifuge tubes while the *Slide Spinner* is designed for drying microarray slides before scanning.

The *Spectrafuge Mini* is supplied with rotors and adapters to accommodate 1.5/2.0 mL, 0.5 mL and 0.4 mL tubes as well as 0.2 mL strips and individual tubes. The rotors are easily interchanged.

The *Slide Spinner* dries two slides in as little as 10 seconds for quality imaging and screening results. The slides are placed into cassettes that catch all liquid without splashing. Cassettes are reusable.

Both centrifuges feature a translucent lid that pivots on a durable stainless steel hinge pin. A switch on the side of the centrifuge starts and stops operation. Alternatively, with the switch in the "on" position, the centrifuge can be started and stopped by closing and opening the lid.

The *Spectrafuge Mini* is available in five different lid colors: gray, blue, teal, purple and red.

- For microtubes and strip tubes or slides
- Convenient on/off switch
- Small footprint
- Starts and stops in seconds



Spectrafuge Mini



Slide Spinner



SPECIFICATIONS

	<u>Spectrafuge Mini</u>	<u>Slide Spinner</u>
Maximum speed/RCF	6,000 rpm/2,000 x g ¹	4,800 rpm
Maximum capacity	2 x 0.2 mL strips or 6 x 1.5/2.0 mL	2 standard slides
Dimensions (W x D x H)	5.9 x 5.9 x 4.6 in/15 x 15 x 11.7 cm	5.9 x 5.9 x 4.6 in/15 x 15 x 11.7 cm
Weight	1 lb/0.45 kg	1 lb/0.45 kg
Electrical	230V~, 50 Hz or 120V~, 60 Hz	230V~, 50 Hz or 120V~, 60 Hz

¹Average max rcf in the strip tube rotor is 1,177 x g

ORDERING INFORMATION

Catalog No.	Description
C1301-230V*	<i>Spectrafuge Mini</i> complete with 1.5/2.0 mL rotor, strip tube rotor, 0.5 mL and 0.4 mL adapters, smoke gray lid, 230V
C1301-B-230V*	<i>Spectrafuge Mini</i> , as above, blue lid, 230V
C1301-T-230V*	<i>Spectrafuge Mini</i> , as above, teal lid, 230V
C1301-P-230V*	<i>Spectrafuge Mini</i> , as above, purple lid, 230V
C1301-R-230V*	<i>Spectrafuge Mini</i> , as above, red lid, 230V
C1300-RT	Replacement rotor for 6 x 1.5/2.0 mL tubes
C1300-RTS	Replacement rotor for 0.2 mL strip tubes or individual 0.2 mL tubes
C1205	Individual adapters for 0.5/0.6 mL tubes, pk of 6
C1206	Individual adapters for 0.4/0.25 mL tubes, pk of 6
C1222	Individual adapters for 0.2 mL thermal cycling tubes, pk of 6
C1303-T-230V*	<i>Slide Spinner</i> with two slide drying cassettes, 230V
C1303-SC	Extra cassettes for <i>Slide Spinner</i> centrifuge, pack of 2

*When ordering 230V units, please specify cord type by adding -EU or -UK to the end of the catalog number. 120V units are available. When ordering, please specify by removing -230V from catalog number.