

# Scanspeed



## Gyrospeed mini



Compact design with small footprint

- Supplied with an autoclavable fixed-angle 12 place rotor
- Separate "Pulse" key for quick run
- Rapid acceleration (< 12 sec) and deceleration (<16 sec)
- Automatic alarm system for over-heat and over-speed
- Automatic door release
- Optional : a rotor for 4 x 8-tube PCR strips



10-001-04-0007



### Optional Accessories



### Technical Specifications

10-001-04-0007	mini
max rpm	13,500 rpm
max rcf	12,300 x g
max capacity	1.5ml / 2.0ml x 12
run time	<30min
dimension mm (hwd)	145 x 208 x 245
weight Kg	4.4*

\* without rotor

  

10-001-04-0037	strip rotor
max rpm	6,000 rpm
max rcf	outer 1,850 / inner 1,542 x g
max capacity	0.2ml x 32

DJB Labcare - The UK's Leading Centrifuge Specialist  
info@djblabcare.co.uk www.djblabcare.co.uk 01908 612598