

Article Specification



(division/multiplication)

Optional accessories:
External rechargeable
battery set
Printer
Protective cover

YRB05Z
YDP20-OCE
YDC06-A

SARTORIUS M-Pact AX6202

The balance for all standard applications in the laboratory. Simple, safe application due to automatic calibration at the touch of a button. Top quality weighing system. Backlit 14-segment display, 2 tare keys, stability indicator and 4 selectable filter levels plus optimized display update rate complete the picture: a balance for the professional.

Weighing capacity	6200	g
Readability	0.01	g
Tare range (subtractive)	6200	g
Measuring time - typical	1.5	s
Repeatability	<=± 0.02	g
Linearity	<=± 0.04	g

Pan size	180x180	mm
Balance housing (WxDxH)	230x310x91	mm
Net weight	approx. 3.6	kg

Standard features/equipment supplied:

- Motorized calibration and adjustment function
- Bidirectional RS-232 data interface port
- ISO/GLP-compliant printing/recording of weighing and calibration data when a Sartorius printer is connected to the balance
- Below-balance weighing equipment
- Overload protection
- Lug for attaching an antitheft locking device
- Leveling feet, level indicator
- AC adapter
- Built-in application programs (mass unit conversion, counting with automatic reference sample updating, weighing in percent, net-total formulation, totalizing, animal weighing, calculation by a factor)